



ENDEAVOR II Pivotal Clinical Trial

Fact Sheet & Clinical Data Summary

Caution: The Endeavor™ Drug Eluting Coronary Stent is an investigational device with an investigational drug (ABT-578) and is exclusively for clinical investigation.

Trial Name: ENDEAVOR II Pivotal Clinical Trial

Start Date: July 2003

Enrollment Completion: January 2004

Purpose: ENDEAVOR II was a randomized (1 to 1), double-blind pivotal clinical trial to assess the safety and efficacy of the investigational Endeavor™ Drug Eluting Coronary Stent compared to the Medtronic Driver™ Cobalt Alloy Coronary Stent in the treatment of restenosis in single de novo native coronary artery lesions.

Size / Location:

- 1,197 patients at 72 sites in Europe, Asia Pacific, Israel, New Zealand and Australia.
- First and largest DES pivotal trial ever conducted outside the United States comparing a drug eluting stent to a bare metal stent.

Stent Diameters: 2.25 mm – 3.5 mm

Stent Lengths: 18 – 30 mm

Lesion Lengths: 14 – 27 mm / Pre Dilatation Required

Randomization: 1 to 1 – Endeavor Stent (600); Driver Stent Control (600)

Primary Endpoint: Target Vessel Failure (TVF; defined as a composite of cardiac death, myocardial infarction and Target Vessel Revascularization) at 9 months.)

Key Secondary Endpoints:

- MACE at 30 days
- Late Lumen Loss at 8 months
- In-Stent / In-Segment Angiographic Binary Restenosis (ABR) at 8 months
- Target Lesion Revascularization (TLR) at 9 months
- Target Vessel Revascularization (TVR) at 9 months

ENDEAVOR II Fact Sheet
Page 2 of 2

Follow Up: 600 patient angiographic follow up
300 patient Intravascular Ultrasound (IVUS)
IVUS for overlapping stents

ENDEAVOR II Key Clinical and Angiographic Results Summary

Parameter	Endeavor Stent (598 patients)	Driver Stent Control (599 patients)
MACE	7.4%	14.7%
TLR	4.6%	12.1%
TVR	5.7%	12.8%
TVF	8.1%	15.4%
%DS In-Segment	32.6%	44.3%
%DS In-Stent	27.9%	42.1%
ABR In-Segment	13.3%	34.2%
ABR In-Stent	9.5%	32.7%
LL In-Segment	0.36	0.71
LL In-Stent	0.62	1.03
LL Index In-Stent	0.34	0.54

###

Caution: The Endeavor™ Drug Eluting Coronary Stent is an investigational device with an investigational drug (ABT-578) and is exclusively for clinical investigation.