



# Cardiac Rhythm Disease Management Product Performance Report

*Important Patient Management Information for Physicians*

2009

First Edition – Issue 60



# A Message from the Vice President

Dear Customer,

At Medtronic, product quality and reliability have been and will continue to be a priority. For over 25 years, Medtronic has compiled and produced product performance reports with one primary goal, to provide you with the product information you need to best care for your patients.

Our commitment to you is best expressed in Medtronic's mission: "To strive without reserve for the greatest possible reliability and quality in our products; to be the unsurpassed standard of comparison and to be recognized as a company of dedication, honesty, integrity, and service." To this end, we continually explore new ways to expand, improve, and learn from our product performance systems and measures.

Our quality goals cannot be reached alone. We welcome your collaboration, insight, and recommendations. Please contact our Technical Services Department at 1 (800) 723-4636 with your feedback comments and any questions.

Your participation and assistance in returning explanted products are also critical. Returned products are tested and evaluated so that we can fully measure the performance of our devices. Please refer to the instructions on page 2 for assistance in returning products to the Medtronic CRDM Returned Product Analysis Laboratory.

As we constantly strive to exceed your expectations, we thank you for your dedication to improving and saving the lives of those suffering from cardiac rhythm disorders.

With appreciation and warm regards,



Tim Samsel  
Vice President, Quality and Regulatory  
Medtronic Cardiac Rhythm Disease Management  
Medtronic, Inc.

# Contact Information

We invite our customers to use these telephone numbers to call with suggestions, inquiries, or specific problems related to our products.

## US Technical Services Department

**Phone:** 1 (800) 723-4636 (Tachy)  
1 (800) 505-4636 (Brady)  
**Fax:** 1 (800) 824-2362  
[www.medtronic.com/corporate/contact.jsp](http://www.medtronic.com/corporate/contact.jsp)

## International Technical Centers

Europe (Heerlen NL) +31-45-566-8844  
Japan (Tokyo) +81-3-5753-4116

For questions related to this CRDM Product Performance Report, please call US Technical Services at the number above, or write to:

Timothy Smith  
Medtronic, Inc.  
8200 Coral Sea Street NE MS MVN61  
Mounds View, MN 55112 USA  
Email: [tim.smith@medtronic.com](mailto:tim.smith@medtronic.com)

For questions related to returning explanted product or returning product that shows signs of malfunction, please contact:

*Outside the United States:*  
Your Medtronic representative or international technical center at the number above.

*Within the United States:*  
Your Medtronic representative or

CRDM Returned Product Analysis Laboratory  
Medtronic, Inc.  
7000 Central Avenue NE MS RCE172  
Minneapolis, MN 55432-3576 USA

Phone: 1 (800) 328-2518, ext. 44800  
Email: [crdm.returnedproduct@medtronic.com](mailto:crdm.returnedproduct@medtronic.com)

## Editorial Staff

### Independent Physician Quality Panel

Angelo Auricchio, MD, Lugano, Switzerland  
Hugh Calkins, MD, Baltimore, MD  
Steven J. Compton, MD, Anchorage, AK  
John P. DiMarco, MD, PhD, Charlottesville, VA  
Kevin Hackett, MD, Columbus, OH  
Mariell Jessup, MD, Philadelphia, PA  
R. Hardwin Mead, MD, Palo Alto, CA

### Editor

Tim Samsel, *Vice President, CRDM Quality and Regulatory*

### Authors

Timothy Smith, *Senior Principal Product Performance Engineer, CRDM, Product Performance Reporting*  
Corrine Buchanan, *Senior Clinical Trial Leader*  
Mike Lenarz, *Senior Statistician, CRDM*  
Scott McRae, *Statistician, CRDM*  
Tim Hamann, *Graphic Designer, CRDM*

### Medtronic Review Board

David Steinhaus, MD, *Vice President and Medical Director, CRDM*  
Lonny Stormo, *Vice President, CRDM, Therapy Delivery*  
Subu Mangipudi, *Director, Product Vigilance and Reliability*

## Trademarks of Medtronic, Inc.

Adapta®	InSync III	Quick Look™
AT500®	Protect™	Relia™
Attain®	InSync Sentry®	Secura™
CapSure®	Intrinsic®	SelectSecure®
CapSure Sense®	Jewel®	Sensia®
CapSureFix®	Kappa®	Sensing Assurance
Capture Management®	Legend®	Sigma®
CareLink®	Marquis®	Spectraflex
Concerto®	Maximo®	Sprint™
Consulta™	Medtronic CareAlert®	Sprint Fidelis®
EnPulse®	Medtronic CareLink®	Sprint Quattro®
EnRhythm®	Micro Jewel	Sprint Quattro Secure®
EnTrust®	Micro Minix	SureFix®
GEM®	Minix	Target Tip®
InSync®	Minuet	Tenax
InSync ICD®	MVP®	Thera®-i
InSync Marquis™	Onyx®	Transvene
InSync II Marquis™	Patient Alert™	Versa®
InSync III Marquis™	Preva	Virtuoso®
InSync Maximo®	Prevail®	
InSync II Protect™	Prodigy	

# CRDM Product Performance Report

2009 First Edition  
Issue 60

Date cutoff for this edition is  
January 31, 2009

This report is available online at  
[www.CRDMPPR.medtronic.com](http://www.CRDMPPR.medtronic.com)

Introduction 4  
Method for Estimating CRT, ICD, and IPG Device Performance 9

**CRT** **Cardiac Resynchronization Therapy** 13  
CRT Survival Summary 20  
CRT Reference Chart 22

**ICD** **Implantable Cardioverter Defibrillators** 23  
ICD Survival Summary 34  
ICD Reference Chart 37  
ICD Connector Styles 39

**IPG** **Implantable Pulse Generators** 40  
IPG Survival Summary 69  
IPG Reference Chart 77

## Leads

Method for Estimating Lead Performance 80

**Left-Heart Leads** 83  
Lead Survival Summary 86  
US Returned Product Analysis Summary 86  
Reference Chart 86

**Defibrillation Leads** 87  
Lead Survival Summary 95  
US Returned Product Analysis Summary 96  
Reference Chart 97

**Pacing Leads** 98  
Lead Survival Summary 130  
US Returned Product Analysis Summary 134  
Reference Chart 136

**Epi/Myocardial Pacing Leads** 138  
Lead Survival Summary 141  
US Returned Product Analysis Summary 142  
Reference Chart 142

**VDD Single Pass Pacing Leads** 143  
Lead Survival Summary 144  
US Returned Product Analysis Summary 144  
Reference Chart 144

**ICD and CRT-D Charge Time Performance** 145

**Advisories** 151

*(in order of communication date, from most recent to oldest)*

**Performance Notes** 160

Clinical Management of VCM near Elective Replacement 160  
Ensuring the Accuracy of Battery Longevity Estimates 161  
Interactions between Cardiac Pacing and Ventricular Arrhythmia Initiation 162  
AT500 Pacing System Follow-Up Protocol 163  
Insertion of the Lead into the Device 164  
GEM II DR/VR and GEM III DR/VR/AT ICD Battery Discharge Behavior 165  
General Follow-Up and Replacement of ICD Leads 166  
Clinical Management of High Voltage Lead System Oversensing 167  
Tests and Observations for Clinical Assessment of Chronic Pacing Leads 168

**Index** 169

CRT

ICD

IPG

Leads

ICD Charge Times

Advisories

Performance Notes

Index

# Introduction

All product performance reports are not created equal. For 26 years, Medtronic has monitored performance via both returned product analysis and multicenter clinical studies.

This Product Performance Report (PPR) presents device survival estimates, advisory summaries, performance notes, and other information pertinent to assessing the performance of Medtronic implantable pulse generators (IPGs), implantable cardioverter defibrillators (ICDs), cardiac resynchronization therapy (CRT) devices, and implantable pacing and defibrillation leads.

This Product Performance Report has been prepared in accordance with International Standard ISO 5841- 2:2000(E).

The survival estimates provided in this report are considered to be representative of worldwide performance.

## Survival Estimates

Medtronic Cardiac Rhythm Disease Management (CRDM) uses both returned product analysis and multicenter clinical studies to monitor performance.

Medtronic, like other companies, monitors CRT, ICD, and IPG device performance using returned product analysis. We also monitor CRT, ICD, and IPG device performance using an active multicenter clinical study. Medtronic CRDM is unique in the industry in that we track CRT, ICD, and IPG device survival using both methods.

Returned product analysis is a passive approach to assessing product performance. This approach provides a suitable measure of product performance only when a significant number of explanted products are returned to the manufacturer. Returned product analysis provides a measure of hardware performance, but not necessarily the total clinical performance (e.g., the incidence of complications such as infection, erosion, muscle stimulation, etc. are not estimated).

The survival estimates provided in this report for CRT, ICD, and IPG devices are based on returned product analysis. This approach is suitable because a significant number of explanted generators are returned for analysis.

Lead performance is monitored differently. In contrast to CRT, ICD, and IPG devices, a very small percentage of leads are returned to the manufacturer due to the difficulty of explanting them. For leads, an active clinical study provides more accurate survival estimates compared to estimates based solely on returned product analysis.

Survival estimates for leads are based on clinical observations recorded via Medtronic CRDM's System Longevity Study. This multicenter clinical study is

designed to record clinical observations representative of the total clinical experience. Therefore, the lead survival estimates include both lead hardware failure and lead-related medical complications, and do not differentiate a lead hardware failure from other clinical events such as exit block, perforation, dislodgement, or concurrent pulse generator failure.

The actuarial life table method is applied to the data collected for CRT, ICD, and IPG devices and leads to provide the survival estimates included in this report. A general introduction to understanding this method of survival analysis is given later in this introduction.

## ICD Charge Times

Since May 2000, Medtronic has provided important information on charge time performance of ICDs. The information provided in this report shows how ICD charge time can vary during the time a device is implanted. The information is presented in graphical format showing charge time as a function of implant time. The data for charge times are collected from devices enrolled in the System Longevity Study.

## Advisory Summaries

This Product Performance Report includes summaries of all advisories applicable to the performance of the products included in the report. An advisory is added to the report when any product affected by the advisory remains in service and at risk of experiencing the behavior described in the advisory. The advisory will remain in the report until Medtronic estimates no product affected by the advisory remains active, or the risk of experiencing the behavior described in the advisory has passed.

For most advisories, the products subject to the advisory retain essentially the same survival probability as the products of the same model(s) not affected by the advisory. For those advisories where the survival probabilities of the affected and non-affected populations do differ significantly, Medtronic will provide separate survival data for each population. The separate survival data will remain in the report until Medtronic estimates no affected product remains in active service.

## Performance Notes

This report concludes with a number of Performance Notes developed by Medtronic to provide additional product performance information relevant to follow-up practice and patient management.

continued

## Medtronic urges all physicians to return explanted products and to notify Medtronic when a product is no longer in use, regardless of reason for explant or removal from use.

### How You Can Help

Medtronic urges all physicians to return explanted products and to notify Medtronic when a product is no longer in use, regardless of the reason for explant or removal from use. The procedures for returning products vary by geographic location.

Mailer kits with prepaid US postage are available for use within the United States to send CRTs, ICDs, IPGs, and leads to Medtronic's CRDM Returned Product Analysis Lab. These mailers are sized to accommodate the devices and leads from a single patient or clinical event and are designed to meet postal regulations for mailing biohazard materials.

If the product being returned is located outside the United States, please contact your local Medtronic representative for instructions.

Medtronic also requests the return of explanted products from non-clinical sources, such as funeral homes, and will assume responsibility for storage and disposal of the product once received. For return of larger quantities of explanted products than the mailer can accommodate, Medtronic has handling and shipping guidelines available upon request.

Both mailers and guidelines can be requested by contacting the Returned Product Lab. For information on how to contact the Lab, refer to Contact Information on page 2 of this report.

We continually strive to improve this CRDM Product Performance Report. In keeping with this philosophy, we ask for your suggestions on the content and format of this report, as well as any information you have regarding the performance of Medtronic products. For information on how to comment on this report, see Contact Information on page 2 of this report.

### Overview of Survival Analysis

Medtronic uses the Cutler-Ederer actuarial life table method to estimate the length of time over which devices and leads will perform within performance limits established by Medtronic. This probability to perform within performance limits over time is called the *survival probability*.

Devices and leads are followed until an *event* occurs where the device or lead ceases to operate within performance limits. The length of time from implant to the event is recorded for each individual device and lead in the *population sample*. The population sample for CRT, ICD, and IPG devices is made up of patients whose devices are registered as implanted in the United States. For leads, the population sample is the patients enrolled in our multicenter, international prospective System Longevity Study.

For IPGs and ICDs, the events can be normal battery depletion or a device malfunction. For leads, the events are complications as defined for the study.

The actuarial life table method allows Medtronic to account for devices and leads removed from service for reasons unrelated to performance. Devices and leads removed for these reasons are said to be *suspended*. Examples include devices and leads:

- still in service at the time the analysis is performed
- removed to upgrade the device or lead
- no longer in service due to the death of the patient for reasons unrelated to the device or leads
- implanted in patients who are lost to follow-up

For each suspension, the device or lead has performed within performance limits for a period of time, after which its performance is unknown.

### An Example

The following example describes the survival analysis method used to establish the survival probability estimates for Medtronic CRDM devices and leads. The example is intended to provide an overview of the analysis process. The definitions of malfunctions and complications, and other details specific to calculating device and lead survival estimates, are provided in the articles *Method for Estimating CRT, ICD, and IPG Device Performance (page 9)* and *Method for Estimating Lead Performance (page 80)*.

continued

This simple example describes the survival analysis method used to establish the survival probability estimates for Medtronic CRDM devices and leads.

**Figure 1**

Implant times for devices of 16 patients. Gray bars with an orange X indicate devices removed from service due to an event. Blue bars indicate suspended devices.




Figure 1 illustrates 16 patients who have implanted devices. The first patient’s device (serial number BH004134) operated within performance limits for 32 months. At that time an event occurred. The fourth patient’s device (serial number BH008926) did not have an event but is suspended, perhaps because it was still in service at the time of the analysis. This patient had 66 months of implant experience. In this example, Figure 1 shows that seven of the 16 devices suffered events, and nine are suspended.

The first step in the life table method is to divide the implant time into intervals of a specific length. This example will use 12-month intervals.

The number of devices entered, suspended, and removed due to an event are counted and summarized, as shown in Table 1. For the first two intervals, all 16 devices survived and none were removed. In the interval (24-36 months), device BH004134 was removed due to an event. Therefore the table entries show that 16 entered the interval, none were suspended, and one was removed due to an event.

For the interval from 36-48 months, only 15 devices entered the interval and one was removed for an event. The remaining intervals are examined and the data entered in columns A, B, and C in like manner. The rest of the columns are filled in using calculations on the data in columns A, B, and C.

The *Effective Sample Size (D)* is the number of devices with full opportunity to experience a qualifying event in the interval. This is computed by subtracting one half the number suspended in the interval from the number that entered the interval. This calculation more accurately reflects the number of devices that could have experienced a qualifying event than simply using the number that entered the interval. Using the number of devices that enter an interval overestimates the sample size because the suspended devices do not complete the interval. Ignoring the suspended devices underestimates the sample size because suspended devices are not credited with their full service time. Using one half the number of suspended devices effectively splits the difference.

The next column in the table is the *Proportion with Event (E)*. This is the proportion of devices that had an event in the interval. It is calculated by dividing the *Number of Events (C)* by the *Effective Sample Size (D)*. The number can be interpreted as the estimated rate at which events occur in the time interval.

The *Interval Survival Probability (F)* is the estimate of probability of surviving to the end of the interval assuming the device was working at the beginning of the interval. It is calculated as 1 minus the *Proportion with Event (E)*. This number can be interpreted as the estimated rate at which events **do not** occur in the time interval.

continued

The Cumulative Survival Probabilities (G) from the last column of the life table can be plotted versus time intervals in the first column to give a survival curve. Figure 2 shows the survival curve for the data shown in Table 1.

**Table 1** Life Table for Figure 1

	A	B	C	D	E	F	G
Interval in Months	Number Entered	Number Suspended	Number of Events	Effective Sample Size	Proportion with Event	Interval Survival Probability	Cumulative Survival Probability
0	16	0	0	16	0.000	1.000	1.000
0-12	16	0	0	16	0.000	1.000	1.000
12-24	16	0	0	16	0.000	1.000	1.000
24-36	16	0	1	16	0.063	0.938	0.938
36-48	15	0	1	15	0.067	0.933	0.875
48-60	14	0	1	14	0.071	0.929	0.813
60-72	13	1	2	12.5	0.160	0.840	0.683
72-84	10	2	0	9	0.000	1.000	0.683
84-96	8	3	0	6.5	0.000	1.000	0.683
96-108	5	2	2	4	0.500	0.500	0.341
108-120	1	0	0	1	0.000	1.000	0.341
120-132	1	1	0	0.5	0.000	1.000	0.341

*Definitions:*

A	B	C	D	E	F	G
Number Entered	Number Suspended	Number of Events	Effective Sample Size	Proportion with Event	Interval Survival Probability	Cumulative Survival Probability
Number of devices active at the start of the interval	Number of devices removed from service for reasons other than an event	Number of units removed from service due to an event	Number of units with full opportunity to experience a qualifying event in the interval. Computed by subtracting one half the Number Suspended from the Number Entered.	Proportion of devices that had an event in the interval. Computed by dividing the Number of Events by the Effective Sample Size.	The probability of surviving to the end of the interval, assuming the device was working at the beginning of the interval. Computed as 1 minus the Proportion With Event.	The overall probability of surviving to the end of the interval. Computed by multiplying the Interval Survival Probability by the previous interval's Cumulative Survival Probability.


*Cumulative Survival Probability (G)* is the estimate of the unconditional probability of surviving to the end of the interval. It is computed by multiplying the *Interval Survival Probability (F)* by the previous interval's Cumulative Survival Probability. The probability of surviving to 132 months in the example is estimated for the table to be 0.341, or 34.1%.

The *Cumulative Survival Probabilities (G)* of the life table can be plotted versus time intervals in the first column to give a survival curve. Figure 2 shows the survival curve for the data in Table 1.

continued



**Figure 2** Survival Curve for Data Given in Table 1



**Confidence Intervals**

Since survival curves are based on a sample of the device and lead population, they are only estimates of survival. The larger the effective sample size, the more confident the estimate. A confidence interval can be calculated to assess the confidence in an estimate. In the Product Performance Report, Medtronic provides a 95% confidence interval. This can be interpreted as meaning that 95% of the time, the true survival of the device will fall somewhere in the interval.

**Survival Curves in the Product Performance Report**

Since the survival estimate can become very imprecise with small effective sample sizes, Medtronic truncates the survival curve when the effective sample size is less than 100 for CRTs, ICDs, and IPGs, and when the number entered is less than 50 for leads. The survival charts in the Product Performance Report show the effective sample size for each year interval where Medtronic has experience. When the effective sample size reaches 100 for CRTs, ICDs, and IPGs or when the number entered reaches 50 for leads, the next data point is added to the survival curve.

Although the report provides tabular data in one-year intervals, the curves are actually computed and plotted using 1-month intervals (for CRT, ICD, and IPG devices) or 3-month intervals (for leads).

A number of references are available for additional information on survival analysis using the Cutler-Ederer life table method.<sup>1</sup>

<sup>1</sup> Lee, Elisa T.(2003) Statistical Methods for Survival Data Analysis – 3rd Edition (Wiley Series in Probability and Statistics).

# Method for Estimating CRT, ICD, and IPG Device Performance

**Medtronic urges all physicians to return explanted products and to notify Medtronic when a product is no longer in use, regardless of reason for explant or removal from use.**

The performance of CRT, ICD, and IPG devices is expressed in terms of device survival estimates, where “survival” refers to the function of the device, not the survival of the patient. These survival estimates are intended to illustrate the probability that a device will survive for a given number of years with neither malfunction nor battery depletion.

The survival estimates are determined from the analysis of Medtronic CRDM’s United States device registration data and US returned product analysis data. These data are presented graphically and numerically.

Because this analysis is based on returned product analysis, the performance data does not reflect any device-related medical complications such as erosion, infection, muscle stimulation, or muscle inhibition.

## **Categorization of Depleted and Malfunctioning Devices for Survival Analysis**

For survival estimation, every device returned to Medtronic CRDM and analyzed in the CRDM Returned Product Analysis laboratory is assigned to one of three categories. The device 1) has functioned normally, 2) has reached normal battery depletion, or 3) has malfunctioned. This categorization is combined with data from our device registry for the total number of implants and the implant durations to create the survival curves presented on the following pages.

## **Definition of Malfunction**

Medtronic CRDM considers a device as having malfunctioned whenever the analysis shows that any parameter was outside the performance limits established by Medtronic while implanted and in service. To be considered a malfunction or battery depletion, the device must have been returned to Medtronic and analyzed.

Devices damaged after explant, damaged due to failure to heed warnings or contraindications in the labeling, or damaged due to interaction with other implanted devices (including leads) are not considered device malfunctions.

A device subject to a safety advisory is not considered to have malfunctioned unless it has been returned to Medtronic CRDM and found, through analysis, to actually have performed outside the performance limits established by Medtronic.

Not all malfunctions expose the patient to a loss of pacing or defibrillation therapy. Some malfunctions included in the following survival estimates may not have been detected at all by the physician or the patient. These malfunctions, however, are included in the survival estimates and provide important feedback to our product development organization.

To provide insight into the nature of malfunctions, each malfunction is categorized as Malfunction with Compromised Therapy Function or Malfunction without Compromised Therapy Function. A summary of these malfunctions is presented for the most recently market-released models.

For this report, Normal Battery Depletion, Malfunction with Compromised Therapy Function, and Malfunction without Compromised Therapy Function are defined as follows:

**Normal Battery Depletion** – The condition when:

- (a) a device is returned with no associated complaint and the device has reached its elective replacement indicator(s) with implant time that meets or exceeds the nominal (50 percentile) predicted longevity at default (labeled) settings, or
- (b) a device is returned and the device has reached its elective replacement indicator(s) with implant time exceeding 80% of the expected longevity calculated using the available device setting information.

Medtronic CRDM establishes expected longevity by statistically characterizing the power consumed by the device and the power available from the device battery. This characterization is applied to a number of parameter configurations to derive a statistical mean longevity value and standard deviation for each parameter configuration. The statistical mean value minus three standard deviations is used as the expected longevity for determining if a battery depleted normally.

continued

The Standard Actuarial Method is used to estimate IPG and ICD survival. This product performance report has been prepared in accordance with International Standard ISO 5841-2:2000(E).

For reference purposes, the following pages include estimated longevity for each model. The actual longevity achieved for any device while implanted will depend on the actual programmed parameters and patient factors, and may differ significantly from these estimates.

### Malfunction with Compromised Therapy Function

The condition when a device is found to have malfunctioned in a manner that compromised pacing or defibrillation therapy (including complete loss or partial degradation), while implanted and in service, as confirmed by returned product analysis.

*Examples:* Sudden loss of battery voltage; accelerated current drain such that low battery was not detected before loss of therapy; sudden malfunction during defibrillation therapy resulting in aborted delivery of therapy, intermittent malfunction where therapy is compromised while in the malfunction state.

### Malfunction without Compromised Therapy Function

The condition when a device is found to have malfunctioned in a manner that *did not* compromise pacing or defibrillation therapy, while implanted and in service, as confirmed by returned product analysis.

*Examples:* Error affecting diagnostic functions, telemetry function, data storage; malfunction of a component that causes battery to lose power quickly enough to cause premature battery depletion, but slowly enough that the condition is detected through normal follow-up before therapy is lost; mechanical problems with connector header that do not affect therapy.

### Expanded Malfunction Detail

The malfunctions are further divided into categories that identify the subject area of the malfunction. The malfunctions are divided into the following subject areas:

Electrical Component – Findings linked to electrical components such as integrated circuits, resistors, capacitors, diodes, etc.

Electrical Interconnect – Findings linked to the connections between electrical components such as wires, solder joints, wire bonds, etc.

Battery – Findings linked to the battery and its components

Software/Firmware – Findings linked to software or firmware function

Possible Early Battery Depletion – Findings where the actual reported implant time is less than 80% of the expected longevity calculated using the available device setting information with no device malfunction observed. There may not be sufficient device setting information to determine conclusively if battery depletion was normal or premature in the absence of a specific root cause finding. However, returned devices meeting the above criteria are conservatively classified as Possible Early Battery Depletion malfunctions.

Other – Findings related to other components such as insulators, grommets, setscrews, and packaging, and findings where analysis is inconclusive

### Returned Product Analysis Process

Analysis of returned product is performed according to written procedures. These procedures determine the minimum analysis required. The analysis required varies depending on the type of device, age of the device, the associated information received with the device, actual experience with models of similar design, and other factors. Additional analysis is performed as necessary to investigate a performance concern from a customer, or to collect specific reliability data.

When a device is returned with a performance concern from a customer, the general analysis process includes a preliminary analysis of the device in its as-received condition, followed by an automated functional test using test equipment equivalent to the equipment used in manufacturing.

When a malfunction is identified, failure analysis is performed to provide the detailed information necessary to investigate possible causes and actions. Medtronic CRDM maintains in-house expertise and performs its failure analysis using facilities it owns and supports. This capability permits detailed failure analysis.

continued

**Medtronic CRDM adjusts all-cause survival estimates to account for underreporting. While this lowers our all-cause survival estimates, we feel it gives a more accurate perspective on real performance.**

## Statistical Methods for Survival Analysis

Of the several different statistical methods available for survival analysis, the Standard Actuarial Method, with suspensions assumed distributed across the intervals (Cutler-Ederer Method), is used to determine estimates of IPG and ICD survival. This method is commonly used by medical researchers and clinicians.

Implant times are calculated from the implant date to the earlier of the explant date or the cutoff date of the report. From this data an estimate of the probability of device survival is calculated at each monthly interval.

On the following pages, each graph includes a survival curve where events include malfunctions and normal battery depletions. This survival curve is a good representation of the probability a device will survive a period of time without malfunction and without battery depletion. For example, if a device survival probability is 95% after 5 years of service, then the device has a 5% chance of being removed due to battery depletion or malfunction in the first 5 years following implant.

In addition, a second curve is included to show survival excluding normal battery depletion. This curve is a good representation of the probability for a device to survive without malfunction. This curve includes only malfunctions as events and excludes normal battery depletion.

Since the survival estimate can become very imprecise with small effective sample sizes, Medtronic truncates the survival curve when the effective sample size is less than 100 for CRT, ICD, and IPG devices. The survival charts in the Product Performance Report show the effective sample size for each year interval where we have experience. When the effective sample size reaches 100, the next data point is added to the survival curve.

Although the report provides tabular data in one-year intervals, the curves are actually computed and plotted using one-month intervals.

The data in the tables are rounded to the nearest tenth of one percent. Occasionally, a graph may show 100% survival, but have one or more malfunctions or battery depletions. This occurs because, even with the malfunctions or battery depletions, the data rounds to 100%.

The survival curves are statistical estimates. As performance experience accumulates, the estimation improves. Confidence intervals are provided as a way

to indicate the degree of certainty of the estimates. Greenwood's formula is used to calculate corresponding 95% confidence intervals for the standard errors, and the complementary log-log method is used to produce the confidence bounds.

## Sample Size and How the Population and Population Samples Are Defined

The population sample from which the survival estimates are derived is comprised of the devices registered as implanted in the United States as of the report cutoff date. The number of registered implants, as well as an estimate of the number that remain in active service, is listed for each model. To be included in the population, the device must have been registered with Medtronic's registration system and implanted for at least one day.

This sample based on US implants is considered to be representative of the worldwide population, and therefore the survival estimates shown in this report should be representative of the performance worldwide of these models.

A CRT, ICD, or IPG model or model family will be included in this report when it has accumulated at least 10,000 implant months and will remain in the report as long as at least 500 devices remain active.

## Methods Used to Adjust for Underreporting of Malfunction and Battery Depletion

The tables on the following pages show the actual number of malfunctions and battery depletions recorded by the analysis lab for US registered devices. Since not all devices are returned to Medtronic CRDM for analysis, these numbers underestimate the true number of malfunctions and battery depletions. To more accurately estimate the all-cause device survival probabilities, the number of malfunctions and battery depletions used to plot each interval of the all-cause survival curves is adjusted (multiplied) by a factor that is based on an estimate of the magnitude of underreporting. The magnitude of underreporting is estimated by analyzing experience in clinical studies (including the System Longevity Study) and the device registration system.

continued

However, at this time, no adjustment for underreporting is applied to the malfunction-free survival curve because a method for estimating malfunction-only underreporting has not been developed.

### **Adjustments to Registered Implants to Compensate for Unreported Devices Removed from Service**


Devices are at times removed from service for reasons other than device malfunction or battery depletion. Examples are devices removed from service due to non-device related patient mortality and devices removed due to changes in the patient's medical condition. Because an accurate estimate of device survival depends on an accurate estimate of the number of devices in service, it is important not to overstate the number of devices in service.

To ensure the number of devices in service is not overstated, the patient mortality rate derived from our device registration system is monitored and compared to published mortality rates for comparable patient populations. If, during calculation of the survival curves, the patient mortality indicated by the data in our device registration is significantly different from published rates, an adjustment is applied to correct the difference.

**7272 InSync ICD**

Product Characteristics


US Market Release	Jul-02	Malfunctions (US)	238	NBD Code	VVED
Registered US Implants	13,000	<b>Therapy Function Not Compromised</b>	223	Serial Number Prefix	PJP
Estimated Active US Implants	1,000	Battery	2	Max Delivered Energy	34 J
Normal Battery Depletions (US)	1,284	Electrical Component	32	Estimated Longevity	<a href="#">See page 22</a>
Advisories	None	Software/Firmware	5		
		Possible Early Battery Depletion	184		
		<b>Therapy Function Compromised</b>	15		
		Battery	1		
		Electrical Component	14		



**7277 InSync Marquis**

Product Characteristics


US Market Release	Mar-03	Malfunctions (US)	74	NBD Code	VVED
Registered US Implants	7,000	<b>Therapy Function Not Compromised</b>	63	Serial Number Prefix	PLT
Estimated Active US Implants	100	Battery	1	Max Delivered Energy	30 J
Normal Battery Depletions (US)	588	Electrical Component	8	Estimated Longevity	<a href="#">See page 22</a>
<b>Advisories:</b> <a href="#">See page 153</a> – 2005 Potential Premature Battery Depletion Due to Battery Short		Software/Firmware	1		
		Possible Early Battery Depletion	53		
		<b>Therapy Function Compromised</b>	11		
		Battery ( <i>10 malfunctions related to advisory</i> )	10		
		Electrical Component	1		



**7289 InSync II Marquis**

Product Characteristics


US Market Release	Jul-03	Malfunctions (US)	284	NBD Code	VVED
Registered US Implants	28,000	<b>Therapy Function Not Compromised</b>	253	Serial Number Prefix	PRJ
Estimated Active US Implants	1,000	Electrical Component	20	Max Delivered Energy	30 J
Normal Battery Depletions (US)	5,001	Software/Firmware	1	Estimated Longevity	<a href="#">See page 22</a>
<b>Advisories:</b> <a href="#">See page 153</a> – 2005 Potential Premature Battery Depletion Due to Battery Short		Possible Early Battery Depletion	232		
		<b>Therapy Function Compromised</b>	31		
		Battery (8 malfunctions related to advisory)	10		
		Electrical Component	21		



**7297 InSync Sentry**

Product Characteristics

US Market Release	Nov-04	Malfunctions (US)	30	NBD Code	VVED
Registered US Implants	9,000	<b>Therapy Function Not Compromised</b>	29	Serial Number Prefix	PRK
Estimated Active US Implants	4,000	Battery	1	Max Delivered Energy	35 J
Normal Battery Depletions (US)	458	Electrical Component	6	Estimated Longevity	<a href="#">See page 22</a>
Advisories	None	Software/Firmware	1		
		Possible Early Battery Depletion	21		
		<b>Therapy Function Compromised</b>	1		
		Electrical Component	1		






**7299 InSync Sentry**

Product Characteristics


US Market Release	Apr-05	Malfunctions (US)	48	NBD Code	VVED
Registered US Implants	31,000	<b>Therapy Function Not Compromised</b>	43	Serial Number Prefix	PRK
Estimated Active US Implants	18,000	Electrical Component	9	Max Delivered Energy	35 J
Normal Battery Depletions (US)	708	Software/Firmware	2	Estimated Longevity	<a href="#">See page 22</a>
Advisories	None	Possible Early Battery Depletion	32		
		<b>Therapy Function Compromised</b>	5		
		Electrical Component	5		



**7303 InSync Maximo**

Product Characteristics


US Market Release	Jun-04	Malfunctions (US)	62	NBD Code	VVED
Registered US Implants	17,000	<b>Therapy Function Not Compromised</b>	57	Serial Number Prefix	PRL
Estimated Active US Implants	6,000	Electrical Component	11	Max Delivered Energy	35 J
Normal Battery Depletions (US)	1,408	Software/Firmware	2	Estimated Longevity	<a href="#">See page 22</a>
Advisories	None	Possible Early Battery Depletion	44		
		<b>Therapy Function Compromised</b>	5		
		Electrical Component	5		



**7304 InSync Maximo**

Product Characteristics


US Market Release	Apr-05	Malfunctions (US)	24	NBD Code	VVED
Registered US Implants	18,000	<b>Therapy Function Not Compromised</b>	22	Serial Number Prefix	PRL
Estimated Active US Implants	11,000	Battery	1	Max Delivered Energy	35 J
Normal Battery Depletions (US)	344	Electrical Component	7	Estimated Longevity	<a href="#">See page 22</a>
Advisories	None	Possible Early Battery Depletion	14		
		<b>Therapy Function Compromised</b>	2		
		Electrical Component	2		



**8040 InSync**

Product Characteristics


US Market Release	Aug-01	Malfunctions (US)	28	NBD Code	DDDR
Registered US Implants	15,000	<b>Therapy Function Not Compromised</b>	7	Serial Number Prefix	PIN
Estimated Active US Implants	3,000	Electrical Component	4	Estimated Longevity	<a href="#">See page 22</a>
Normal Battery Depletions (US)	533	Possible Early Battery Depletion	3		
Advisories	None	<b>Therapy Function Compromised</b>	21		
		Electrical Interconnect	21		



**8042 InSync III**

Product Characteristics


US Market Release	Feb-03	Malfunctions (US)	6	NBD Code	DDDR
Registered US Implants	29,000	<b>Therapy Function Not Compromised</b>	3	Serial Number Prefix	PKF
Estimated Active US Implants	16,000	Electrical Component	2	Estimated Longevity	<a href="#">See page 22</a>
Normal Battery Depletions (US)	243	Possible Early Battery Depletion	1		
Advisories	None	<b>Therapy Function Compromised</b>	3		
		Electrical Interconnect	3		



**C154DWK, C164AWK, C174AWK Concerto**

Product Characteristics


US Market Release	May-06	Malfunctions (US)	50	NBD Code	VVED
Registered US Implants	68,000	<b>Therapy Function Not Compromised</b>	36	Serial Number Prefix	PVU, PVT, PVR
Estimated Active US Implants	57,000	Electrical Component	7	Max Delivered Energy	35 J
Normal Battery Depletions (US)	168	Possible Early Battery Depletion	29	Estimated Longevity	<a href="#">See page 22</a>
Advisories	None	<b>Therapy Function Compromised</b>	14		
		Electrical Component	13		
		Electrical Interconnect	1		



**D224TRK Consulta CRT-D**

Product Characteristics


US Market Release	Aug-08	Malfunctions (US)	0	NBD Code	DDED
Registered US Implants	2,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PUD
Estimated Active US Implants	2,000	<b>Therapy Function Compromised</b>	0	Max Delivered Energy	35 J
Normal Battery Depletions (US)	0			Estimated Longevity	<a href="#">See page 22</a>
Advisories	None				



**D284TRK Maximo II CRT-D**

Product Characteristics

US Market Release	Mar-08	Malfunctions (US)	0	NBD Code	VVED
Registered US Implants	1,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PZP
Estimated Active US Implants	1,000	<b>Therapy Function Compromised</b>	0	Max Delivered Energy	35 J
Normal Battery Depletions (US)	0			Estimated Longevity	<a href="#">See page 22</a>
Advisories	None				



**Device Survival Summary** (95% Confidence Interval)  
 The following table shows CRT device survival estimates with 95% confidence intervals. Estimates are shown both with and without normal battery depletions included.

Model Number	Family	US Market Release	Registered US Implants	Estimated Active US Implants	Normal Battery Depletions (US)	Malfunctions (US)		Device Survival Probability (%)									
						Therapy Function Compromised	Therapy Function Not Compromised	Years After Implant									
						1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr			
<b>7272</b>	InSync ICD	Jul-02	13,000	1,000	1,284	15 + 223 = 238	Total = 238	99.9 +0.0/-0.1	99.8 +0.1/-0.1	99.1 +0.2/-0.2	97.8 +0.3/-0.4	96.8 +0.4/-0.5	95.5 +0.6/-0.7 at 71 mo				
								99.7 +0.1/-0.1	99.2 +0.2/-0.2	95.4 +0.4/-0.5	87.6 +0.8/-0.8	66.7 +1.3/-1.3	22.4 +2.3/-2.2 at 71 mo				
<b>7277</b>	InSync Marquis	Mar-03	7,000	100	588	11 + 63 = 74	Total = 74	100.0 +0.0/-0.1	98.9 +0.3/-0.4	98.2 +0.2/-0.2	97.7 +0.5/-0.6	97.6 +0.5/-0.7					
								99.8 +0.1/-0.1	91.9 +0.8/-0.9	86.4 +1.1/-1.2	78.7 +1.6/-1.7	18.0 +3.0/-2.8					
						(10) + (0) = (10) (advisory-related subset)	(10)	99.9 +0.0/-0.0	99.7 +0.1/-0.1	98.7 +0.2/-0.2	98.0 +0.2/-0.3	97.9 +0.2/-0.3					
<b>7289</b>	InSync II Marquis	Jul-03	28,000	1,000	5,001	31 + 253 = 284	Total = 284	99.7 +0.1/-0.1	96.9 +0.2/-0.3	82.9 +0.6/-0.6	21.0 +1.0/-1.0	6.6 +0.8/-0.7 at 49 mo					
						(8) + (0) = (8) (advisory-related subset)	(8)	100.0 +0.0/-0.1	99.9 +0.1/-0.1	99.6 +0.1/-0.2	99.5 +0.2/-0.2 at 47 mo						
<b>7297</b>	InSync Sentry	Nov-04	9,000	4,000	458	1 + 29 = 30	Total = 30	99.8 +0.1/-0.1	97.8 +0.3/-0.4	89.8 +0.7/-0.8	73.6 +3.1/-3.4 at 47 mo						
								100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1 at 42 mo						
<b>7299</b>	InSync Sentry	Apr-05	31,000	18,000	708	5 + 43 = 48	Total = 48	99.8 +0.0/-0.1	97.9 +0.2/-0.2	90.5 +0.5/-0.6	82.9 +2.1/-2.3 at 42 mo						
								100.0 +0.0/-0.0	99.8 +0.1/-0.1	99.6 +0.1/-0.1	99.4 +0.1/-0.2 at 52 mo						
<b>7303</b>	InSync Maximo	Jun-04	17,000	6,000	1,408	5 + 57 = 62	Total = 62	99.9 +0.1/-0.1	97.7 +0.2/-0.3	89.0 +0.6/-0.6	65.1 +1.2/-1.3	42.9 +3.8/-3.9 at 52 mo					
								100.0 +0.0/-0.0	99.9 +0.2/-0.3	99.7 +0.1/-0.1	99.7 +0.1/-0.1 at 42 mo						
<b>7304</b>	InSync Maximo	Apr-05	18,000	11,000	344	2 + 22 = 24	Total = 24	99.8 +0.1/-0.1	98.1 +0.2/-0.3	92.2 +0.6/-0.7	86.6 +1.3/-1.5 at 42 mo						
								100.0 +0.0/-0.0	99.9 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1 at 42 mo						
								99.8 +0.1/-0.1	98.1 +0.2/-0.3	92.2 +0.6/-0.7	86.6 +1.3/-1.5 at 42 mo						

continued

**Device Survival Summary** continued

Model Number	Family	US Market Release	Registered US Implants	Estimated Active US Implants	Normal Battery Depletions (US)	Malfunctions (US)		Device Survival Probability (%)									
						Therapy Function Compromised	Therapy Function Not Compromised	Total	1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr
8040	InSync	Aug-01	15,000	3,000	533	21	7	28	100.0 +0.0/-0.0	100.0 +0.0/-0.1	99.9 +0.0/-0.1	99.8 +0.1/-0.1	99.7 +0.1/-0.2	99.6 +0.1/-0.2	99.6 +0.1/-0.2	99.6 +0.1/-0.2 at 86 mo	99.6 +0.1/-0.2 at 86 mo
						21	7	28	99.9 +0.0/-0.1	99.7 +0.1/-0.1	98.3 +0.2/-0.3	96.2 +0.4/-0.4	91.3 +0.7/-0.7	84.6 +0.9/-1.0	73.5 +1.9/-2.0	69.5 +2.7/-2.9 at 86 mo	69.5 +2.7/-2.9 at 86 mo
8042	InSync III	Feb-03	29,000	16,000	243	3	3	6	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1 at 73 mo		
						3	3	6	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.2 +0.1/-0.2	97.2 +0.3/-0.4	93.9 +0.6/-0.7	82.1 +4.1/-5.1	79.0 +4.8/-5.9 at 73 mo		
C154DWK, C164AWK, C174AWK	Concerto	May-06	68,000	57,000	168	14	36	50	99.9 +0.0/-0.0	99.8 +0.1/-0.1	99.4 +0.3/-0.5 at 29 mo						
						14	36	50	99.8 +0.0/-0.1	97.5 +0.3/-0.3	93.8 +1.3/-1.7 at 29 mo						
D224TRK	Consulta CRT-D	Aug-08	2,000	2,000	0	0	0	0	100.0 +0.0/-0.0 at 3 mo								
						0	0	0	100.00 +0.0/-0.0 at 3 mo								
D284TRK	Maximo II CRT-D	Mar-08	1,000	1,000	0	0	0	0	100.0 +0.0/-0.0 at 2 mo								
						0	0	0	100.00 +0.0/-0.0 at 2 mo								

## Reference Chart

The longevity estimates provided are mean values calculated for the parameters given. The actual longevity achieved for any device while implanted will depend on the actual programmed parameters and patient factors, and may differ significantly from these estimates.

Model Number	Family	Connector Style	Volume/ Mass*	Delivered Energy	Estimated Longevity					Elective Replacement (ERI)***		End of Life (EOL) Battery Voltage
					Charging Frequency**	100% Pacing#	50% Pacing#	15% Pacing#	100% Sensing	Battery Voltage	Charge Time	
7272	InSync ICD	DR+LV	66 cc 117 g	34 J	Monthly	5.4	6.3	7.3	7.8	≤ 4.91 V	—	≤ 4.57 V
					Quarterly	6.5	8.0	9.4	10.3			
					Semiannual	6.9	8.5	10.3	11.2			
7277	InSync Marquis	DR+LV split	38 cc 77 g	30 J	Monthly	3.7	4.3	4.7	4.9	≤ 2.62 V	> 16 second charge time	3 months after ERI
					Quarterly	5.0	6.0	7.0	7.5			
					Semiannual	5.5	6.7	8.0	8.6			
7289	InSync II Marquis	DR+LV true	38 cc 76 g	30 J	Monthly	3.3	3.6	4.0	4.2	≤ 2.62 V	> 16 second charge time	3 months after ERI
					Quarterly	4.2	4.9	5.5	5.8			
					Semiannual	4.5	5.4	6.1	6.6			
7295	InSync II Protect	DR+LV true	38 cc 77 g	30 J	Monthly	3.3	3.7	4.0	4.2	≤ 2.62 V	> 16 second charge time	3 months after ERI
					Quarterly	4.2	4.9	5.5	5.9			
					Semiannual	4.5	5.4	6.2	6.6			
7297	InSync Sentry	DR+LV true	40 cc 78 g	35 J	Monthly	3.3	3.7	4.1	4.3	≤ 2.62 V	> 16 second charge time	3 months after ERI
					Quarterly	4.5	5.3	6.2	6.6			
					Semiannual	5.0	6.0	7.1	7.7			
7299	InSync Sentry	DR+LV true	40 cc 78 g	35 J	Monthly	3.3	3.7	4.1	4.3	≤ 2.62 V	> 16 second charge time	3 months after ERI
					Quarterly	4.5	5.3	6.2	6.6			
					Semiannual	5.0	6.0	7.1	7.7			
7303	InSync Maximo	DR+LV true	40 cc 78 g	35 J	Monthly	3.3	3.7	4.1	4.3	≤ 2.62 V	> 16 second charge time	3 months after ERI
					Quarterly	4.5	5.3	6.2	6.6			
					Semiannual	5.0	6.0	7.1	7.7			
7304	InSync Maximo	DR+LV true	40 cc 78 g	35 J	Monthly	3.3	3.7	4.1	4.3	≤ 2.62 V	> 16 second charge time	3 months after ERI
					Quarterly	4.5	5.3	6.2	6.6			
					Semiannual	5.0	6.0	7.1	7.7			

Family	Model Number	Amplitude Setting	Estimated Longevity		Elective Replacement Time Indicators
			500 Lead Ω	1000 Lead Ω	
InSync	8040	Low 2.5 V (A, RV)	11.9	13.7	**
		Nominal 3.5 V (A, RV)	8.9	11.4	
		High 5.0 V (A, RV)	6.6	9.1	
InSync III	8042	Low 2.5 V (A, RV, LV)	8.3	9.9	**
		Nominal 3.5 V (A, RV, LV)	5.9	7.8	
		High 5.0 V (A, RV, LV)	4.1	6.0	

Model Number	Family	Connector Style	Volume/ Mass*	Delivered Energy	Estimated Longevity					Recommended Replacement (RRT)***		End of Service (EOS)
					Charging Frequency**	100% Pacing#	50% Pacing#	15% Pacing#	100% Sensing	Battery Voltage	Charge Time	
C154DWK, C164AWK, C174AWK	Concerto	DR+LV true	38 cc 68 g	35 J	Monthly	3.8	4.3	4.8	5.0	≤ 2.62 V	—	3 month after RRT or > 16-second charge time
					Quarterly	5.5	6.8	8.0	8.8			
					Semiannual	6.3	8.0	9.8	11.0			
D224TRK	Consulta	DR+LV true	38 cc/ 68 g	35 J	Monthly	3.2	3.8	4.4	4.7	≤ 2.63 V	—	3 month after RRT or > 16-second charge time
					Quarterly	4.4	5.5	6.8	7.5			
					Semiannual	4.8	6.2	7.9	9.0			
D284TRK	Maximo II	DR+LV true	38 cc/ 68 g	35 J	Monthly	3.2	3.8	4.4	4.7	≤ 2.63 V	—	3 month after RRT or > 16-second charge time
					Quarterly	4.4	5.5	6.8	7.5			
					Semiannual	4.8	6.2	7.9	9.0			

\* Volume and mass differ by connector style.

\*\* A full charge is a full energy therapeutic shock or capacitor reformation.

\*\*\* The minimum time between ERI and EOL (or RRT and EOS) is 3 months (100% pacing, two charges per month).

‡ Pacing mode is DDD for CRT models. Parameter settings; lower rate at 60 ppm, sensing rate at 70 bpm, (A, RV, LV) 3.0 V amplitude, 0.4 ms pulse width, and 510-ohm pace load per applicable channel. CRT models with shared biventricular pacing; InSync Marquis 7277 (LV impedance set to 510 ohms), InSync ICD 7272 (RV amplitude set to 4.0 V).

**7227 GEM**

Product Characteristics

US Market Release	Oct-98
Registered US Implants	22,000
Estimated Active US Implants	3,000
Normal Battery Depletions (US)	1,824


Malfunctions (US) 148

NBD Code	VVEV
Serial Number Prefix	PIP, PLN, PLP, PLR

Max Delivered Energy 35 J

Estimated Longevity [See page 37](#)

**Advisories:** [See page 155](#) – 1999 Potential Circuit Overload



	1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	9 yr	at 112 mo
% ———	99.7	99.6	99.5	99.4	99.2	99.2	99.1	99.1	99.0	99.0
% - - - -	99.3	98.9	98.7	98.1	97.2	94.8	82.7	64.6	33.4	13.1
#	20,000	17,000	15,000	13,000	11,000	8,000	5,000	3,000	1,000	200
Effective Sample Size										

**7229 GEM II VR**

Product Characteristics

US Market Release	Jul-99
Registered US Implants	11,000
Estimated Active US Implants	40
Normal Battery Depletions (US)	1,932

Malfunctions (US) 27


NBD Code	VVEV
Serial Number Prefix	PJJ

Max Delivered Energy 30 J

Estimated Longevity [See page 37](#)

**Advisories:** [See page 155](#) – 1999 Potential Circuit Overload

[Also see page 165](#) – Performance note on ICD Battery Discharge Behavior



	1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	at 89 mo
% ———	99.9	99.8	99.8	99.7	99.7	99.7	99.7	99.7
% - - - -	99.8	99.4	99.2	98.4	94.1	82.8	43.2	8.3
#	9,000	8,000	7,000	6,000	5,000	4,000	2,000	300
Effective Sample Size								






**7230 Marquis VR**

Product Characteristics


US Market Release	Dec-02	Malfunctions (US)	31	NBD Code	VVEV
Registered US Implants	19,000	<b>Therapy Function Not Compromised</b>	21	Serial Number Prefix	PKD, PLW, PLY
Estimated Active US Implants	9,000	Electrical Component	11	Max Delivered Energy	30 J
Normal Battery Depletions (US)	52	Software/Firmware	1	Estimated Longevity	<a href="#">See page 37</a>
<b>Advisories:</b> <a href="#">See page 153</a> – 2005 Potential Premature Battery Depletion Due to Battery Short		Possible Early Battery Depletion	8		
		Other	1		
		<b>Therapy Function Compromised</b>	10		
		Battery (3 malfunction related to advisory)	5		
		Electrical Component	5		



**7231 GEM III VR**

Product Characteristics


US Market Release	Dec-00	Malfunctions (US)	37	NBD Code	VVEV
Registered US Implants	17,000	<b>Therapy Function Not Compromised</b>	27	Serial Number Prefix	PJL
Estimated Active US Implants	7,000	Battery	1	Max Delivered Energy	30 J
Normal Battery Depletions (US)	795	Electrical Component	22	Estimated Longevity	<a href="#">See page 37</a>
<b>Performance Note:</b> <a href="#">See page 165</a> – Performance note on ICD Battery Discharge Behavior		Possible Early Battery Depletion	4		
		<b>Therapy Function Compromised</b>	10		
		Battery	1		
		Electrical Component	9		



**7232 Maximo VR**

Product Characteristics


US Market Release	Oct-03	Malfunctions (US)	35	NBD Code	VVED
Registered US Implants	43,000	<b>Therapy Function Not Compromised</b>	25	Serial Number Prefix	PRN
Estimated Active US Implants	30,000	Electrical Component	13	Max Delivered Energy	35 J
Normal Battery Depletions (US)	47	Possible Early Battery Depletion	12	Estimated Longevity	<a href="#">See page 37</a>
<b>Advisories:</b> <a href="#">See page 153</a> – 2005 Potential Premature Battery Depletion Due to Battery Short		<b>Therapy Function Compromised</b>	10		
		Electrical Component	8		
		Electrical Interconnect	1		
		Possible Early Battery Depletion	1		



**7271 GEM DR**

Product Characteristics


US Market Release	Oct-98	Malfunctions (US)	96	NBD Code	VVED
Registered US Implants	15,000			Serial Number Prefix	PIM
Estimated Active US Implants	2,000			Max Delivered Energy	27 J
Normal Battery Depletions (US)	1,370			Estimated Longevity	<a href="#">See page 37</a>
Advisories	None				



**7274 Marquis DR**

Product Characteristics


US Market Release	Mar-02	Malfunctions (US)	150	NBD Code	VVED
Registered US Implants	48,000	<b>Therapy Function Not Compromised</b>	75	Serial Number Prefix	PKC
Estimated Active US Implants	15,000	Battery (3 malfunctions related to advisory)	5	Max Delivered Energy	30 J
Normal Battery Depletions (US)	1,495	Electrical Component	23	Estimated Longevity	<a href="#">See page 37</a>
<b>Advisories:</b> <a href="#">See page 153</a> – 2005 Potential Premature Battery Depletion Due to Battery Short		Possible Early Battery Depletion	47		
		<b>Therapy Function Compromised</b>	75		
		Battery (42 malfunctions related to advisory)	52		
		Electrical Component	23		



**7275 GEM III DR**

Product Characteristics


US Market Release	Nov-00	Malfunctions (US)	41	NBD Code	VVED
Registered US Implants	20,000	<b>Therapy Function Not Compromised</b>	30	Serial Number Prefix	PJM
Estimated Active US Implants	2,000	Battery	1	Max Delivered Energy	30 J
Normal Battery Depletions (US)	3,471	Electrical Component	11	Estimated Longevity	<a href="#">See page 37</a>
<b>Performance Note:</b> <a href="#">See page 165</a> – Performance note on ICD Battery Discharge Behavior		Software/Firmware	1		
		Possible Early Battery Depletion	17		
		<b>Therapy Function Compromised</b>	11		
		Battery	2		
		Electrical Component	8		
		Electrical Interconnect	1		



**7276 GEM III AT**

Product Characteristics


US Market Release	Feb-01	Malfunctions (US)	46	NBD Code	DDED
Registered US Implants	14,000	<b>Therapy Function Not Compromised</b>	39	Serial Number Prefix	PKE
Estimated Active US Implants	1,000	Electrical Component	7	Max Delivered Energy	30 J
Normal Battery Depletions (US)	2,700	Software/Firmware	1	Estimated Longevity	<a href="#">See page 37</a>
<b>Performance Note:</b> <a href="#">See page 165</a> –		Possible Early Battery Depletion	31		
Performance note on ICD Battery Discharge Behavior		<b>Therapy Function Compromised</b>	7		
		Electrical Component	7		



**7278 Maximo DR**

Product Characteristics


US Market Release	Oct-03	Malfunctions (US)	29	NBD Code	VVED
Registered US Implants	37,000	<b>Therapy Function Not Compromised</b>	21	Serial Number Prefix	PRM
Estimated Active US Implants	25,000	Electrical Component	12	Max Delivered Energy	35 J
Normal Battery Depletions (US)	194	Possible Early Battery Depletion	9	Estimated Longevity	<a href="#">See page 37</a>
<b>Advisories:</b> <a href="#">See page 153</a> – 2005 Potential Premature Battery Depletion Due to Battery Short		<b>Therapy Function Compromised</b>	8		
		Electrical Component	7		
		Possible Early Battery Depletion	1		



**7288 Intrinsic**

Product Characteristics


US Market Release	Aug-04	Malfunctions (US)	33	NBD Code	VVED
Registered US Implants	31,000	<b>Therapy Function Not Compromised</b>	27	Serial Number Prefix	PUB
Estimated Active US Implants	22,000	Battery	2	Max Delivered Energy	35 J
Normal Battery Depletions (US)	85	Electrical Component	10	Estimated Longevity	<a href="#">See page 37</a>
Advisories	None	Software/Firmware	1		
		Possible Early Battery Depletion	14		
		<b>Therapy Function Compromised</b>	6		
		Electrical Component	6		



**7290 Onyx**

Product Characteristics


US Market Release	Mar-04	Malfunctions (US)	4	NBD Code	VVEV
Registered US Implants	1,000	<b>Therapy Function Not Compromised</b>	3	Serial Number Prefix	PRP
Estimated Active US Implants	1,000	Electrical Component	2	Max Delivered Energy	30 J
Normal Battery Depletions (US)	6	Possible Early Battery Depletion	1	Estimated Longevity	<a href="#">See page 37</a>
Advisories	None	<b>Therapy Function Compromised</b>	1		
		Electrical Component	1		



**D153ATG, D153DRG EnTrust**

Product Characteristics


US Market Release	Jun-05	Malfunctions (US)	3	NBD Code	DDED, VVED
Registered US Implants	400	<b>Therapy Function Not Compromised</b>	2	Serial Number Prefix	PNR
Estimated Active US Implants	300	Possible Early Battery Depletion	2	Max Delivered Energy	30 J
Normal Battery Depletions (US)	7	<b>Therapy Function Compromised</b>	1	Estimated Longevity	<a href="#">See page 38</a>
Advisories	None	Electrical Component	1		



**D154ATG, D154DRG EnTrust**

Product Characteristics

US Market Release	Jun-05	Malfunctions (US)	20	NBD Code	DDED, VVED
Registered US Implants	28,000	<b>Therapy Function Not Compromised</b>	14	Serial Number Prefix	PNR
Estimated Active US Implants	22,000	Electrical Component	5	Max Delivered Energy	35 J
Normal Battery Depletions (US)	33	Possible Early Battery Depletion	9	Estimated Longevity	<a href="#">See page 38</a>
Advisories	None	<b>Therapy Function Compromised</b>	6		
		Electrical Component	6		




**D154AWG, D164AWG Virtuoso**

Product Characteristics

US Market Release	May-06
Registered US Implants	62,000
Estimated Active US Implants	56,000
Normal Battery Depletions (US)	14
Advisories	None

Malfunctions (US)	25
<b>Therapy Function Not Compromised</b>	10
Electrical Component	8
Electrical Interconnect	1
Possible Early Battery Depletion	1
<b>Therapy Function Compromised</b>	15
Electrical Component	15

NBD Code	VVED
Serial Number Prefix	PVV, PUL
Max Delivered Energy	35 J
Estimated Longevity	<a href="#">See page 38</a>




**D154VRC EnTrust**

Product Characteristics

US Market Release	Jun-05
Registered US Implants	14,000
Estimated Active US Implants	11,000
Normal Battery Depletions (US)	7
Advisories	None

Malfunctions (US)	15
<b>Therapy Function Not Compromised</b>	11
Electrical Component	5
Possible Early Battery Depletion	6
<b>Therapy Function Compromised</b>	4
Electrical Component	4


NBD Code	VVEV
Serial Number Prefix	PNT
Max Delivered Energy	35 J
Estimated Longevity	<a href="#">See page 38</a>



**D154VWC, D164VWC Virtuoso**

Product Characteristics


US Market Release	May-06	Malfunctions (US)	8	NBD Code	VVEV
Registered US Implants	26,000	<b>Therapy Function Not Compromised</b>	3	Serial Number Prefix	PUN
Estimated Active US Implants	24,000	Electrical Component	2	Max Delivered Energy	35 J
Normal Battery Depletions (US)	7	Electrical Interconnect	1	Estimated Longevity	<a href="#">See page 38</a>
Advisories	None	<b>Therapy Function Compromised</b>	5		
		Electrical Component	5		



**D224DRG Secura DR**

Product Characteristics

US Market Release	Aug-08	Malfunctions (US)	0	NBD Code	DDED
Registered US Implants	1,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PUG
Estimated Active US Implants	1,000	<b>Therapy Function Compromised</b>	0	Max Delivered Energy	35 J
Normal Battery Depletions (US)	0			Estimated Longevity	<a href="#">See page 38</a>
Advisories	None				






**D224VRC Secura VR**

Product Characteristics

US Market Release	Aug-08
Registered US Implants	1,000
Estimated Active US Implants	1,000
Normal Battery Depletions (US)	0
Advisories	None

Malfunctions (US)	0
<b>Therapy Function Not Compromised</b>	0
<b>Therapy Function Compromised</b>	0

NBD Code	VVEV
Serial Number Prefix	PUX
Max Delivered Energy	35 J
Estimated Longevity	<a href="#">See page 38</a>




**D284DRG Maximo II DR**

Product Characteristics

US Market Release	Mar-08
Registered US Implants	1,000
Estimated Active US Implants	1,000
Normal Battery Depletions (US)	0
Advisories	None

Malfunctions (US)	0
<b>Therapy Function Not Compromised</b>	0
<b>Therapy Function Compromised</b>	0


NBD Code	VVED
Serial Number Prefix	PZM
Max Delivered Energy	35 J
Estimated Longevity	<a href="#">See page 38</a>



**D284VRC Maximo II VR**

Product Characteristics

US Market Release	Mar-08	Malfunctions (US)	0	NBD Code	VVEV
Registered US Implants	400	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PZN
Estimated Active US Implants	400	<b>Therapy Function Compromised</b>	0	Max Delivered Energy	35 J
Normal Battery Depletions (US)	0			Estimated Longevity	<a href="#">See page 38</a>
Advisories	None				



**Device Survival Summary** (95% Confidence Interval)  
 The following table shows ICD device survival estimates with 95% confidence intervals. Estimates are shown both with and without normal battery depletions included.

Model Number	Family	US Market Release	Registered US Implants	Estimated Active US Implants	Normal Battery Depletions (US)	Malfunctions (US)		Device Survival Probability (%)										
						Therapy Function Compromised	Therapy Function Not Compromised	Years After Implant										
						1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr				
7227	GEM	Oct-98	22,000	3,000	1,824	—	—	148	99.7 +0.1/-0.1	99.6 +0.1/-0.1	99.5 +0.1/-0.1	99.4 +0.1/-0.1	99.2 +0.1/-0.1	99.2 +0.1/-0.1	99.2 +0.1/-0.1	99.1 +0.1/-0.2	99.1 +0.2/-0.2	99.0 +0.2/-0.3 at 112 mo
	Advisories: See page 155 – 1999 Potential Circuit Overload						—	—	—	99.3 +0.1/-0.1	98.9 +0.1/-0.2	98.7 +0.2/-0.2	98.1 +0.2/-0.2	97.2 +0.2/-0.3	94.8 +0.4/-0.4	82.7 +0.8/-0.8	64.6 +1.2/-1.2	13.1 +1.9/-1.8 at 112 mo
7229	GEM IIVR	Jul-99	11,000	40	1,932	—	—	27	99.9 +0.0/-0.1	99.8 +0.1/-0.1	99.8 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.2	99.7 +0.1/-0.2	99.7 +0.1/-0.2 at 89 mo	
	Advisories: See page 155 – 1999 Potential Circuit Overload; Also see page 165 – Performance note on ICD Battery Discharge Behavior						—	—	—	99.8 +0.1/-0.1	99.4 +0.1/-0.2	99.2 +0.2/-0.2	98.4 +0.3/-0.3	94.1 +0.6/-0.6	82.8 +1.0/-1.1	43.2 +1.6/-1.6	8.3 +1.4/-1.3 at 89 mo	
7230	Marquis VR	Dec-02	19,000	9,000	52	10 + 21 = 31	10 + 21 = 31	31	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.1/-0.1	99.7 +0.1/-0.1	99.6 +0.1/-0.2 at 71 mo				
	Advisories: See page 153 – 2005 Potential Premature Battery Depletion Due to Battery Short						(1) (0) (1) (advisory-related subset)	(1) (0) (1) (advisory-related subset)	(1) (0) (1) (advisory-related subset)	99.8 +0.0/-0.1	99.7 +0.1/-0.1	99.5 +0.1/-0.1	99.2 +0.1/-0.2	98.7 +0.2/-0.3	97.3 +0.5/-0.6 at 71 mo			
7231	GEM III VR	Dec-00	17,000	7,000	795	10 + 27 = 37	10 + 27 = 37	37	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.8 +0.1/-0.1	99.8 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.2	99.7 +0.1/-0.2	99.7 +0.1/-0.2 at 91 mo	
	See page 165 – Performance note on ICD Battery Discharge Behavior									99.8 +0.1/-0.1	99.4 +0.1/-0.1	99.0 +0.1/-0.2	98.5 +0.2/-0.2	97.0 +0.3/-0.4	87.4 +0.7/-0.8	71.4 +1.4/-1.5	55.8 +2.9/-3.1 at 91 mo	
7232	Maximo VR	Oct-03	43,000	30,000	47	10 + 25 = 35	10 + 25 = 35	35	100.0 +0.0/-0.0	99.9 +0.0/-0.0	99.9 +0.0/-0.0	99.9 +0.0/-0.0	99.8 +0.1/-0.1					
	Advisories: See page 153 – 2005 Potential Premature Battery Depletion Due to Battery Short						(0) + (0) = (0) (advisory-related subset)	(0) + (0) = (0) (advisory-related subset)	(0) + (0) = (0) (advisory-related subset)	99.9 +0.0/-0.0	99.8 +0.0/-0.1	99.5 +0.1/-0.1	99.3 +0.1/-0.1	99.1 +0.2/-0.2				
7271	GEM DR	Oct-98	15,000	2,000	1,370	—	—	96	99.9 +0.0/-0.1	99.8 +0.1/-0.1	99.7 +0.1/-0.1	99.5 +0.1/-0.2	99.3 +0.1/-0.2	99.1 +0.2/-0.2	99.0 +0.2/-0.2	98.8 +0.2/-0.3	98.8 +0.2/-0.3 at 105 mo	
	Advisories: See page 153 – 2005 Potential Premature Battery Depletion Due to Battery Short									99.7 +0.1/-0.1	99.4 +0.1/-0.1	98.8 +0.2/-0.2	97.8 +0.3/-0.3	95.5 +0.4/-0.4	89.3 +0.7/-0.7	75.7 +1.2/-1.2	44.5 +1.8/-1.8 at 105 mo	
7274	Marquis DR	Mar-02	48,000	15,000	1,495	75 + 75 = 150	75 + 75 = 150	150	100.0 +0.0/-0.0	99.9 +0.0/-0.0	99.8 +0.0/-0.1	99.6 +0.1/-0.1	99.5 +0.1/-0.1	99.3 +0.1/-0.2	99.3 +0.1/-0.2	99.3 +0.1/-0.2 at 79 mo		
	Advisories: See page 153 – 2005 Potential Premature Battery Depletion Due to Battery Short						(36) + (3) = (39) (advisory-related subset)	(36) + (3) = (39) (advisory-related subset)	(39) (advisory-related subset)	99.8 +0.0/-0.0	99.5 +0.1/-0.1	98.6 +0.1/-0.1	97.4 +0.2/-0.2	91.9 +0.4/-0.4	68.4 +1.2/-1.2	42.9 +3.0/-3.1 at 79 mo		

continued

Device Survival Summary continued

Model Number	Family	US Market Release	Registered US Implants	Estimated Active US Implants	Normal Battery Depletions	Malfunctions		Device Survival Probability (%)												
						Therapy Function Compromised	Therapy Function Not Compromised	Years After Implant												
								1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr				
7275	GEM III DR	Nov-00	20,000	2,000	3,471	11	30	41	99.9	99.9	99.8	99.7	99.7	99.7	99.7	99.7	99.7	99.7	99.7	99.7
	See page 165 – Performance note on ICD Battery Discharge Behavior									99.6	99.0	96.9	90.1	64.5	8.8					
7276	GEM III AT	Feb-01	14,000	1,000	2,700	7	39	46	99.9	99.9	99.7	99.5	99.5	99.5	99.5	99.5	99.5	99.5	99.5	99.5
	See page 165 – Performance note on ICD Battery Discharge Behavior									99.6	98.8	96.1	84.6	50.5	3.7					
7278	Maximo DR	Oct-03	37,000	25,000	194	8	21	29	100.0	100.0	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9
	Advisories: See page 153 – 2005 Potential Premature Battery Depletion Due to Battery Short					(0)	(0)	(0)	99.9	99.7	99.3	97.7	93.8	11.1						
7288	Intrinsic	Aug-04	31,000	22,000	85	6	27	33	100.0	99.9	99.9	99.8	99.8	99.8	99.8	99.8	99.8	99.8	99.8	99.8
	Onyx	Mar-04	1,000	1,000	6	1	3	4	99.9	99.7	99.4	98.2	97.8	97.8	97.8	97.8	97.8	97.8	97.8	97.8
D153ATG, D153DRG	EnTrust DR	Jun-05	400	300	7	1	2	3	99.8	99.8	99.1	99.1	99.1	99.1	99.1	99.1	99.1	99.1	99.1	99.1
	See page 165 – Performance note on ICD Battery Discharge Behavior									99.5	99.5	97.3	91.1	3.4						
D154ATG, D154DRG	EnTrust DR	Jun-05	28,000	22,000	33	6	14	20	100.0	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9
	See page 165 – Performance note on ICD Battery Discharge Behavior									99.9	99.7	99.3	99.2	99.2	99.2	99.2	99.2	99.2	99.2	99.2

continued

Device Survival Summary continued

Model Number	Family	US Market Release	Registered US Implants	Estimated Active US Implants	Normal Battery Depletions	Malfunctions			Device Survival Probability (%)								
						Therapy Function Compromised	Therapy Function Not Compromised	Total	Years After Implant								
									1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	
D154AWG D164AWG	Virtuoso DR	May-06	62,000	56,000	14	15 + 10 = 25	Excluding Normal Battery Depletion	100.0 +0.0/-0.0	99.9 +0.0/-0.0	99.7 +0.2/-0.5 at 29 mo							
							Including Normal Battery Depletion	99.9 +0.0/-0.0	99.8 +0.0/-0.1	99.4 +0.3/-0.5 at 29 mo							
D154VRC	EnTrust VR	Jun-05	14,000	11,000	7	4 + 11 = 15	Excluding Normal Battery Depletion	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.1/-0.1 at 39 mo							
							Including Normal Battery Depletion	99.9 +0.0/-0.1	99.7 +0.1/-0.1	99.6 +0.1/-0.2 at 39 mo							
D154VWC, D164VWC	Virtuoso VR	May-06	26,000	24,000	7	5 + 3 = 8	Excluding Normal Battery Depletion	100.0 +0.0/-0.0	99.9 +0.0/-0.1								
							Including Normal Battery Depletion	99.9 +0.0/-0.1	99.7 +0.1/-0.2								
D224DRG	Secura DR	Aug-08	1,000	1,000	0	0 + 0 = 0	Excluding Normal Battery Depletion	100.0 +0.0/-0.0 at 2 mo									
							Including Normal Battery Depletion	100.0 +0.0/-0.0 at 2 mo									
D224VRC	Secura VR	Aug-08	1,000	1,000	0	0 + 0 = 0	Excluding Normal Battery Depletion	100.0 +0.0/-0.0 at 2 mo									
							Including Normal Battery Depletion	100.0 +0.0/-0.0 at 2 mo									
D284DRG	Maximo II DR	Mar-08	1,000	1,000	0	0 + 0 = 0	Excluding Normal Battery Depletion	100.0 +0.0/-0.0 at 2 mo									
							Including Normal Battery Depletion	100.0 +0.0/-0.0 at 2 mo									
D284VRC	Maximo II VR	Mar-08	400	400	0	0 + 0 = 0	Excluding Normal Battery Depletion	100.0 +0.0/-0.0 at 2 mo									
							Including Normal Battery Depletion	100.0 +0.0/-0.0 at 2 mo									

### Reference Chart

The longevity estimates provided are mean values calculated for the parameters given. The actual longevity achieved for any device while implanted will depend on the actual programmed parameters and patient factors, and may differ significantly from these estimates.

Model Number	Family	Connector Style	Volume/ Mass*	Delivered Energy	Estimated Longevity					Elective Replacement (ERI)***		End of Life (EOL) Battery Voltage
					Charging Frequency**	100% Pacing#	50% Pacing#	15% Pacing#	100% Sensing	Battery Voltage	Charge Time	
7227	GEM	B, Cx, D, E	49 cc* 90 g	35 J	Monthly	5.3	5.7	6.0	6.1	≤ 2.55 V	—	≤ 2.40 V <sup>§</sup>
					Quarterly	7.7	8.5	9.3	9.6			
					Semiannual	8.8	10.0	11.0	11.5			
7229	GEM II VR	Cx	39 cc 77 g	30 J	Monthly	3.6	3.9	4.1	4.2	≤ 2.55 V	—	≤ 2.40 V
					Quarterly	5.0	5.5	6.0	6.2			
					Semiannual	5.6	6.3	6.9	7.1			
7230	Marquis VR	B, Cx, E	36 cc 75 g	30 J	Monthly	4.9	5.2	5.4	5.5	≤ 2.62 V	> 16-second charge time	3 months after ERI
					Quarterly	7.3	8.0	8.5	8.7			
					Semiannual	8.5	9.3	10.0	10.4			
7231	GEM III VR	Cx	39 cc 77 g	30 J	Monthly	4.3	4.7	5.0	5.2	≤ 2.55 V	—	≤ 2.40 V
					Quarterly	6.0	6.8	7.4	7.8			
					Semiannual	6.6	7.5	8.5	8.9			
7232	Maximo VR	B, Cx, E	39 cc 76 g	35 J	Monthly	4.4	4.7	4.8	4.9	≤ 2.62 V	> 16-second charge time	3 months after ERI
					Quarterly	7.0	7.5	8.0	8.3			
					Semiannual	8.2	9.0	9.7	10.0			
7271	GEM DR	DR	62 cc 115 g	35 J	Monthly	6.0	6.9	7.5	7.8	≤ 4.91 V	—	≤ 4.57 V <sup>§</sup>
					Quarterly	7.4	8.4	9.3	9.8			
					Semiannual	7.9	9.0	10.0	10.6			
7274	Marquis DR	DR+LV	36 cc 75 g	30 J	Monthly	4.0	4.4	4.8	4.9	≤ 2.62 V	> 16-second charge time	3 months after ERI
					Quarterly	5.6	6.4	7.1	7.5			
					Semiannual	6.2	7.2	8.1	8.6			
7275	GEM III DR	DR	39.5 cc 78 g	30 J	Monthly	3.3	3.8	4.3	4.4	≤ 2.55 V	—	≤ 2.40 V
					Quarterly	4.2	5.0	5.8	6.3			
					Semiannual	4.5	5.5	6.5	7.0			
7276	GEM III AT	DR	39 cc 77 g	30 J	Monthly	3.3	3.8	4.3	4.5	≤ 2.55 V	—	≤ 2.40 V
					Quarterly	4.3	5.1	5.9	6.3			
					Semiannual	4.5	5.5	6.5	7.0			
7278	Maximo DR	DR	39 cc 77 g	35 J	Monthly	3.7	4.1	4.3	4.5	≤ 2.62 V	> 16-second charge time	3 months after ERI
					Quarterly	5.3	6.1	6.8	7.1			
					Semiannual	6.0	7.0	8.0	8.5			
7287	Intrinsic 30	DR	36 cc 75 g	30 J	Monthly	4.0	4.3	4.7	4.8	≤ 2.62 V	> 16-second charge time	3 months after ERI
					Quarterly	5.5	6.3	7.0	7.4			
					Semiannual	6.2	7.2	8.2	8.6			
7288	Intrinsic	DR	38 cc 76 g	35 J	Monthly	3.7	4.1	4.3	4.5	≤ 2.62 V	> 16-second charge time	3 months after ERI
					Quarterly	5.4	6.1	6.8	7.1			
					Semiannual	6.1	7.0	8.0	8.5			
7290	Onyx	Cx	39 cc 77 g	30 J	Monthly	3.8	4.1	4.3	4.5	≤ 2.55 V	> 16-second charge time	≤ 2.40 V
					Quarterly	5.0	5.6	6.2	6.4			
					Semiannual	5.4	6.1	6.7	7.0			

\* Volume and mass differ by connector style.

\*\* A full charge is a full energy therapeutic shock or capacitor reformation.

\*\*\* The minimum time between ERI and EOL is 3 months (100% pacing, two charges per month).

‡ Pacing mode is VVI for single chamber models and DDD for dual chamber and CRT models. Parameter settings; lower rate at 60 ppm, sensing rate at 70 bpm, (A, RV, LV) 3.0 V amplitude, 0.4 ms pulse width, and 510-ohm pace load per applicable channel. CRT models with shared biventricular pacing; InSync Marquis 7277 (LV impedance set to 510 ohms), InSync ICD 7272 (RV amplitude set to 4.0 V).

‡‡ For Model 7223 devices, if charge time exceeds 60 seconds, the devices are at EOL. If two consecutive charge cycles exceed 60 seconds, the “charge circuit inactive” indicator is tripped and all therapies except emergency output VVI pacing are disabled.

§ For Model 7271 and 7227 devices, if charge time exceeds 30 seconds, the device is at EOL. Immediate replacement is recommended. If three consecutive charge cycles exceed 30 seconds, the “charge circuit inactive” indicator is tripped and all therapies except emergency VVI pacing are disabled.

continued

Reference Chart continued

Model Number	Family	Connector Style	Volume/ Mass*	Delivered Energy	Charging Frequency**	Estimated Longevity				Recommended Replacement (RRT)***		End of Service (EOS)
						100% Pacing†	50% Pacing†	15% Pacing†	100% Sensing	Battery Voltage	Charge Time	
D153ATG, D153DRG	EnTrust	DR	33 cc 63 g	30 J	Monthly	3.5	3.8	4.1	4.2	≤ 2.61 V	—	3 months after RRT or > 16-second charge time
					Quarterly	4.8	5.4	6.0	6.3			
					Semiannual	5.3	6.1	6.9	7.2			
D153VRC	EnTrust	Cx	32 cc 63 g	30 J	Monthly	4.4	4.7	4.9	5.0	≤ 2.61 V	—	3 months after RRT or > 16-second charge time
					Quarterly	6.8	7.4	7.9	8.1			
					Semiannual	7.9	8.7	9.5	9.8			
D154ATG, D154DRG	EnTrust	DR	35 cc 68 g	35 J	Monthly	3.8	4.2	4.4	4.6	≤ 2.61 V	—	3 months after RRT or > 16-second charge time
					Quarterly	5.5	6.1	6.8	7.0			
					Semiannual	6.1	7.0	7.9	8.3			
D154AWG, D164AWG	Virtuoso	DR	37 cc 68 g	35 J	Monthly	4.1	4.5	4.8	5.0	≤ 2.62 V	—	3 months after RRT or > 16-second charge time
					Quarterly	6.3	7.3	8.3	8.8			
					Semiannual	7.3	8.7	10.1	11.0			
D154VRC	EnTrust	Cx	35 cc 68 g	35 J	Monthly	4.8	5.0	5.2	5.3	≤ 2.61 V	—	3 months after RRT or > 16-second charge time
					Quarterly	7.5	8.3	8.8	9.0			
					Semiannual	9.0	10.0	10.7	11.0			
D154VWC, D164VWC	Virtuoso	Cx	37 cc 68 g	35 J	Monthly	4.8	5.1	5.3	5.4	≤ 2.62 V	—	3 months after RRT or > 16-second charge time
					Quarterly	8.1	9.0	9.6	10.0			
					Semiannual	10.0	11.2	12.3	12.9			
D224DRG	Secura DR	DR	37 cc 68 g	35J	Monthly	3.60	4.08	4.50	4.67	≤ 2.63 V	—	3 months after RRT or > 16-second charge time
					Quarterly	5.07	6.05	7.00	7.50			
					Semiannual	5.70	7.00	8.27	9.00			
D224VRC	Secura VR	Cx	37 cc 68 g	35J	Monthly	4.33	4.67	4.92	5.00	≤ 2.63 V	—	3 months after RRT or > 19-second charge time
					Quarterly	6.67	7.45	8.05	8.41			
					Semiannual	7.76	8.85	9.79	10.25			
D284DRG	Maximo II DR	DR	37 cc 68 g	35J	Monthly	3.60	4.08	4.50	4.67	≤ 2.63 V	—	3 months after RRT or > 16-second charge time
					Quarterly	5.07	6.05	7.00	7.50			
					Semiannual	5.70	7.00	8.27	9.00			
D284VRC	Maximo II VR	Cx	37 cc 68 g	35J	Monthly	4.33	4.67	4.92	5.00	≤ 2.63 V	—	3 months after RRT or > 19-second charge time
					Quarterly	6.67	7.45	8.05	8.41			
					Semiannual	7.76	8.85	9.79	10.25			










\* Volume and mass differ by connector style.

\*\* A full charge is a full energy therapeutic shock or capacitor reformation.

\*\*\* The minimum time between RRT and EOS is 3 months (100% pacing, two charges per month).

† Pacing mode is VVI for single chamber models and DDD for dual chamber models. Parameter settings; lower rate at 60 ppm, sensing rate at 70 bpm, (A, RV, LV) 3.0 V amplitude, 0.4 ms pulse width, and 510-ohm pace load per applicable channel.

ICD Connector Styles






**Adapta DR ADDR01, ADDR03, ADDR06, ADD01**

Product Characteristics


US Market Release	Jul-06	Malfunctions (US)	10	NBG Code	DDDR, DDD
Registered US Implants	116,000	<b>Therapy Function Not Compromised</b>	6	Serial Number Prefix	PWB, PWD, PWC, PWF, NWB, NWC, NWD
Estimated Active US Implants	104,000	Electrical Component	6		
Normal Battery Depletions (US)	1	<b>Therapy Function Compromised</b>	4		
Advisories	None	Electrical Component	4	Estimated Longevity	<a href="#">See page 77</a>



**Adapta DR ADDR1**

Product Characteristics


US Market Release	Jul-06	Malfunctions (US)	0	NBG Code	DDDR
Registered US Implants	11,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PWE, NWE
Estimated Active US Implants	10,000	<b>Therapy Function Compromised</b>	0		
Normal Battery Depletions (US)	0			Estimated Longevity	<a href="#">See page 77</a>
Advisories	None				



**Adapta DR ADDR1**

Product Characteristics


US Market Release	Jul-06	Malfunctions (US)	1	NBG Code	SSIR
Registered US Implants	10,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PWA
Estimated Active US Implants	9,000	<b>Therapy Function Compromised</b>	1		
Normal Battery Depletions (US)	1	Electrical Component	1	Estimated Longevity	<a href="#">See page 77</a>
Advisories	None				



**Adapta SR ADSR01, ADSR03, ADSR06**

Product Characteristics


US Market Release	Jul-06	Malfunctions (US)	0	NBG Code	SSIR
Registered US Implants	22,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	NWN, NWM, NWP
Estimated Active US Implants	19,000	<b>Therapy Function Compromised</b>	0		
Normal Battery Depletions (US)	1			Estimated Longevity	<a href="#">See page 77</a>
Advisories	None				



**Adapta VDD ADVDD01**

Product Characteristics


US Market Release	Jul-06	Malfunctions (US)	0	NBG Code	VDO
Registered US Implants	400	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PWG, NWG
Estimated Active US Implants	400	<b>Therapy Function Compromised</b>	0	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	0				
Performance Note: <a href="#">See page 163</a> – Performance note on AT500 Pacing System Follow-Up Protocol					



**AT500 AT501, 7253**

Product Characteristics


US Market Release	Mar-03	Malfunctions (US)	10	NBG Code	DDDRP
Registered US Implants	11,000	<b>Therapy Function Not Compromised</b>	5	Serial Number Prefix	IJF
Estimated Active US Implants	5,000	Electrical Component	2	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	481	Possible Early Battery Depletion	3		
Performance Note: <a href="#">See page 163</a> – Performance note on AT500 Pacing System Follow-Up Protocol					
		<b>Therapy Function Compromised</b>	5		
		Electrical Component	3		
		Electrical Interconnect	1		
		Possible Early Battery Depletion	1		



**EnPulse DR E1DR01, E1DR03, E1DR06**

Product Characteristics


US Market Release	Dec-03	Malfunctions (US)	1	NBG Code	DDDR
Registered US Implants	7,000	<b>Therapy Function Not Compromised</b>	1	Serial Number Prefix	PRA
Estimated Active US Implants	4,000	Electrical Component	1	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	26	<b>Therapy Function Compromised</b>	0		
Advisories	None				



**EnPulse DR E1DR21**

Product Characteristics


US Market Release	Dec-03	Malfunctions (US)	0	NBG Code	DDDR
Registered US Implants	2,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PPT
Estimated Active US Implants	1,000	<b>Therapy Function Compromised</b>	0	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	23				
Advisories	None				



**EnPulse 2 DR E2DR01, E2DR03, E2DR06**

Product Characteristics


US Market Release	Feb-04	Malfunctions (US)	13	NBG Code	DDDR
Registered US Implants	101,000	<b>Therapy Function Not Compromised</b>	10	Serial Number Prefix	PNB, PNC, PNH
Estimated Active US Implants	70,000	Electrical Component	9	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	81	Possible Early Battery Depletion	1		
Advisories	None	<b>Therapy Function Compromised</b>	3		
		Battery	1		
		Electrical Component	2		



**EnPulse 2 DR E2DR21**

Product Characteristics


US Market Release	Feb-04	Malfunctions (US)	1	NBG Code	DDDR
Registered US Implants	12,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PMU
Estimated Active US Implants	8,000	<b>Therapy Function Compromised</b>	1	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	48	Electrical Component	1		
Advisories	None				



**EnPulse 2 DR E2DR31, E2DR33**

Product Characteristics


US Market Release	Feb-04	Malfunctions (US)	0	NBG Code	DDDR
Registered US Implants	1,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PNL
Estimated Active US Implants	400	<b>Therapy Function Compromised</b>	0	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	0				
Advisories	None				



**EnPulse 2 SR E2SR01, E2SR03, E2SR06**

Product Characteristics


US Market Release	Dec-03	Malfunctions (US)	3	NBG Code	SSIR
Registered US Implants	25,000	<b>Therapy Function Not Compromised</b>	2	Serial Number Prefix	PMW, PMY, PNA
Estimated Active US Implants	15,000	Electrical Component	1	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	42	Possible Early Battery Depletion	1		
Advisories	None	<b>Therapy Function Compromised</b>	1		
		Other	1		



**EnPulse 2 VDD E2VDD01**

Product Characteristics


US Market Release	Dec-03	Malfunctions (US)	0	NBG Code	VDD
Registered US Implants	1,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PMV
Estimated Active US Implants	500	<b>Therapy Function Compromised</b>	0	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	1				
Advisories	None				



**EnRhythm DR P1501DR**

Product Characteristics


US Market Release	May-05	Malfunctions (US)	30	NBG Code	DDDRP
Registered US Implants	78,000	<b>Therapy Function Not Compromised</b>	5	Serial Number Prefix	PNP
Estimated Active US Implants	62,000	Electrical Component	5	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	4	<b>Therapy Function Compromised</b>	25		
Advisories	None	Electrical Component	23		
		Possible Early Battery Depletion	1		



**Kappa 400 DR KDR401, KDR403**

Product Characteristics

US Market Release	Jan-98	Malfunctions (US)	22	NBG Code	DDD/RO
Registered US Implants	47,000	<b>Therapy Function Not Compromised</b>	13	Serial Number Prefix	PER, PET
Estimated Active US Implants	5,000	Electrical Component	9	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	5,271	Electrical Interconnect	1		
Advisories	None	Possible Early Battery Depletion	2		
		Other	1		
		<b>Therapy Function Compromised</b>	9		
		Electrical Component	7		
		Electrical Interconnect	2		






**Kappa 400 SR** KSR401, KSR403

Product Characteristics


US Market Release	Feb-98	Malfunctions (US)	5	NBG Code	SSI/R
Registered US Implants	15,000	<b>Therapy Function Not Compromised</b>	4	Serial Number Prefix	PEU, PGD
Estimated Active US Implants	2,000	Electrical Component	3	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	870	Possible Early Battery Depletion	1		
Advisories	None	<b>Therapy Function Compromised</b>	1		
		Electrical Interconnect	1		



**Kappa 600 DR** KDR601, KDR603, KDR606

Product Characteristics


US Market Release	Jan-99	Malfunctions (US)	33	NBG Code	DDD/RO
Registered US Implants	24,000	<b>Therapy Function Not Compromised</b>	3	Serial Number Prefix	PHF, PHH, PHG
Estimated Active US Implants	2,000	Electrical Component	3	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	2,417	<b>Therapy Function Compromised</b>	30		
Advisories: <a href="#">See page 154</a> – 2002 Potential Fractured Power Supply Wires		Electrical Component	2		
		Electrical Interconnect	28		
		<i>(15 malfunctions related to advisory)</i>			



**Kappa 600 DR** KDR651, KDR653

Product Characteristics


US Market Release	Mar-01	Malfunctions (US)	19	NBG Code	DDD/RO
Registered US Implants	14,000	<b>Therapy Function Not Compromised</b>	2	Serial Number Prefix	PLJ, PLK
Estimated Active US Implants	5,000	Electrical Component	1	Estimated Longevity	<a href="#">See page 77</a>
Normal Battery Depletions (US)	569	Possible Early Battery Depletion	1		
<b>Advisories:</b> <a href="#">See page 154</a> – 2002 Potential Fractured Power Supply Wires		<b>Therapy Function Compromised</b>	17		
		Electrical Component	1		
		Electrical Interconnect	16		
		<i>(1 malfunction related to advisory)</i>			



**Kappa 700 D** KD701, KD703, KD706

Product Characteristics


US Market Release	Jan-99	Malfunctions (US)	0	NBG Code	DDD
Registered US Implants	300	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PHK
Estimated Active US Implants	100	<b>Therapy Function Compromised</b>	0	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	11				
<b>Advisories:</b> <a href="#">See page 154</a> – 2002 Potential Fractured Power Supply Wires					



**Kappa 700 DR** KDR701, KDR703, KDR706

Product Characteristics


US Market Release	Feb-99	Malfunctions (US)	351	NBG Code	DDD/RO
Registered US Implants	192,000	<b>Therapy Function Not Compromised</b>	28	Serial Number Prefix	PGU, PGY, PGW
Estimated Active US Implants	54,000	Battery	1	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	12,539	Electrical Component	21		
<b>Advisories:</b> <a href="#">See page 154</a> – 2002 Potential Fractured Power Supply Wires		Electrical Interconnect	1		
		Possible Early Battery Depletion	3		
		Other	2		
		<b>Therapy Function Compromised</b>	323		
		Electrical Component	15		
		Electrical Interconnect	308		
		<i>(136 malfunctions related to advisory)</i>			



**Kappa 700 DR KDR721**

Product Characteristics

US Market Release	Feb-99	Malfunctions (US)	5	NBG Code	DDD/RO
Registered US Implants	10,000	<b>Therapy Function Not Compromised</b>	1	Serial Number Prefix	PGR
Estimated Active US Implants	100	Electrical Component	1	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	1,237	<b>Therapy Function Compromised</b>	4		
<b>Advisories:</b> <a href="#">See page 154</a> – 2002 Potential Fractured Power Supply Wires		Electrical Interconnect	4		
		<i>(4 malfunctions related to advisory)</i>			




	1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr
% ———	100.0	100.0	100.0	100.0	99.9	99.9	99.9
% - - - - -	99.9	99.7	98.9	96.8	91.1	68.8	20.0
#	8,000	7,000	6,000	6,000	4,000	2,000	200
Effective Sample Size							



**Kappa 700 SR** KSR701, KSR703, KSR706

Product Characteristics


US Market Release	Feb-99	Malfunctions (US)	17	NBG Code	SSI/R
Registered US Implants	55,000	<b>Therapy Function Not Compromised</b>	3	Serial Number Prefix	PHT, PHW, PHU
Estimated Active US Implants	13,000	Electrical Component	2	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	2,265	Possible Early Battery Depletion	1		
Advisories	None	<b>Therapy Function Compromised</b>	14		
		Electrical Component	4		
		Electrical Interconnect	10		



**Kappa 700 VDD** KVDD701

Product Characteristics


US Market Release	Jan-99	Malfunctions (US)	3	NBG Code	VDD/RO
Registered US Implants	2,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PHP
Estimated Active US Implants	100	<b>Therapy Function Compromised</b>	3	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	149	Electrical Interconnect	3		
<b>Advisories:</b> <a href="#">See page 154</a> – 2002 Potential Fractured Power Supply Wires		<i>(3 malfunctions related to advisory)</i>			



**Kappa 800 DR** KDR801, KDR803

Product Characteristics


US Market Release	Jan-02	Malfunctions (US)	1	NBG Code	DDD/RO
Registered US Implants	4,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PKW, PKY
Estimated Active US Implants	2,000	<b>Therapy Function Compromised</b>	1	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	27	Electrical Interconnect	1		
Advisories	None				



**Kappa 900 DR** KDR901, KDR903, KDR906

Product Characteristics


US Market Release	Jan-02	Malfunctions (US)	36	NBG Code	DDD/RO
Registered US Implants	125,000	<b>Therapy Function Not Compromised</b>	12	Serial Number Prefix	PKM, PKN, PKP
Estimated Active US Implants	70,000	Electrical Component	12	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	907	<b>Therapy Function Compromised</b>	24		
Advisories	None	Electrical Component	8		
		Electrical Interconnect	16		



**Kappa 920 DR** **KDR921**

Product Characteristics


US Market Release	Jan-02	Malfunctions (US)	3	NBG Code	DDD/RO
Registered US Implants	16,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PKR
Estimated Active US Implants	6,000	<b>Therapy Function Compromised</b>	3	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	902	Electrical Interconnect	3		
Advisories	None				



**Kappa 900 SR** **KSR901, KSR903, KSR906**

Product Characteristics


US Market Release	Jan-02	Malfunctions (US)	11	NBG Code	VVEV
Registered US Implants	37,000	<b>Therapy Function Not Compromised</b>	8	Serial Number Prefix	PLF, PLG, PLH
Estimated Active US Implants	17,000	Electrical Component	7	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	240	Possible Early Battery Depletion	1		
Advisories	None	<b>Therapy Function Compromised</b>	3		
		Electrical Interconnect	3		



**Kappa 900 VDD KVDD901**

Product Characteristics


US Market Release	Jan-02	Malfunctions (US)	0	NBG Code	VDD
Registered US Implants	1,000	Therapy Function Not Compromised	0	Serial Number Prefix	PLE
Estimated Active US Implants	300	Therapy Function Compromised	0	Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	9				
Advisories	None				



**Legend 8416, 8417, 8417M, 8418, 8419**

Product Characteristics

US Market Release	Aug-89	Malfunctions (US)	143	NBG Code	SSIRO
Registered US Implants	57,000			Serial Number Prefix	XT, WJ, WN, ZT
Estimated Active US Implants	2,000			Estimated Longevity	<a href="#">See page 78</a>
Normal Battery Depletions (US)	2,900				
Advisories	None				






**Legend II 8424, 8426, 8427**

Product Characteristics

US Market Release	Nov-91
Registered US Implants	59,000
Estimated Active US Implants	4,000
Normal Battery Depletions (US)	2,203
Advisories	None

Malfunctions (US) 36

NBG Code	SSIRO
Serial Number Prefix	2P, 2T, 2U
Estimated Longevity	<a href="#">See page 78</a>



**Minix/Minix ST 8330, 8331, 8331M, 8340, 8341, 8341M, 8342**


Product Characteristics

US Market Release	Dec-89
Registered US Implants	58,000
Estimated Active US Implants	4,000
Normal Battery Depletions (US)	1,594

Malfunctions (US) 49

NBG Code	SSIRO
Serial Number Prefix	2P, 2T, 2U
Estimated Longevity	<a href="#">See page 78</a>

Advisories: [See page 159](#) – 1991 Potential Delayed Restoration of Permanent Settings




**Minuet 7107, 7108**

Product Characteristics

US Market Release	Mar-92
Registered US Implants	17,000
Estimated Active US Implants	1,000
Normal Battery Depletions (US)	762
Advisories	None

Malfunctions (US) 4

NBG Code	DDDCO
Serial Number Prefix	1Z1, 2G1
Estimated Longevity	<a href="#">See page 78</a>




**Preva DR 7088, 7089**

Product Characteristics

US Market Release	Jul-96
Registered US Implants	26,000
Estimated Active US Implants	3,000
Normal Battery Depletions (US)	1,828
Advisories	None

Malfunctions (US) 4

NBG Code	DDD/RO
Serial Number Prefix	PGJ, PGK
Estimated Longevity	<a href="#">See page 78</a>




**Preva SR 8088, 8089**

Product Characteristics

US Market Release	Jul-96
Registered US Implants	18,000
Estimated Active US Implants	2,000
Normal Battery Depletions (US)	687
Advisories	None

Malfunctions (US) 1

NBG Code	SSI/R
Serial Number Prefix	PGL, PGM
Estimated Longevity	<a href="#">See page 78</a>




**Prevail S 8085, 8086**

Product Characteristics

US Market Release	Oct-95
Registered US Implants	4,000
Estimated Active US Implants	1,000
Normal Battery Depletions (US)	38
Advisories	None

Malfunctions (US) 1

NBG Code	SSI
Serial Number Prefix	PGL, PGM
Estimated Longevity	<a href="#">See page 78</a>




**Prodigy D 7864, 7865, 7866**

Product Characteristics

US Market Release	Oct-95
Registered US Implants	3,000
Estimated Active US Implants	300
Normal Battery Depletions (US)	160
Advisories	None

Malfunctions (US) 0

NBG Code	DDDCO
Serial Number Prefix	PDL, PDM, PDN
Estimated Longevity	<a href="#">See page 79</a>




**Prodigy DR 7860, 7861, 7862**

Product Characteristics

US Market Release	Oct-95
Registered US Implants	38,000
Estimated Active US Implants	4,000
Normal Battery Depletions (US)	2,522
Advisories	None

Malfunctions (US) 11

NBG Code	DDD/RO
Serial Number Prefix	PDH, PDJ, PDK
Estimated Longevity	<a href="#">See page 79</a>



**Prodigy S 8164, 8165, 8166**


Product Characteristics

US Market Release	Oct-95
Registered US Implants	2,000
Estimated Active US Implants	300
Normal Battery Depletions (US)	41
Advisories	None

Malfunctions (US) 0

NBG Code	SSIC
Serial Number Prefix	PEG, PEH, PEJ

Estimated Longevity [See page 79](#)



**Prodigy SR 8158, 8160, 8161, 8162**


Product Characteristics

US Market Release	Oct-95
Registered US Implants	22,000
Estimated Active US Implants	3,000
Normal Battery Depletions (US)	836
Advisories	None

Malfunctions (US) 4

NBG Code	SSI/R
Serial Number Prefix	PEM, PED, PEE, PEF


Estimated Longevity [See page 79](#)



**Sensia DR SEDR01, SED01**

Product Characteristics


US Market Release	Jul-06	Malfunctions (US)	1	NBG Code	DDD, DDR
Registered US Implants	42,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PWL, PWK, NWL
Estimated Active US Implants	37,000	<b>Therapy Function Compromised</b>	1		
Normal Battery Depletions (US)	0	Electrical Component	1	Estimated Longevity	<a href="#">See page 79</a>
Advisories	None				



**Sensia SR SESR01, SES01**

Product Characteristics


US Market Release	Jul-06	Malfunctions (US)	0	NBG Code	SSIR, SSI
Registered US Implants	26,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PWR, PWS, NWR
Estimated Active US Implants	22,000	<b>Therapy Function Compromised</b>	0		
Normal Battery Depletions (US)	1			Estimated Longevity	<a href="#">See page 79</a>
Advisories	None				



**Sigma 100 S SS103, SS106**

Product Characteristics


US Market Release	Aug-99	Malfunctions (US)	0	NBG Code	SSI
Registered US Implants	1,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PJG, PJH
Estimated Active US Implants	200	<b>Therapy Function Compromised</b>	0	Estimated Longevity	<a href="#">See page 79</a>
Normal Battery Depletions (US)	5				
<b>Advisories:</b> <a href="#">See page 152</a> – 2005 Potential Separation of Interconnect Wires					



**Sigma 200 DR SDR203**

Product Characteristics

US Market Release	Aug-99	Malfunctions (US)	6	NBG Code	DDD/RO
Registered US Implants	16,000	<b>Therapy Function Not Compromised</b>	1	Serial Number Prefix	PJD
Estimated Active US Implants	6,000	Electrical Component	1	Estimated Longevity	<a href="#">See page 79</a>
Normal Battery Depletions (US)	117	<b>Therapy Function Compromised</b>	5		
<b>Advisories:</b> <a href="#">See page 152</a> – 2005 Potential Separation of Interconnect Wires		Electrical Component	1		
		Electrical Interconnect	4		
		<i>(2 malfunction related to advisory)</i>			




**Sigma 200 SR SSR203**

Product Characteristics

US Market Release	Sep-99	Malfunctions (US)	6	NBG Code	SSI/R
Registered US Implants	12,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PJG
Estimated Active US Implants	4,000	<b>Therapy Function Compromised</b>	6	Estimated Longevity	<a href="#">See page 79</a>
Normal Battery Depletions (US)	65	Electrical Interconnect	6		
		<i>(4 malfunctions related to advisory)</i>			

**Advisories:** [See page 152](#) – 2005 Potential Separation of Interconnect Wires




**Sigma 300 DR SDR303, SDR306**

Product Characteristics

US Market Release	Aug-99	Malfunctions (US)	76	NBG Code	DDD/RO
Registered US Implants	107,000	<b>Therapy Function Not Compromised</b>	5	Serial Number Prefix	PJD, PJE
Estimated Active US Implants	53,000	Electrical Component	4	Estimated Longevity	<a href="#">See page 79</a>
Normal Battery Depletions (US)	394	Possible Early Battery Depletion	1		
		<b>Therapy Function Compromised</b>	71		
		Electrical Component	6		
		Electrical Interconnect	65		
		<i>(28 malfunctions related to advisory)</i>			

**Advisories:** [See page 152](#) – 2005 Potential Separation of Interconnect Wires






**Sigma 300 SR** SSR303, SSR306

Product Characteristics


US Market Release	Sep-99	Malfunctions (US)	15	NBG Code	SSI/R
Registered US Implants	54,000	<b>Therapy Function Not Compromised</b>	2	Serial Number Prefix	PJG, PJH
Estimated Active US Implants	20,000	Electrical Component	1	Estimated Longevity	<a href="#">See page 79</a>
Normal Battery Depletions (US)	191	Electrical Interconnect	1		
<b>Advisories:</b> <a href="#">See page 152</a> – 2005 Potential Separation of Interconnect Wires		<b>Therapy Function Compromised</b>	13		
		Electrical Component	3		
		Electrical Interconnect	10		
		<i>(5 malfunctions related to advisory)</i>			



**Sigma 300 VDD** SVDD303

Product Characteristics

US Market Release	Sep-99	Malfunctions (US)	0	NBG Code	VDDD
Registered US Implants	1,000	<b>Therapy Function Not Compromised</b>	0	Serial Number Prefix	PJD
Estimated Active US Implants	200	<b>Therapy Function Compromised</b>	0	Estimated Longevity	<a href="#">See page 79</a>
Normal Battery Depletions (US)	5				
<b>Advisories:</b> <a href="#">See page 152</a> – 2005 Potential Separation of Interconnect Wires					



**Thera-i D 7964i, 7965i, 7966i**


Product Characteristics

US Market Release	Oct-95
Registered US Implants	3,000
Estimated Active US Implants	300
Normal Battery Depletions (US)	195
Advisories	None

Malfunctions (US) 1

NBG Code	DDDCO
Serial Number Prefix	PDE, PDF, PDG

Estimated Longevity [See page 79](#)



**Thera-i DR 7960i, 7961i, 7962i**


Product Characteristics

US Market Release	Oct-95
Registered US Implants	122,000
Estimated Active US Implants	7,000
Normal Battery Depletions (US)	10,066
Advisories	None

Malfunctions (US) 50

NBG Code	DDD/RO
Serial Number Prefix	PDB, PDC, PDD

Estimated Longevity [See page 79](#)




**Thera-i DR 7968i**

Product Characteristics

US Market Release	Jul-96
Registered US Implants	4,000
Estimated Active US Implants	100
Normal Battery Depletions (US)	295
Advisories	None

Malfunctions (US) 3

NBG Code	DDD/RO
Serial Number Prefix	PGH
Estimated Longevity	<a href="#">See page 79</a>




**Thera-i S 8964i, 8965i, 8966i**

Product Characteristics

US Market Release	Oct-95
Registered US Implants	4,000
Estimated Active US Implants	500
Normal Battery Depletions (US)	84
Advisories	None

Malfunctions (US) 1

NBG Code	SSIR
Serial Number Prefix	PDY, PEA, PEB
Estimated Longevity	<a href="#">See page 79</a>




**Thera-i SR 8960i, 8961i, 8962i**

Product Characteristics

US Market Release	Oct-95
Registered US Implants	50,000
Estimated Active US Implants	5,000
Normal Battery Depletions (US)	2,012
Advisories	None

Malfunctions (US) 7

NBG Code	SSIR
Serial Number Prefix	PDU, PDV, PDW
Estimated Longevity	<a href="#">See page 79</a>




**Thera-i VDD 8968i**

Product Characteristics

US Market Release	Mar-96
Registered US Implants	5,000
Estimated Active US Implants	1,000
Normal Battery Depletions (US)	196
Advisories	None

Malfunctions (US) 0


NBG Code	VDD
Serial Number Prefix	PEC
Estimated Longevity	<a href="#">See page 79</a>



Versa DR VEDR01

Product Characteristics

US Market Release	Jul-06	Malfunctions (US)	2	NBG Code	DDDR
Registered US Implants	38,000	<b>Therapy Function Not Compromised</b>	2	Serial Number Prefix	PWH, NWH
Estimated Active US Implants	34,000	Electrical Component	2	Estimated Longevity	<a href="#">See page 79</a>
Normal Battery Depletions (US)	1	<b>Therapy Function Compromised</b>	0		
Advisories	None				



**Device Survival Summary** (95% Confidence Interval)

The following table shows IPG device survival estimates with 95% confidence intervals. Estimates are shown both with and without normal battery depletions included.

Family	Model Number	US Market Release	Registered US Implants	Estimated Active US Implants	Normal Battery Depletions (US)	Malfunctions (US)		Device Survival Probability (%)												
						Therapy Function Compromised	Therapy Function Not Compromised	Years After Implant												
						4	6	10 yr	12 yr	14 yr	16 yr									
Adapta DR	ADDR01, ADDR03, ADDR06, ADD01	Jul-06	116,000	104,000	1	4 + 6 = 10	Excluding Normal Battery Depletion	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0		
							Including Normal Battery Depletion	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0		
							Excluding Normal Battery Depletion	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0		
Adapta DR	ADDR11	Jul-06	11,000	10,000	0	0 + 0 = 0	Excluding Normal Battery Depletion	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0		
Adapta DR	ADDR51	Jul-06	10,000	9,000	1	1 + 0 = 1	Excluding Normal Battery Depletion	100.0	99.9	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0		
Adapta SR	ADSR01, ADSR03, ADSR06	Jul-06	22,000	19,000	1	0 + 0 = 0	Excluding Normal Battery Depletion	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0		
Adapta VDD	ADVDD01	Jul-06	400	400	0	0 + 0 = 0	Excluding Normal Battery Depletion	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0		
AT500	AT501, 7253	Mar-03	11,000	5,000	481	5 + 5 = 10	Excluding Normal Battery Depletion	100.0	99.9	99.9	100.0	99.9	99.9	100.0	99.9	99.9	99.9	100.0		
							Including Normal Battery Depletion	100.0	99.9	99.9	100.0	99.9	99.9	100.0	99.9	99.9	100.0	99.9	99.9	
EnPulse DR	E1DR01, E1DR03, E1DR06	Dec-03	7,000	4,000	26	0 + 1 = 1	Excluding Normal Battery Depletion	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0		
EnPulse DR	E1DR21	Dec-03	2,000	1,000	23	0 + 0 = 0	Excluding Normal Battery Depletion	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0		

**Device Survival Summary** continued

		Device Survival Probability (%)													
		Years After Implant													
		1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr	12 yr	14 yr	16 yr		
EnPulse 2 DR	E2DR01, E2DR03, E2DR06	US Market Release	Feb-04	101,000	70,000	81	Total	3 + 10 = 13	Therapy Function Compromised	Therapy Function Not Compromised	Normal Battery Depletions (US)	81	Estimated Active US Implants	70,000	81
		Model Number													
		Family	EnPulse 2 DR												
EnPulse 2 DR	E2DR21	US Market Release	Feb-04	12,000	8,000	48	Total	1 + 0 = 1	Therapy Function Compromised	Therapy Function Not Compromised	Normal Battery Depletions (US)	48	Estimated Active US Implants	8,000	48
		Model Number													
		Family	EnPulse 2 DR												
EnPulse 2 DR	E2DR31, E2DR33	US Market Release	Feb-04	1,000	400	0	Total	0 + 0 = 0	Therapy Function Compromised	Therapy Function Not Compromised	Normal Battery Depletions (US)	0	Estimated Active US Implants	400	0
		Model Number													
		Family	EnPulse 2 DR												
EnPulse SR	E2SR01, E2SR03, E2SR06	US Market Release	Dec-03	25,000	15,000	42	Total	1 + 2 = 3	Therapy Function Compromised	Therapy Function Not Compromised	Normal Battery Depletions (US)	42	Estimated Active US Implants	15,000	42
		Model Number													
		Family	EnPulse SR												
EnPulse 2 VDD	E2VDD01	US Market Release	Dec-03	1,000	500	1	Total	0 + 0 = 0	Therapy Function Compromised	Therapy Function Not Compromised	Normal Battery Depletions (US)	1	Estimated Active US Implants	500	1
		Model Number													
		Family	EnPulse 2 VDD												
EnRhythm DR	P1501DR	US Market Release	May-05	78,000	62,000	4	Total	25 + 5 = 30	Therapy Function Compromised	Therapy Function Not Compromised	Normal Battery Depletions (US)	4	Estimated Active US Implants	62,000	4
		Model Number													
		Family	EnRhythm DR												
Kappa 400 DR	KDR401, KDR403	US Market Release	Jan-98	47,000	5,000	5,271	Total	9 + 13 = 22	Therapy Function Compromised	Therapy Function Not Compromised	Normal Battery Depletions (US)	5,271	Estimated Active US Implants	5,000	5,271
		Model Number													
		Family	Kappa 400 DR												
Kappa 400 SR	KSR401, KSR403	US Market Release	Feb-98	15,000	2,000	870	Total	1 + 4 = 5	Therapy Function Compromised	Therapy Function Not Compromised	Normal Battery Depletions (US)	870	Estimated Active US Implants	2,000	870
		Model Number													
		Family	Kappa 400 SR												

continued

Family	Model Number	US Market Release	Registered US Implants	Estimated Active US Implants	Normal Battery Depletions (US)	Malfunctions (US)		Device Survival Probability (%)	Years After Implant																																
						Therapy Function Compromised	Therapy Function Not Compromised		1 yr		2 yr		3 yr		4 yr		5 yr		6 yr		7 yr		8 yr		10 yr		12 yr		14 yr		16 yr										
									100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.1	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.0/-0.1	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.1	100.0 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.0/-0.1	100.0 +0.0/-0.0	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1			
Kappa 600 DR	KDR601, KDR603, KDR606 Advisories: See page 154 – 2002 Potential Fractured Power Supply Wires	Jan-99	24,000	2,000	2,417	30	3 = 33	Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.1	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.0/-0.1	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.0/-0.1	100.0 +0.0/-0.0	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1			
						(12) + (0) = (12) (advisory-related subset)	Including Normal Battery Depletion	99.9 +0.0/-0.1	99.8 +0.1/-0.1	99.5 +0.1/-0.1	98.8 +0.2/-0.2	96.8 +0.3/-0.3	87.9 +0.6/-0.6	15.5 +1.1/-1.1	58.9 +1.1/-1.1																										
Kappa 600 DR	KDR651, KDR653 Advisories: See page 154 – 2002 Potential Fractured Power Supply Wires	Mar-01	14,000	5,000	569	17	2 = 19	Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.1	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.8 +0.1/-0.1	99.8 +0.1/-0.1	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1	99.7 +0.1/-0.1		
						(1) + (0) = (1) (advisory-related subset)	Including Normal Battery Depletion	99.9 +0.0/-0.1	99.9 +0.1/-0.1	99.4 +0.1/-0.2	98.1 +0.3/-0.3	81.4 +1.1/-1.2	66.8 +2.8/-3.0																												
Kappa 700 D	KD701, KD703, KD706 Advisories: See page 154 – 2002 Potential Fractured Power Supply Wires	Jan-99	300	100	11	0	0 = 0	Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0			
						(0) + (0) = (0) (advisory-related subset)	Including Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.0 +0.8/-3.1	97.8 +1.4/-3.5	93.8 +2.3/-4.5																													
Kappa 700 DR	KDR701, KDR703, KDR706 Advisories: See page 154 – 2002 Potential Fractured Power Supply Wires	Feb-99	192,000	54,000	12,539	323	28 = 351	Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.0/-0.0	99.9 +0.0/-0.1	99.8 +0.1/-0.1	99.9 +0.0/-0.0	99.9 +0.0/-0.0	99.9 +0.0/-0.1	99.8 +0.0/-0.0	99.8 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.1/-0.1	99.5 +0.8/-0.8	9.4 +0.8/-0.8	54.6 +0.5/-0.6	84.6 +0.3/-0.3	95.3 +0.1/-0.1	95.3 +0.1/-0.1	99.9 +0.0/-0.0	99.9 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	
						(133) + (0) = (133) (advisory-related subset)	Including Normal Battery Depletion	99.9 +0.0/-0.0	99.8 +0.1/-0.1	99.2 +0.0/-0.0	98.2 +0.1/-0.1	84.6 +0.3/-0.3																													
Kappa 700 DR	KDR721 Advisories: See page 154 – 2002 Potential Fractured Power Supply Wires	Feb-99	10,000	100	1,237	4	1 = 5	Excluding Normal Battery Depletion	100.0 +0.0/-0.1	100.0 +0.0/-0.1	100.0 +0.0/-0.1	100.0 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	100.0 +0.0/-0.1	100.0 +0.0/-0.1	100.0 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1		
						(4) + (0) = (4) (advisory-related subset)	Including Normal Battery Depletion	99.9 +0.0/-0.1	99.7 +0.1/-0.2	96.8 +0.4/-0.5	91.1 +0.7/-0.8	20.0 +2.2/-2.1																													
Kappa 700 SR	KSR701, KSR703, KSR706 Advisories: See page 154 – 2002 Potential Fractured Power Supply Wires	Feb-99	55,000	13,000	2,265	14	3 = 17	Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.8 +0.1/-0.1	99.8 +0.1/-0.1	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	
						(3) + (0) = (3) (advisory-related subset)	Including Normal Battery Depletion	99.9 +0.0/-0.0	99.8 +0.1/-0.1	98.8 +0.1/-0.1	97.4 +0.2/-0.2	85.2 +0.6/-0.6	62.0 +1.1/-1.2	24.0 +2.6/-2.5																											
Kappa 700 VDD	KVDD701 Advisories: See page 154 – 2002 Potential Fractured Power Supply Wires	Jan-99	2,000	100	149	3	0 = 3	Excluding Normal Battery Depletion	99.9 +0.1/-0.4	99.9 +0.1/-0.4	99.8 +0.1/-0.5	99.8 +0.1/-0.5	99.8 +0.1/-0.4	99.8 +0.1/-0.4	99.8 +0.1/-0.5	99.8 +0.1/-0.5	99.7 +0.2/-0.4	99.7 +0.2/-0.4	99.8 +0.1/-0.4	99.8 +0.1/-0.4	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6	99.7 +0.2/-0.6
						(3) + (0) = (3) (advisory-related subset)	Including Normal Battery Depletion	99.7 +0.2/-0.4	99.7 +0.2/-0.4	98.9 +0.5/-0.8	98.5 +0.6/-0.9	69.7 +3.4/-3.7																													



Device Survival Summary continued

Family		Model Number	US Market Release	Registered US Implants	Estimated Active US Implants	Normal Battery Depletions	Malfunctions		Device Survival Probability (%)											
							Therapy Function Compromised	Therapy Function Not Compromised	Total	Years After Implant										
							1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr	12 yr	14 yr	16 yr		
Kappa 800 DR	KDR801, KDR803	27	Jan-02	4,000	2,000	Normal Battery Depletions	1	0	=	1	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.1/-0.6 at 75 mo	99.9 +0.1/-0.6 at 75 mo					
							Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.1/-0.6 at 75 mo	99.9 +0.1/-0.6 at 75 mo								
Kappa 900 DR	KDR901, KDR903, KDR906	907	Jan-02	125,000	70,000	Normal Battery Depletions	24	12	=	36	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.1/-0.6 at 82 mo	99.9 +0.1/-0.6 at 82 mo					
							Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.1/-0.6 at 82 mo	99.9 +0.1/-0.6 at 82 mo								
Kappa 900 SR	KSR901, KSR903, KSR906	240	Jan-02	37,000	17,000	Normal Battery Depletions	3	8	=	11	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.8 +0.1/-0.6 at 81 mo	99.8 +0.1/-0.6 at 81 mo					
							Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.8 +0.1/-0.6 at 81 mo	99.8 +0.1/-0.6 at 81 mo								
Kappa 900 VDD	KVDD901	9	Jan-02	1,000	300	Normal Battery Depletions	0	0	=	0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0					
							Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0								
Kappa 920 DR	KDR921	902	Jan-02	16,000	6,000	Normal Battery Depletions	3	0	=	3	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.1/-0.6 at 77 mo	99.9 +0.1/-0.6 at 77 mo					
							Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.1/-0.6 at 77 mo	99.9 +0.1/-0.6 at 77 mo								
Legend	8416, 8417, 8417M, 8418, 8419	2,900	Aug-89	57,000	2,000	Normal Battery Depletions	—	—	=	143	99.9 +0.0/-0.0	99.9 +0.0/-0.0	99.9 +0.0/-0.0	99.6 +0.1/-0.6 at 217 mo	99.6 +0.1/-0.6 at 217 mo					
							Excluding Normal Battery Depletion	99.9 +0.0/-0.0	99.9 +0.0/-0.0	99.9 +0.0/-0.0	99.6 +0.1/-0.6 at 217 mo	99.6 +0.1/-0.6 at 217 mo								
Legend II	8424, 8426, 8427	2,203	Nov-91	59,000	4,000	Normal Battery Depletions	—	—	=	36	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.0	99.9 +0.0/-0.0					
							Excluding Normal Battery Depletion	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	99.9 +0.0/-0.0	99.9 +0.0/-0.0								

continued

Device Survival Summary continued

		Malfunctions					Device Survival Probability (%)															
Family	Model Number	US Market Release	Registered US Implanters	Estimated Active US Implants	Normal Battery Depletions	Therapy Function Compromised	Therapy Function Not Compromised	Total	Years After Implant													
									1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr	12 yr	14 yr	16 yr		
Minix/ Minix ST	8330, 8331, 8331M, 8340, 8341, 8341M, 8342	Dec-89	58,000	4,000	1,594	—	—	49	Excluding Normal Battery Depletion	100.0	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9	99.9
										+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0
									Including Normal Battery Depletion	99.9	99.7	99.5	99.2	98.6	97.6	95.1	91.8	87.0	83.6	80.4	71.7	
										+0.0/-0.0	+0.0/-0.0	+0.1/-0.1	+0.1/-0.1	+0.1/-0.1	+0.2/-0.2	+0.3/-0.3	+0.4/-0.4	+0.5/-0.5	+0.6/-0.6	+0.7/-0.7	+1.6/-1.7	
										100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	
										+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.1	+0.0/-0.1	+0.0/-0.1	+0.0/-0.1	+0.0/-0.1	+0.0/-0.1	+0.0/-0.1	+0.0/-0.1	+0.0/-0.1	
										100.0	100.0	99.8	99.5	98.9	98.2	96.9	94.7	89.2	81.9	69.0	51.8	
										+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.1/-0.1	+0.2/-0.2	+0.3/-0.3	+0.4/-0.4	+0.5/-0.5	+0.7/-0.8	+1.0/-1.1	+1.5/-1.6	+3.0/-3.1	
										100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	
										+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.1	+0.0/-0.1	+0.0/-0.1	
										99.9	99.8	99.6	99.4	98.8	97.9	96.3	92.5	67.4	25.8			
										+0.0/-0.0	+0.0/-0.0	+0.1/-0.1	+0.1/-0.1	+0.2/-0.2	+0.3/-0.3	+0.4/-0.4	+0.5/-0.5	+1.1/-1.1	+2.3/-2.3			
										100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	
										+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.1	+0.0/-0.1	+0.0/-0.1	
										99.9	99.8	99.5	99.1	98.4	97.3	94.9	92.2	75.9	46.4			
										+0.0/-0.1	+0.1/-0.1	+0.1/-0.1	+0.2/-0.2	+0.3/-0.3	+0.4/-0.4	+0.5/-0.5	+0.6/-0.7	+1.4/-1.4	+3.4/-3.5			
										100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	
										+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.1/-0.4	+0.1/-0.4	+0.1/-0.4	+0.1/-0.4	+0.1/-0.4	+0.1/-0.4	+0.1/-0.4	+0.1/-0.4	
										99.9	99.8	99.8	99.8	99.1	98.9	98.1	97.2	93.5	89.7	86.4	82.9	
										+0.1/-0.2	+0.1/-0.2	+0.1/-0.2	+0.1/-0.2	+0.4/-0.6	+0.4/-0.7	+0.6/-0.9	+0.8/-1.2	+1.5/-2.0	+2.2/-2.8	+4.6/-4.9		
										100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	
										+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	
										99.9	99.7	99.4	98.8	98.6	97.6	96.8	95.1	75.2	52.9			
										+0.1/-0.2	+0.2/-0.3	+0.3/-0.4	+0.4/-0.6	+0.4/-0.6	+0.6/-0.8	+0.8/-1.0	+1.0/-1.3	+2.9/-3.2	+4.6/-4.9			
										100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	100.0	
										+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	+0.0/-0.0	
										99.9	99.9	99.7	99.4	98.8	97.6	96.4	92.5	67.9	23.8			
										+0.0/-0.0	+0.0/-0.0	+0.1/-0.1	+0.1/-0.1	+0.1/-0.1	+0.2/-0.2	+0.3/-0.3	+0.4/-0.4	+0.9/-0.9	+2.3/-2.3			

continued



Device Survival Summary continued

		Device Survival Probability (%)												
		Years After Implant												
		1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr	12 yr	14 yr	16 yr	
Sigma 300 DR	Model Number	SDR303, SDR306	71 + 5 = 76		394		394		394		394		394	
	US Market Release	Aug-99	107,000		53,000		53,000		53,000		53,000		53,000	
	Registered US Implants		(15) + (0) = (15)		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion	
Sigma 300 SR	Model Number	SSR303, SSR306	13 + 2 = 15		191		191		191		191		191	
	US Market Release	Sep-99	54,000		20,000		20,000		20,000		20,000		20,000	
	Registered US Implants		(4) + (0) = (4)		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion	
Sigma 300 VDD	Model Number	SVDD303	0 + 0 = 0		5		5		5		5		5	
	US Market Release	Sep-99	1,000		200		200		200		200		200	
	Registered US Implants		(0) + (0) = (0)		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion	
Thera-iD	Model Number	7964i, 7965i, 7966i	—		1		1		1		1		1	
	US Market Release	Oct-95	3,000		300		300		300		300		300	
	Registered US Implants		—		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion	
Thera-iDR	Model Number	7960i, 7961i, 7962i	—		50		50		50		50		50	
	US Market Release	Oct-95	122,000		7,000		7,000		7,000		7,000		7,000	
	Registered US Implants		—		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion	
Thera-iDR	Model Number	7968i	—		3		3		3		3		3	
	US Market Release	Jul-96	4,000		100		100		100		100		100	
	Registered US Implants		—		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion	
Thera-iS	Model Number	8964i, 8965i, 8966i	—		1		1		1		1		1	
	US Market Release	Oct-95	4,000		500		500		500		500		500	
	Registered US Implants		—		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion		Including Normal Battery Depletion		Excluding Normal Battery Depletion	

Device Survival Summary continued

		Malfunctions								Device Survival Probability (%)												
Family	Model Number	US Market Release	Registered US Implants	Estimated Active US Implants	Normal Battery Depletions	Therapy Function		Total	Years After Implant													
						Complimented	Not Complimented		1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr	12 yr	14 yr	16 yr		
Thera-i SR	8960i, 8961i, 8962i	Oct-95	50,000	5,000	2,012	—	—	7	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0
									100.0 +0.0/-0.0	99.9 +0.0/-0.0	99.7 +0.1/-0.1	99.4 +0.1/-0.1	98.9 +0.1/-0.1	98.2 +0.2/-0.2	96.9 +0.2/-0.2	94.5 +0.3/-0.3	80.6 +0.7/-0.7	52.5 +1.5/-1.5	43.2 +2.3/-2.4	at 153 mo		
Thera-i VDD	8968i	Mar-96	5,000	1,000	196	—	—	0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0
									99.9 +0.1/-0.1	99.8 +0.1/-0.2	99.6 +0.2/-0.3	99.5 +0.2/-0.3	99.2 +0.3/-0.4	99.1 +0.3/-0.4	99.0 +0.3/-0.4	98.3 +0.4/-0.6	88.7 +1.5/-1.7	53.4 +4.7/-5.0				
Versa DR	VEDR01	Jul-06	38,000	34,000	1	0	2	2	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0	100.0 +0.0/-0.0
									100.0 +0.0/-0.0	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1	99.9 +0.0/-0.1

## Reference Chart

The longevity estimates provided are mean values calculated for the parameters given. The longevity estimates shown here assume a lower rate of 60 ppm, 100% pacing, and pulse width of 0.4 ms unless noted otherwise. The actual longevity achieved for any device while implanted will depend on the actual programmed parameters and patient factors, and may differ significantly from these estimates. **The elective replacement time is indicated via telemetry indication, and rate and mode change to 65 ppm and VVI respectively (VOO/65 with magnet), unless noted otherwise.**

Family	Model Number	Amplitude Setting	Estimated Longevity		Elective Replacement Indicators
			500 Lead $\Omega$	1000 Lead $\Omega$	
Adapta DR	ADDR01, ADDR03, ADDR06, ADD01	Low 2.5 V (A, RV)	7.4	8.2	**
		Nominal 3.5 V (A, RV)	6.0	7.3	
		High 5.0 V (A, RV)	4.5	6.0	
Adapta DR	ADDRS1	Low 2.5 V (A, RV)	5.5	6.1	**
		Nominal 3.5 V (A, RV)	4.3	5.4	
		High 5.0 V (A, RV)	3.2	4.4	
Adapta DR	ADDRL1	Low 2.5 V (A, RV)	9.1	10.1	**
		Nominal 3.5 V (A, RV)	7.4	9.0	
		High 5.0 V (A, RV)	5.4	7.3	
Adapta SR	ADSR01, ADSR03, ADSR06	Low 2.5 V (RV)	7.3	7.8	**
		Nominal 3.5 V (RV)	6.4	7.4	
		High 5.0 V (RV)	5.0	6.2	
Adapta VDD	ADVDD01	Low 2.5 V (RV)	6.2	6.5	**
		Nominal 3.5 V (RV)	5.5	6.2	
		High 5.0 V (RV)	4.4	5.4	
AT500	AT501, 7253	Low 2.0 V (A, RV)	7.7	8.3	Telemetry indication. Pacing mode and rate (magnet and non-magnet) as programmed.
		Nominal 3.0 V (A, RV)	5.8	7.0	
		High 5.0 V (A, RV)	3.7	5.2	
EnPulse DR	E1DR01, E1DR03, E1DR06	Low 2.5 V (A, RV)	7.5	8.5	**
		Nominal 3.5 V (A, RV)	6.2	7.6	
		High 5.0 V (A, RV)	4.4	5.9	
EnPulse DR	E1DR21	Low 2.5 V (A, RV)	5.4	6.0	**
		Nominal 3.5 V (A, RV)	4.3	5.4	
		High 5.0 V (A, RV)	3.0	4.2	
EnPulse 2 DR	E2DR01, E2DR03, E2DR06	Low 2.5 V (A, RV)	7.5	8.5	**
		Nominal 3.5 V (A, RV)	6.2	7.6	
		High 5.0 V (A, RV)	4.4	5.9	
EnPulse 2 DR	E2DR21	Low 2.5 V (A, RV)	5.4	6.0	**
		Nominal 3.5 V (A, RV)	4.3	5.4	
		High 5.0 V (A, RV)	3.0	4.2	
EnPulse 2 DR	E2DR31, E2DR33	Low 2.5 V (A, RV)	9.0	10.1	**
		Nominal 3.5 V (A, RV)	7.4	9.1	
		High 5.0 V (A, RV)	5.2	7.1	
EnPulse 2 SR	E2SR01, E2SR03, E2SR06	Low 2.5 V (A, RV)	7.2	7.7	**
		Nominal 3.5 V (A, RV)	6.3	7.3	
		High 5.0 V (A, RV)	4.8	6.1	
EnPulse 2 VDD	E2VDD01	Low 2.5 V (RV)	6.1	6.5	**
		Nominal 3.5 V (RV)	5.5	6.2	
		High 5.0 V (RV)	4.3	5.4	
EnRhythm DR	P1501DR	Low 2.5 V (A, RV)	10.6	12.3	**
		Nominal 3.5 V (A, RV)	8.0	10.3	
		High 5.0 V (A, RV)	5.4	7.8	
Kappa 400 DR	KDR401, KDR403	Low 2.5 V (A, RV)	7.8	8.5	**
		Nominal 3.5 V (A, RV)	6.4	7.5	
		High 5.0 V (A, RV)	5.1	6.5	
Kappa 400 SR	KSR401, KSR403	Low 2.5 V (RV)	7.9	8.4	**
		Nominal 3.5 V (RV)	6.9	7.7	
		High 5.0 V (RV)	5.8	7.0	
Kappa 600 DR	KDR601, KDR603, KDR606	Low 2.5 V (A, RV)	7.7	8.6	**
		Nominal 3.5 V (A, RV)	6.3	7.7	
		High 5.0 V (A, RV)	4.4	6.0	
Kappa 600 DR	KDR651, KDR653	Low 2.5 V (A, RV)	7.7	8.6	**
		Nominal 3.5 V (A, RV)	6.3	7.7	
		High 5.0 V (A, RV)	4.4	6.0	

\*\*Telemetry indication. Rate and mode change to 65 ppm and VVI respectively (VOO/65 with magnet).

continued



Reference Chart continued

Family	Model Number	Amplitude Setting	Estimated Longevity		Elective Replacement Indicators
			500 Lead $\Omega$	1000 Lead $\Omega$	
Kappa 700 D	KDR701, KD703, KD706	Low 2.5 V (A, RV)	7.7	8.6	**
		Nominal 3.5 V (A, RV)	6.3	7.7	
		High 5.0 V (A, RV)	4.4	6.0	
Kappa 700 DR	KDR701, KDR703, KDR706	Low 2.5 V (A, RV)	7.7	8.6	**
		Nominal 3.5 V (A, RV)	6.3	7.7	
		High 5.0 V (A, RV)	4.4	6.0	
Kappa 700 DR	KDR721	Low 2.5 V (A, RV)	5.5	6.1	**
		Nominal 3.5 V (A, RV)	4.4	5.5	
		High 5.0 V (A, RV)	3.0	4.2	
Kappa 700 SR	KSR701, KSR703, KSR706	Low 2.5 V (RV)	7.4	7.9	**
		Nominal 3.5 V (RV)	6.5	7.5	
		High 5.0 V (RV)	4.9	6.2	
Kappa 700 VDD	KVDD701	Low 2.5 V (RV)	6.2	6.6	**
		Nominal 3.5 V (RV)	5.6	6.3	
		High 5.0 V (RV)	4.4	5.3	
Kappa 800 DR	KDR801, KDR803	Low 2.5 V (A, RV)	7.7	8.6	**
		Nominal 3.5 V (A, RV)	6.3	7.7	
		High 5.0 V (A, RV)	4.4	6.0	
Kappa 900 DR	KDR901, KDR903, KDR906	Low 2.5 V (A, RV)	7.7	8.6	**
		Nominal 3.5 V (A, RV)	6.3	7.7	
		High 5.0 V (A, RV)	4.4	6.0	
Kappa 920 DR	KDR921	Low 2.5 V (A, RV)	5.5	6.1	**
		Nominal 3.5 V (A, RV)	4.4	5.5	
		High 5.0 V (A, RV)	3.0	4.3	
Kappa 900 SR	KSR901, KSR903, KSR906	Low 2.5 V (RV)	7.3	7.9	**
		Nominal 3.5 V (RV)	6.4	7.4	
		High 5.0 V (RV)	4.9	6.1	
Kappa 900 VDD	KVDD901	Low 2.5 V (RV)	6.2	6.6	**
		Nominal 3.5 V (RV)	5.6	6.3	
		High 5.0 V (RV)	4.4	5.4	
Legend	8416, 8417, 8417M, 8418, 8419	Low 2.5 V (RV)	15.6	17.7	If programmed to non-rate responsive mode (e.g., VVI), rate decrease of 10% from programmed rate. Telemetry indication. If programmed to rate responsive mode (e.g., VVIR), rate change to 65 ppm and mode change to VVI. Telemetry indication.
		Nominal 3.3 V (RV)	11.3	14.5	
		High 5.0 V (RV)	9.0	12.5	
Legend II	8424, 8426, 8427	Low 2.5 V, 0.36 ms (RV)	12.9	14.5	If programmed to non-rate responsive mode (e.g., VVI), rate decrease of 10% from programmed rate. Telemetry indication. If programmed to rate responsive mode (e.g., VVIR), rate change to 65 ppm and mode change to VVI. Telemetry indication.
		Nominal 3.3 V, 0.36 ms (RV)	9.4	11.8	
		High 5.0 V, 0.36 ms (RV)	7.8	10.5	
Minix	8340, 8341, 8341M, 8342	Low 2.5 V (RV)	14.9	17.3	Telemetry indication. Rate decrease of 10% from programmed rate.
		Nominal 3.3 V (RV)	10.2	13.6	
		High 5.0 V (RV)	7.9	11.3	
Minix ST	8330, 8331, 8331M	Low 2.5 V (RV)	14.9	17.3	Telemetry indication. Rate decrease of 10% from programmed rate.
		Nominal 5.0 V (RV)	7.9	11.4	
		High 8.0 V (RV)	4.0	7.0	
Minuet	7107, 7108	Low 2.5 V, 0.36 ms (A, RV)	12.5	15.6	**
		Nominal 4.0 V, 0.36 ms (A, RV)	7.7	10.9	
		High 5.0 V, 0.36 ms (A, RV)	4.7	7.6	
Preva DR	7088, 7089	Low 2.5 V (A, RV)	9.9	11.3	**
		Nominal 3.5 V (A, RV)	7.4	9.4	
		High 5.0 V (A, RV)	5.4	7.5	
Preva SR	8088, 8089	Low 2.5 V (RV)	9.8	10.7	**
		Nominal 3.5 V (RV)	8.0	9.5	
		High 5.0 V (RV)	6.4	8.1	
Prevail S	8085, 8086	Low 2.5 V, 0.42 ms (RV)	16.4	19.4	Telemetry indication. Rate decrease of 10% from programmed rate.
		Nominal 3.3 V, 0.42 ms (RV)	10.8	14.4	
		High 5.0 V, 0.42 ms (RV)	8.6	12.4	

\*\*Telemetry indication. Rate and mode change to 65 ppm and VVI respectively (VOO/65 with magnet).

continued

Reference Chart continued

Family	Model Number	Amplitude Setting	Estimated Longevity		Elective Replacement Indicators
			500 Lead $\Omega$	1000 Lead $\Omega$	
Prodigy D	7864, 7865, 7866	Low 2.5 V (A, RV)	10.0	11.4	**
		Nominal 3.5 V (A, RV)	7.4	9.5	
		High 5.0 V (A, RV)	5.4	7.6	
Prodigy DR	7860, 7861, 7862	Low 2.5 V (A, RV)	9.9	11.3	**
		Nominal 3.5 V (A, RV)	7.4	9.4	
		High 5.0 V (A, RV)	5.4	7.5	
Prodigy S	8164, 8165, 8166	Low 2.5 V (RV)	10.0	10.9	**
		Nominal 3.5 V (RV)	8.1	9.6	
		High 5.0 V (RV)	6.4	8.2	
Prodigy SR	8158, 8160, 8161, 8162	Low 2.5 V (RV)	9.8	10.7	**
		Nominal 3.5 V (RV)	8.0	9.5	
		High 5.0 V (RV)	6.4	8.1	
Sensia DR	SEDR01, SED01	Low 2.5 V (A, RV)	7.5	8.3	**
		Nominal 3.5 V (A, RV)	6.1	7.4	
		High 5.0 V (A, RV)	4.5	6.0	
Sensia DR	SEDR1	Low 2.5 V (A, RV)	9.1	10.1	**
		Nominal 3.5 V (A, RV)	7.4	9.0	
		High 5.0 V (A, RV)	5.4	7.3	
Sensia SR	SESR01, SES01	Low 2.5 V (RV)	7.3	7.8	**
		Nominal 3.5 V (RV)	6.4	7.4	
		High 5.0 V (RV)	5.0	6.2	
Sigma 100 S	SS103, SS106	Low 2.5 V (RV)	10.1	11.1	**
		Nominal 3.5 V (RV)	8.2	9.8	
		High 5.0 V (RV)	6.4	8.4	
Sigma 200 DR	SDR203	Low 2.5 V (A, RV)	10.1	11.7	**
		Nominal 3.5 V (A, RV)	7.5	9.6	
		High 5.0 V (A, RV)	5.5	7.8	
Sigma 200 SR	SSR203	Low 2.5 V (RV)	10.1	11.1	**
		Nominal 3.5 V (RV)	8.2	9.8	
		High 5.0 V (RV)	6.4	8.4	
Sigma 300 DR	SDR303, SDR306	Low 2.5 V (A, RV)	10.1	11.7	**
		Nominal 3.5 V (A, RV)	7.5	9.6	
		High 5.0 V (A, RV)	5.5	7.8	
Sigma 300 SR	SSR303, SSR306	Low 2.5 V (RV)	10.1	11.1	**
		Nominal 3.5 V (RV)	8.2	9.8	
		High 5.0 V (RV)	6.4	8.4	
Sigma 300 VDD	SVDD303	Low 2.5 V (RV)	8.9	9.7	**
		Nominal 3.5 V (RV)	7.3	8.6	
		High 5.0 V (RV)	5.8	7.4	
Thera-i D	7964i, 7965i, 7966i	Low 2.5 V (A, RV)	10.0	11.4	**
		Nominal 3.5 V (A, RV)	7.4	9.5	
		High 5.0 V (A, RV)	5.4	7.6	
Thera-i DR	7960i, 7961i, 7962i	Low 2.5 V (A, RV)	9.9	11.3	**
		Nominal 3.5 V (A, RV)	7.4	9.4	
		High 5.0 V (A, RV)	5.4	7.5	
Thera-i DR	7968i	Low 2.5 V (A, RV)	7.2	8.3	**
		Nominal 3.5 V (A, RV)	5.4	6.9	
		High 5.0 V (A, RV)	3.9	5.5	
Thera-i S	8964i, 8965i, 8966i	Low 2.5 V (RV)	10.0	10.9	**
		Nominal 3.5 V (RV)	8.1	9.6	
		High 5.0 V (RV)	6.4	8.2	
Thera-i SR	8960i, 8961i, 8962i	Low 2.5 V (RV)	9.8	10.7	**
		Nominal 3.5 V (RV)	8.0	9.5	
		High 5.0 V (RV)	6.4	8.1	
Thera-i VDD	8968i	Low 2.5 V (RV)	11.5	12.4	**
		Nominal 3.5 V (RV)	9.6	11.1	
		High 5.0 V (RV)	7.7	9.7	
Versa DR	VEDR01	Low 2.5 V (A, RV)	7.5	8.3	**
		Nominal 3.5 V (A, RV)	6.1	7.4	
		High 5.0 V (A, RV)	4.5	6.0	

\*\*Telemetry indication. Rate and mode change to 65 ppm and VVI respectively (VOO/65 with magnet).





# Method for Estimating Lead Performance

Medtronic CRDM has tracked lead survival for over 20 years with its multicenter, global chronic lead studies.

## Leads Performance Analysis

Implanted leads operate in the challenging biochemical environment of the human body and the body's response to foreign objects. Implanted leads are also subject to mechanical stresses associated with heart motion, body motion, and patient anatomy.

In this environment, pacemaker and defibrillation leads cannot be expected to last forever. While IPGs and ICDs have a battery that will deplete after a predictable length of time, a lead's longevity cannot be predicted, nor are there simple indicators that a lead is approaching the end of its service life. Therefore, regular monitoring while implanted, and evaluation of lead integrity upon IPG or ICD replacement, is necessary to determine if a lead may be approaching the end of its service life.

## Returned Product Analysis Shortfalls

All leads and lead segments returned to Medtronic are analyzed to determine whether or not they meet performance limits established by Medtronic. Although returned product analyses are valuable for gaining insight into lead failure mechanisms, this data cannot be used by itself for determining the survival probability of leads because only a small fraction of leads are explanted and returned for analysis. Additionally, those leads that are returned cannot be assumed to be statistically representative of the performance of the total population for a given lead model. Partial or total lead extraction can result in significant damage to a lead, making a definitive analysis of a suspected failure and its cause impossible. Thus, lead survival probabilities are more appropriately determined through a clinical surveillance study. Although returned product analysis results are presented in this report, Medtronic tracks lead survival through its System Longevity Study.

## System Longevity Study (SLS)

The SLS is a prospective, multicenter, global study designed to monitor the performance of market-released cardiac therapy products. Medtronic has been monitoring the performance of its cardiac therapy products for 26 years and has evaluated the performance of more than 75,000 leads, with data reported from 14 countries on four continents.

Patients are eligible for enrollment in the study if:

- 1 They are within 6 months post-implant of a Medtronic market-released lead connected to a market-released CRT, ICD, or IPG device, and the lead is used for a pacing, sensing, or defibrillation application, or
- 2 They participated in a qualifying study of a market-released Medtronic cardiac therapy product; complete implant and follow-up data are available; and the data is appropriately and legally released for use in the study.

## Lead Complications

The SLS complication criteria are defined below. These criteria do not, however, enable a lead integrity or "hardware" failure to be conclusively differentiated from other clinical events such as an undetected lead dislodgement, exit block, or concurrent pulse generator failure manifested as a sensing or capture problem.

A lead-related complication is considered to have occurred if at least one of the following clinical observations is reported and at least one of the following clinical actions is made 30 days or more after the implant.

## Clinical Observations

- Failure to capture
- Failure to sense/undersensing
- Oversensing
- Abnormal pacing impedance (based on lead model, but normal range is typically 200-3,000 ohms)
- Abnormal defibrillation impedance (based on lead model, but normal range is typically 20-200 ohms)
- Insulation breach, observed visually, that has degraded system performance
- Conductor fracture, observed visually or radiographically
- Extracardiac stimulation
- Cardiac perforation
- Lead dislodgement

continued

The Standard Actuarial Method is used to determine estimates of lead survival.

## Clinical Actions

- Lead surgically abandoned/capped
- Lead electrically abandoned/capped
- Lead explanted
- Lead replaced
- Polarity reprogrammed (i.e., bipolar to unipolar; unipolar to bipolar)
- Lead use continued based on medical judgment despite a known clinical performance issue
- Other lead-related surgery performed (e.g., lead mechanical alteration or unsuccessful repositioning)

**Note:** Successful lead repositioning is not a qualifying action.

## Methods

The performance of leads is expressed in terms of lead survival estimates, where “survival” refers to the function of the lead, not the survival of the patient. These survival estimates are intended to illustrate the probability that a lead will survive for a given number of years without a lead-related complication.

The survival estimates are determined from the analysis of the data collected through the SLS. These data are presented graphically and numerically.

The SLS protocol requires regular follow-up reporting on all leads actively followed in the study. Each study center must inform Medtronic whenever a lead complication has occurred or when a patient is no longer participating in the study. Under the study protocol, each lead is assumed to be normally active unless a lead-related complication is confirmed, the lead is abandoned or explanted, or the patient is no longer available for follow-up. The data analyses assume that the patient is still part of the study and there are no lead complications at the time of the report cutoff date unless specifically reported by the center.

Medtronic evaluates center compliance with study protocol through, at a minimum, annual clinical monitoring at each study site. Additionally, study center personnel must be trained in the study procedures prior to participating, and they must adhere to the policies and procedures of their local ethics boards.

Implant times are calculated from the implant date to the earlier of the complication date, out-of-service date (for example, patient leaves the study or the lead is no longer being used), or the cutoff date of the report. If a lead experiences more than one complication, the first is used to calculate survival time; although all complications associated with a lead are reported in PPR tables.

Of the several different statistical methods available for survival analysis, the Standard Actuarial Method, with suspensions assumed distributed across the intervals (Cutler-Ederer Method), is used to determine estimates of lead survival. This method is commonly used by medical researchers and clinicians.

On the following pages, each graph includes a survival curve where events include qualifying lead-related complications. This survival estimate is a good representation of the probability a lead will survive a period of time without a lead-related complication. For example, if a survival probability is 95% after 5 years of service, then the lead has a 5% chance of experiencing a lead-related complication in the first 5 years following implant.

Since the survival estimate can become very imprecise with small effective sample sizes, Medtronic truncates the survival curve when the number of leads entering an interval is less than 50 leads. When the number of leads entering an interval reaches 50, the next data point is added to the survival curve.

Although the report provides tabular data in 1-year intervals, the curves are actually computed and plotted using 3-month intervals.

The data in the tables is rounded to the nearest tenth of one percent. Occasionally, a graph may show 100% survival, but have one or more complications. This occurs because even with the complications, the data rounds to 100%.

The survival curves are statistical estimates. As performance experience accumulates, the estimation improves. Confidence intervals are provided as a way to indicate the degree of certainty of the estimates. Greenwood’s formula is used to calculate corresponding 95% confidence intervals for the standard errors, and the complementary log-log method is used to produce the confidence bounds.

continued

**Medtronic urges all physicians to return explanted products and to notify Medtronic when a product is no longer in use, regardless of reason for explant or removal from use.**

## **Sample Size and How the Population and Population Samples Are Defined**

The population sample from which the survival estimates are derived is comprised of the patients successfully enrolled in the SLS as of the report cutoff date. The number of enrolled implants is listed for each model.

This sample based on SLS enrollments is considered to be representative of the worldwide population, including data from 14 countries on four continents, and therefore the survival estimates shown in this report should be representative of the performance worldwide of these models.

In general, a model or model family will be included in this report when more than 100 leads have been enrolled and no fewer than 50 leads followed for at least 6 months. Models will remain in the report as long as Medtronic estimates at least 500 leads remain active in the United States, based on estimated US implants.

## **Returned Product Analysis Results**

Every lead or lead portion returned to Medtronic receives an analysis. Although the returned product analysis data is not used to generate the survival estimates, the data provides valuable insight into the causes of lead malfunction.

For reporting returned product analysis results, Medtronic CRDM considers a lead as having malfunctioned whenever the analysis shows that any parameter was outside the performance limits established by Medtronic while implanted and in service. To be considered a malfunction for returned product analysis reporting, the lead must have been returned to Medtronic and analyzed.

Lead malfunctions are divided into three categories: Implant Damage, Electrical, and Other. Typical examples of implant damage are stylet perforation, cut or torn insulation, bent or distorted conductors, over-retracted helices, and conductor fracture due to over-torquing. An electrical malfunction is defined as a hardware malfunction resulting in a break in the insulation or a break in the conductor that could affect the electrical performance of the lead. A break in the insulation is defined as a breach allowing entry of body fluids or inappropriate current flow between the conductors, or between the conductor and the body. Examples include cuts, tears, depressions, environmental stress cracking (ESC), and metal ion

oxidation (MIO). A break in the conductor of a lead is defined as the loss of continuity in the metallic components that could interrupt the current flow or voltage. Examples include fractured conductors and defective crimps.

Leads damaged after explant or damaged due to failure to heed warnings or contraindications in the labeling are not considered device malfunctions.

A lead subject to a safety advisory is not considered to have malfunctioned unless it has been returned to Medtronic CRDM and found, through analysis, to actually have performed outside the performance limits established by Medtronic.

For leads designed for either ventricular or atrial use, the numbers listed in the Returned Product Analysis tables include both.

The numbers of malfunctions listed in the Returned Product Analysis tables are the actual numbers confirmed in the returned product analysis from the United States. The numbers of complications listed in the complications tables are the actual numbers observed in the SLS centers around the world.

## **Estimated Number of Implanted and Active Leads in the United States**

In addition to providing the number of leads enrolled in the SLS, this report also provides the number of leads registered as implanted and the number remaining active based on the status recorded in the Medtronic Device Registration System.

Some lead models do not have a survival curve presented in this report. These lead models do not have a survival curve because they have insufficient sample size in the System Longevity Study. Returned Product Analysis results for these models are included here for reference and comparison.

# Left-Heart Leads

## 2187 Attain


### Product Characteristics

US Market Release	Aug-01	Serial Number Prefix	LEY	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	12,000	Type and/or Fixation	Transvenous, Left Ventricular Cardiac Vein, Distal Continuous Curve	Implant Damage	7
Estimated Active US Implants	4,300	Polarity	Unipolar	Electrical Malfunction	0
Advisories	None	Steroid	No	Other	16

### System Longevity Study Results

### Qualifying Complications 1 Total

Number of Leads Enrolled in Study	134	Failure to Capture	1
Cumulative Months of Follow-Up	5,978		
Number of Leads Active in Study	38		



## 2188 Attain


### Product Characteristics

US Market Release	Aug-01	Serial Number Prefix	LEB	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	1,800	Type and/or Fixation	Transvenous, Coronary Sinus/ Cardiac Vein, Canted	Implant Damage	1
Estimated Active US Implants	400	Polarity	Bipolar	Electrical Malfunction	1
Advisories	None	Steroid	No	Other	0

### System Longevity Study Results

### Qualifying Complications 1 Total

Number of Leads Enrolled in Study	14	Extra Cardiac Stimulation	1
Cumulative Months of Follow-Up	383		
Number of Leads Active in Study	1		



Leads

# Left-Heart Leads continued

## 4193 Attain


### Product Characteristics

US Market Release	May-02	Serial Number Prefix	BAA	<u>US Returned Product Analysis</u>	
Registered US Implants	98,100	Type and/or Fixation	Transvenous, Left Ventricular Cardiac Vein, Distal Double Curve	Implant Damage	65
Estimated Active US Implants	54,100	Polarity	Unipolar	Electrical Malfunction	23
Advisories	None	Steroid	Yes	Other	64

### System Longevity Study Results

### Qualifying Complications 34 Total

Number of Leads Enrolled in Study	673	Conductor Fracture	1	Lead Dislodgement	13
Cumulative Months of Follow-Up	22,161	Extra Cardiac Stimulation	6	Unspecified Clinical Failure	3
Number of Leads Active in Study	259	Failure to Capture	11		



## 4194 Attain


### Product Characteristics

US Market Release	Aug-04	Serial Number Prefix	LFG	<u>US Returned Product Analysis</u>	
Registered US Implants	75,000	Type and/or Fixation	Transvenous, Left Ventricular Cardiac Vein, Distal Double Curve	Implant Damage	96
Estimated Active US Implants	57,700	Polarity	Bipolar	Electrical Malfunction	14
Advisories	None	Steroid	Yes	Other	7

### System Longevity Study Results

### Qualifying Complications 7 Total

Number of Leads Enrolled in Study	517	Failure to Capture	1
Cumulative Months of Follow-Up	10,982	Lead Dislodgement	6
Number of Leads Active in Study	401		



## 4195 Attain


### Product Characteristics

				<u>US Returned Product Analysis</u>	
US Market Release	Aug-08	Serial Number Prefix	AAD		
Registered US Implants	2,420	Type and/or Fixation	Transvenous, Left Ventricular Cardiac Vein, Deployable Lobe Fixation	Implant Damage	11
Estimated Active US Implants	2,400	Polarity	Unipolar	Electrical Malfunction	2
Advisories	None	Steroid	Yes	Other	8

### System Longevity Study Results

Qualifying Complications 0 Total

Number of Leads Enrolled in Study	24
Cumulative Months of Follow-Up	373
Number of Leads Active in Study	23



## Lead Survival Summary (95% Confidence Interval)

Model Number	Family	US Market Release	Leads Enrolled	Leads Active in Study	Qualifying Complications	Cumulative Months of Follow-Up in Study	Device Survival Probability (%)									
							Years After Implant									
							1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	9 yr	10 yr
2187	Attain	Aug-01	134	38	1	5,978	99.1 +0.8/-5.2	99.1 +0.8/-5.2	99.1 +0.8/-5.2	99.1 +0.8/-5.2	99.1 +0.8/-5.2 at 54 mo					
2188	Attain	Aug-01	14	1	1	383	Survival estimate not available due to insufficient sample size									
4193	Attain	May-02	673	259	34	22,161	95.8 +1.3/-2.0	94.7 +1.6/-2.2	93.9 +1.7/-2.5	93.9 +1.7/-2.5	93.9 +1.7/-2.5	93.9 +1.7/-2.5 at 69 mo				
4194	Attain	Aug-04	517	401	7	10,982	99.1 +0.6/-1.6	98.0 +1.1/-2.2	98.0 +1.1/-2.2	98.0 +1.1/-2.2 at 42 mo						
4195	Attain	Aug-08	24	23	0	373	Survival estimate not available due to insufficient sample size									

Source: System Longevity Study  
Data as of January 31, 2009

## US Returned Product Analysis Summary

Model Number	Family	US Market Release	Estimated US Implants	Estimated US Active	Implant Damage	Electrical Malfunction	Other
2187	Attain	Aug-01	12,000	4,300	7	0	16
2188	Attain	Aug-01	1,800	400	1	1	0
4193	Attain	May-02	98,100	54,100	64	19	65
4194	Attain	Aug-04	75,000	57,700	93	6	7
4195	Attain	Jun-08	2,420	2,400	11	2	8

Source: Returned Product Analysis  
Data as of January 31, 2009

## Reference Chart

Model Number	Family	Type	Insulation	Conductor Material	Tip Electrode	Connector Type
2187	Attain	Transvenous Cardiac Vein Preformed Body	Polyurethane (55D)	MP35N	Platinized Platinum	IS-1 UNI
2188	Attain	Transvenous Cardiac Vein Preformed Body	Polyurethane (55D)	MP35N	Platinized Platinum	IS-1 BI
4193	Attain	Transvenous Cardiac Vein Preformed Body	Polyurethane (55D)	MP35N	Platinized Platinum	IS-1 UNI
4194	Attain	Transvenous Cardiac Vein Preformed Body	Polyurethane (55D)/ Silicone (4719)	MP35N	Platinum Alloy	IS-1 BI
4195	Attain	Transvenous Cardiac Vein Deployable Lobes	Polyurethane (55D)	MP35N	Platinum Alloy	IS-1 Uni

# Defibrillation Leads

## 6721, 6921 Epicardial Patch


### Product Characteristics

US Market Release	Feb-93	Serial Number Prefix	TBH, TBG, TBB, TAD, TAC, or TAB	<u>US Returned Product Analysis</u>	
Registered US Implants	8,300	Type and/or Fixation	Epicardial Defib Patch, Suture	Implant Damage	5
Estimated Active US Implants	1,300	Polarity	Defib Electrode only	Electrical Malfunction	80
Advisories	None	Steroid	No	Other	0

### System Longevity Study Results

### Qualifying Complications 28 Total

Number of Leads Enrolled in Study	407	Conductor Fracture	20	Insulation (not further defined)	3
Cumulative Months of Follow-Up	18,283	Failure to Capture	2		
Number of Leads Active in Study	12	Impedance Out of Range	3		



Leads

## 6930 Sprint Fidelis


### Product Characteristics

US Market Release	Sep-04	Serial Number Prefix	LFK	<u>US Returned Product Analysis</u>	
Registered US Implants	400	Type and/or Fixation	Transvenous, Vent, Defib and Pace/ Sense, Tines	Implant Damage	0
Estimated Active US Implants	300	Polarity	True Bipolar/One Coil	Electrical Malfunction	2
Advisories	1	Steroid	Yes	Other	0
See page 151 – 2007 Potential Conductor Wire Fracture					

### System Longevity Study Results

### Qualifying Complications 0 Total

Number of Leads Enrolled in Study	4
Cumulative Months of Follow-Up	83
Number of Leads Active in Study	4





# Defibrillation Leads continued

## 6931 Sprint Fidelis


### Product Characteristics

US Market Release	Sep-04	Serial Number Prefix	LFL	<u>US Returned Product Analysis</u>	
Registered US Implants	8,100	Type and/or Fixation	Transvenous, Vent, Defib and Pace/ Sense, Screw-in	Implant Damage	28
Estimated Active US Implants	6,000	Polarity	True Bipolar/One Coil	Electrical Malfunction	102
Advisories	1	Steroid	Yes	Other	0
See page 151 – 2007 Potential Conductor Wire Fracture					

### System Longevity Study Results

### Qualifying Complications 11 Total

Number of Leads Enrolled in Study	292	Failure to Capture	3	Lead Dislodgement	2
Cumulative Months of Follow-Up	6,804	Failure to Sense	1	Oversensing	3
Number of Leads Active in Study	228	Impedance Out of Range	2		



## 6932 Sprint


### Product Characteristics

US Market Release	Aug-96	Serial Number Prefix	TCA	<u>US Returned Product Analysis</u>	
Registered US Implants	15,000	Type and/or Fixation	Transvenous, Vent, Defib and Pace/ Sense, Tines	Implant Damage	16
Estimated Active US Implants	5,700	Polarity	True Bipolar/One Coil	Electrical Malfunction	39
Advisories	None	Steroid	Yes	Other	7

### System Longevity Study Results

### Qualifying Complications 8 Total

Number of Leads Enrolled in Study	410	Extra Cardiac Stimulation	1
Cumulative Months of Follow-Up	20,997	Failure to Capture	2
Number of Leads Active in Study	72	Failure to Sense	2
		Oversensing	3



# Defibrillation Leads continued

## 6933, 6937, 6937A, 6963 SVC/CS Product Characteristics


US Market Release	Dec-93	Serial Number Prefix	TAT, TBU, or TAF	<u>US Returned Product Analysis</u>	
Registered US Implants	15,700	Type and/or Fixation	Transvenous CS or SVC Defib	Implant Damage	31
Estimated Active US Implants	2,700	Polarity	One Defib Coil	Electrical Malfunction	200
Advisories	None	Steroid	No	Other	13

### System Longevity Study Results

Number of Leads Enrolled in Study	966
Cumulative Months of Follow-Up	47,888
Number of Leads Active in Study	32

### Qualifying Complications 25 Total

Conductor Fracture	15	Lead Dislodgement	1
Failure to Capture	1	Unspecified Clinicial Failure	4
Impedance Out of Range	2		
Insulation (not further defined)	2		




## 6935 Sprint Quattro Secure Product Characteristics

US Market Release	Nov-08	Serial Number Prefix	TAU	<u>US Returned Product Analysis</u>	
Registered US Implants	590	Type and/or Fixation	Transvenous, Vent, Defib and Pace/ Sense, Screw-in	Implant Damage	0
Estimated Active US Implants	580	Polarity	True Bipolar/One Coil	Electrical Malfunction	0
Advisories	None	Steroid	Yes	Other	0

### System Longevity Study Results

Number of Leads Enrolled in Study	2
Cumulative Months of Follow-Up	0
Number of Leads Active in Study	2

### Qualifying Complications 0 Total



# Defibrillation Leads continued

## 6936, 6966 Transvene

### Product Characteristics

US Market Release	Dec-93	Serial Number Prefix	TAV or TAL
Registered US Implants	23,700	Type and/or Fixation	Transvenous, Vent, Defib and Pace/Sense, Screw-in
Estimated Active US Implants	3,200	Polarity	True Bipolar/One Coil
Advisories	None	Steroid	No

### US Returned Product Analysis


Implant Damage	90
Electrical Malfunction	487
Other	19

### System Longevity Study Results

Number of Leads Enrolled in Study	1,350
Cumulative Months of Follow-Up	67,840
Number of Leads Active in Study	40

### Qualifying Complications 150 Total

Conductor Fracture	18	Impedance Out of Range	5
Extra Cardiac Stimulation	2	Insulation (not further defined)	14
Failure to Capture	9	Oversensing	93
Failure to Sense	4	Unspecified Clinical Failure	5



## 6939, 6999 Sub-Q Patch

### Product Characteristics

US Market Release	Dec-93	Serial Number Prefix	TBA or TAP
Registered US Implants	3,600	Type and/or Fixation	Subcutaneous Defib Patch, Suture
Estimated Active US Implants	300	Polarity	Defib Electrode Only
Advisories	None	Steroid	No

### US Returned Product Analysis


Implant Damage	4
Electrical Malfunction	33
Other	1

### System Longevity Study Results

Number of Leads Enrolled in Study	384
Cumulative Months of Follow-Up	17,726
Number of Leads Active in Study	5

### Qualifying Complications 20 Total

Conductor Fracture	10	Unspecified Clinical Failure	2
Failure to Capture	2		
Insulation (not further defined)	6		



# Defibrillation Leads continued

## 6942 Sprint


### Product Characteristics

US Market Release	Jul-97	Serial Number Prefix	TCB	<u>US Returned Product Analysis</u>	
Registered US Implants	17,700	Type and/or Fixation	Transvenous, Vent, Defib and Pace/Sense, Tines	Implant Damage	31
Estimated Active US Implants	7,400	Polarity	Integrated Bipolar/Two Coils	Electrical Malfunction	37
Advisories	None	Steroid	Yes	Other	5

### System Longevity Study Results

### Qualifying Complications 7 Total

Number of Leads Enrolled in Study	351	Conductor Fracture	1	Oversensing	3
Cumulative Months of Follow-Up	15,650	Failure to Sense	1	Unspecified Clinical Failure	1
Number of Leads Active in Study	57	Lead Dislodgement	1		



Leads

## 6943 Sprint


### Product Characteristics

US Market Release	Oct-97	Serial Number Prefix	TCE	<u>US Returned Product Analysis</u>	
Registered US Implants	20,800	Type and/or Fixation	Transvenous, Vent, Defib and Pace/Sense, Screw-in	Implant Damage	51
Estimated Active US Implants	8,800	Polarity	True Bipolar/One Coil	Electrical Malfunction	71
Advisories	None	Steroid	Yes	Other	8

### System Longevity Study Results

### Qualifying Complications 67 Total

Number of Leads Enrolled in Study	1,311	Conductor Fracture	15	Insulation (not further defined)	1
Cumulative Months of Follow-Up	65,157	Failure to Capture	7	Lead Dislodgement	1
Number of Leads Active in Study	381	Failure to Sense	5	Oversensing	31
		Impedance Out of Range	4	Unspecified Clinical Failure	3



# Defibrillation Leads continued

## 6944 Sprint Quattro


### Product Characteristics

US Market Release	Dec-00	Serial Number Prefix	TDC	<u>US Returned Product Analysis</u>	
Registered US Implants	31,900	Type and/or Fixation	Transvenous, Vent, Defib and Pace/Sense, Tines	Implant Damage	31
Estimated Active US Implants	19,200	Polarity	True Bipolar/Two Coils	Electrical Malfunction	35
Advisories	None	Steroid	Yes	Other	8

### System Longevity Study Results

### Qualifying Complications 3 Total

Number of Leads Enrolled in Study	188	Oversensing	2
Cumulative Months of Follow-Up	7,482	Unspecified Clinical Failure	1
Number of Leads Active in Study	75		



## 6945 Sprint


### Product Characteristics

US Market Release	Sep-97	Serial Number Prefix	TDA	<u>US Returned Product Analysis</u>	
Registered US Implants	42,800	Type and/or Fixation	Transvenous, Vent, Defib and Pace/Sense, Screw-in	Implant Damage	198
Estimated Active US Implants	18,100	Polarity	Integrated Bipolar/Two Coils	Electrical Malfunction	118
Advisories	None	Steroid	Yes	Other	11

### System Longevity Study Results

### Qualifying Complications 27 Total

Number of Leads Enrolled in Study	1,157	Conductor Fracture	3	Impedance Out of Range	5
Cumulative Months of Follow-Up	54,899	Extra Cardiac Stimulation	1	Oversensing	12
Number of Leads Active in Study	229	Failure to Capture	1	Unspecified Clinical Failure	1
		Failure to Sense	4		



# Defibrillation Leads continued

## 6947 Sprint Quattro Secure

### Product Characteristics


US Market Release	Nov-01	Serial Number Prefix	TDG	<u>US Returned Product Analysis</u>	
Registered US Implants	206,100	Type and/or Fixation	Transvenous, Vent, Defib and Pace/Sense, Screw-in	Implant Damage	476
Estimated Active US Implants	151,200	Polarity	True Bipolar/Two Coils	Electrical Malfunction	137
Advisories	None	Steroid	Yes	Other	25

### System Longevity Study Results

Number of Leads Enrolled in Study	1,397
Cumulative Months of Follow-Up	54,528
Number of Leads Active in Study	597

### Qualifying Complications 23 Total

Conductor Fracture	3	Lead Dislodgement	3
Failure to Sense	2	Oversensing	7
Impedance Out of Range	4	Unspecified Clinical Failure	2
Insulation (not further defined)	2		



## 6948 Sprint Fidelis


### Product Characteristics

US Market Release	Sep-04	Serial Number Prefix	LFH	<u>US Returned Product Analysis</u>	
Registered US Implants	10,400	Type and/or Fixation	Transvenous, Vent, Defib and Pace/Sense, Tines	Implant Damage	9
Estimated Active US Implants	7,900	Polarity	True Bipolar/Two Coils	Electrical Malfunction	17
Advisories	1	Steroid	Yes	Other	4
<a href="#">See page 151</a> – 2007 Potential Conductor Wire Fracture					

### System Longevity Study Results

Number of Leads Enrolled in Study	30
Cumulative Months of Follow-Up	797
Number of Leads Active in Study	26

### Qualifying Complications 0 Total



# Defibrillation Leads continued

## 6949 Sprint Fidelis


### Product Characteristics

US Market Release	Sep-04	Serial Number Prefix	LFJ	<u>US Returned Product Analysis</u>	
Registered US Implants	186,700	Type and/or Fixation	Transvenous, Vent, Defib and Pace/ Sense, Screw-in	Implant Damage	461
Estimated Active US Implants	135,900	Polarity	True Bipolar/Two Coils	Electrical Malfunction	1,435
Advisories	1	Steroid	Yes	Other	45
See page 151 – 2007 Potential Conductor Wire Fracture					

### System Longevity Study Results

### Qualifying Complications 30 Total

Number of Leads Enrolled in Study	789	Conductor Fracture	9	Insulation (not further defined)	1
Cumulative Months of Follow-Up	21,439	Failure to Capture	2	Lead Dislodgement	1
Number of Leads Active in Study	576	Failure to Sense	2	Oversensing	11
		Impedance Out of Range	4		



## 6996 Sub-Q Lead


### Product Characteristics

US Market Release	Jun-01	Serial Number Prefix	TCR	<u>US Returned Product Analysis</u>	
Registered US Implants	2,300	Type and/or Fixation	Subcutaneous Defib Coil, Suture	Implant Damage	1
Estimated Active US Implants	1,600	Polarity	One Defib Coil	Electrical Malfunction	3
Advisories	None	Steroid	No	Other	0

### System Longevity Study Results

### Qualifying Complications 0 Total

Number of Leads Enrolled in Study	17
Cumulative Months of Follow-Up	445
Number of Leads Active in Study	14



# Defibrillation Leads continued

## Lead Survival Summary (95% Confidence Interval)

Model Number	Family	US Market Release	Leads Enrolled	Leads Active in Study	Qualifying Complications	Cumulative Months of Follow-Up in Study	Device Survival Probability (%)																
							Years After Implant																
							1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr	12 yr	14 yr	16 yr	18 yr	20 yr			
<b>6721, 6921</b>	<b>Epicardial Patch</b>	Feb-93	407	12	28	18,283	100.0	98.6 +0.9/-2.3	96.0 +1.9/-3.3	95.0 +2.1/-3.8	91.6 +3.2/-4.9	84.4 +4.9/-7.0	82.1 +5.5/-7.7	82.1 +5.5/-7.7									
<b>6930</b>	<b>Sprint Fidelis</b>	Sep-04	4	4	—	83	Survival estimate not available due to insufficient sample size																
Advisories: See page 151 – 2007 Potential Conductor Wire Fracture																							
<b>6931</b>	<b>Sprint Fidelis</b>	Sep-04	292	228	11	6,804	98.2 +1.0/-2.5	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo	96.4 +1.7/-3.1 at 27 mo		
Advisories: See page 151 – 2007 Potential Conductor Wire Fracture																							
<b>6932</b>	<b>Sprint</b>	Aug-96	410	72	8	20,997	99.4 +0.4/-1.8	98.4 +0.9/-2.3	98.4 +0.9/-2.3	98.4 +0.9/-2.3	97.8 +1.2/-2.9	97.8 +1.2/-2.9	96.7 +1.9/-4.1 at 114 mo	96.7 +1.9/-4.1 at 114 mo									
<b>6933, 6937, 6937A, 6963</b>	<b>SVC/CS</b>	Dec-93	966	32	25	47,888	99.6 +0.3/-0.8	99.2 +0.4/-1.0	99.2 +0.4/-1.0	98.6 +0.7/-1.4	97.0 +1.2/-2.0	96.4 +1.4/-2.2	95.2 +1.8/-2.8	94.6 +2.0/-3.1	92.2 +2.9/-4.6 at 141 mo								
<b>6935</b>	<b>Sprint Quattro Secure</b>	Nov-08	2	2	0	0	Survival estimate not available due to insufficient sample size																
<b>6936, 6966</b>	<b>Transvene</b>	Dec-93	1,350	40	150	67,840	99.2 +0.4/-0.7	98.4 +0.6/-1.0	97.0 +1.0/-1.3	96.0 +1.1/-1.6	92.7 +1.7/-2.3	87.4 +2.4/-3.1	79.5 +3.4/-4.0	75.5 +3.9/-4.5	62.5 +2.9/-6.2 at 153 mo	60.3 +6.1/-6.7 at 153 mo							
<b>6939, 6999</b>	<b>Sub-Q Patch</b>	Dec-93	384	5	20	17,726	99.0 +0.7/-1.9	98.7 +0.8/-2.3	98.2 +1.0/-2.6	98.2 +1.0/-2.6	94.1 +2.7/-4.7	91.1 +3.6/-5.7	87.7 +4.7/-7.1	86.1 +5.2/-7.9 at 99 mo									
<b>6942</b>	<b>Sprint</b>	Jul-97	351	57	7	15,650	98.9 +0.8/-2.2	98.9 +0.8/-2.2	97.8 +1.3/-3.1	97.2 +1.5/-3.6	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	96.3 +2.0/-4.1	
<b>6943</b>	<b>Sprint</b>	Oct-97	1,311	381	67	65,157	98.8 +0.5/-0.8	98.0 +0.6/-1.1	96.8 +1.0/-1.2	95.8 +1.1/-1.5	93.6 +1.5/-2.0	91.8 +1.9/-2.3	91.2 +2.0/-2.5	90.8 +2.1/-2.6									
<b>6944</b>	<b>Sprint Quattro</b>	Dec-00	188	75	3	7,482	100.0	100.0	98.9 +0.9/-6.7	97.6 +1.8/-7.0	95.9 +2.8/-8.1 at 69 mo	95.9 +2.8/-8.1 at 69 mo											
<b>6945</b>	<b>Sprint</b>	Sep-97	1,157	229	27	54,899	99.6 +0.2/-0.7	99.1 +0.5/-0.9	98.8 +0.6/-1.0	98.1 +0.8/-1.2	97.0 +1.1/-1.8	95.9 +1.4/-2.2	95.5 +1.6/-2.3	94.0 +2.1/-3.1 at 111 mo									
<b>6947</b>	<b>Sprint Quattro Secure</b>	Nov-01	1,397	597	23	54,528	99.2 +0.4/-0.7	98.9 +0.5/-0.8	98.5 +0.6/-0.9	98.1 +0.7/-1.2	97.4 +1.0/-1.4	97.0 +1.1/-1.8 at 78 mo	97.0 +1.1/-1.8 at 78 mo										
<b>6948</b>	<b>Sprint Fidelis</b>	Sep-04	30	26	0	797	Survival estimate not available due to insufficient sample size																
Advisories: See page 151 – 2007 Potential Conductor Wire Fracture																							
<b>6949</b>	<b>Sprint Fidelis</b>	Sep-04	789	576	30	21,439	98.6 +0.7/-1.2	96.6 +1.2/-1.7	94.3 +1.8/-2.7	93.6 +2.1/-3.1													
Advisories: See page 151 – 2007 Potential Conductor Wire Fracture																							
<b>6996</b>	<b>Sub-Q Lead</b>	Jun-01	17	14	0	445	Survival estimate not available due to insufficient sample size																

## Leads



## US Returned Product Analysis Summary

Model Number	Family	US Market Release	Estimated US Implants	Estimated US Active	Implant Damage	Electrical Malfunction	Other
6721, 6921	Epicardial Patch	Feb-93	8,300	1,300	5	80	0
6930	Sprint Fidelis	Sep-04	400	300	0	2	0
6931	Sprint Fidelis	Sep-04	8,100	6,000	28	102	0
6932	Sprint	Aug-96	15,000	5,700	16	39	7
6933, 6937, 6937A, 6963	SVC/CS	Dec-93	15,700	2,700	31	200	13
6935	Sprint Quattro Secure	Nov-08	590	580	0	0	0
6936, 6966	Transvene	Dec-93	23,700	3,200	90	487	19
6939, 6999	Sub-Q Patch	Dec-93	3,600	300	4	33	1
6942	Sprint	Jul-97	17,700	7,400	31	37	5
6943	Sprint	Oct-97	20,800	8,800	51	71	8
6944	Sprint Quattro	Dec-00	31,900	19,200	31	35	8
6945	Sprint	Sep-97	42,800	18,100	198	118	11
6947	Sprint Quattro Secure	Nov-01	206,100	151,200	476	137	25
6948	Sprint Fidelis	Sep-04	10,400	7,900	9	17	4
6949	Sprint Fidelis	Sep-04	186,700	135,900	461	1,435	45
6996	Sub-Q Lead	Jun-01	2,300	1,600	1	3	0

## Reference Chart

Model Number	Family	Type	Pin Configuration		Lead Body Diameter	Insulation, Lead Body	Fixation, Steroid
			Pace/Sense	High Voltage			
6721	Epicardial Patch	Epi Patch	—	DF-1	S, M, L	Silicone, Single Lumen	Suture
6921	Epicardial Patch	Epi Patch	—	6.5 mm	S, M, L	Silicone, Single Lumen	Suture
6930	Sprint Fidelis	Endo RV True Bipolar Sensing	IS-1	DF-1	6.6 Fr	Silicone with Polyurethane Overlay, Multilumen	Passive, Steroid
6931	Sprint Fidelis	Endo RV True Bipolar Sensing	IS-1	DF-1	6.6 Fr	Silicone with Polyurethane Overlay, Multilumen	Active, Steroid
6932	Sprint	Endo RV True Bipolar Sensing	IS-1	DF-1	7.8 Fr	Silicone, Multilumen	Passive, Steroid
6933	SVC/CS	Endo SVC/CS Coil	—	DF-1	7 Fr	Silicone, Single Lumen	Passive
6934S	Transvene	Endo RV True Bipolar Sensing	IS-1	DF-1	12 Fr	Silicone, Coaxial	Passive, Steroid
6935	Sprint Quattro Secure	Endo RV True Bipolar Sensing	IS-1	DF-1	8.2 Fr	Silicone with Polyurethane Overlay, Multilumen	Active, Steroid
6936	Transvene	Endo RV True Bipolar Sensing	IS-1	DF-1	10 Fr	Polyurethane, Coaxial	Active
6937	SVC/CS	Endo SVC Coil	—	DF-1	5.5 Fr	Silicone, Single Lumen	Passive
6937A	SVC/CS	Endo SVC Coil	—	DF-1	7.5 Fr	Silicone with Polyurethane Overlay, Single Lumen	Passive
6939	Sub-Q Patch	SQ Patch	—	DF-1	One Size	Silicone, Single Lumen	Suture
6942	Sprint	Endo RV/SVC Integrated Bipolar Sensing	IS-1	2 DF-1	7.8 Fr	Silicone, Multilumen	Passive, Steroid
6943	Sprint	Endo RV True Bipolar Sensing	IS-1	DF-1	7.8 Fr	Silicone, Multilumen	Active, Steroid
6944	Sprint Quattro	Endo RV/SVC True Bipolar Sensing	IS-1	2 DF-1	8.2 Fr	Silicone with Polyurethane Overlay, Multilumen	Passive, Steroid
6945	Sprint	Endo RV/SVC Integrated Bipolar Sensing	IS-1	2 DF-1	7.8 Fr	Silicone, Multilumen	Active, Steroid
6947	Sprint Quattro Secure	Endo RV/SVC True Bipolar Sensing	IS-1	2 DF-1	8.2 Fr	Silicone with Polyurethane Overlay, Multilumen	Active, Steroid
6948	Sprint Fidelis	Endo RV/SVC True Bipolar Sensing	IS-1	2 DF-1	6.6 Fr	Silicone with Polyurethane Overlay, Multilumen	Passive, Steroid
6949	Sprint Fidelis	Endo RV/SVC True Bipolar Sensing	IS-1	2 DF-1	6.6 Fr	Silicone with Polyurethane Overlay, Multilumen	Active, Steroid
6963	SVC/CS	Endo SVC/CS Coil	—	6.5 mm	7 Fr	Silicone, Single Lumen	Passive
6966	Transvene	Endo RV True Bipolar Sensing	3.2 mm L.P.	6.5 mm	10 Fr	Polyurethane, Coaxial	Active
6996	Sub-Q Lead	SQ Coil	—	DF-1	7.5 Fr	Silicone, Single Lumen	Passive
6999	Sub-Q Patch	SQ Patch	—	6.5 mm	One Size	Silicone, Single Lumen	Suture

# Pacing Leads

## 3830 SelectSecure

### Product Characteristics


US Market Release	Aug-05	Serial Number Prefix	LFF	<u>US Returned Product Analysis</u>	
Registered US Implants	12,200	Type and/or Fixation	Transvenous, V or A, Screw-in	Implant Damage	21
Estimated Active US Implants	10,600	Polarity	Bipolar	Electrical Malfunction	5
Advisories	None	Steroid	Yes	Other	1

## Atrial Placement

### System Longevity Study Results

### Qualifying Complications 1 Total

Number of Leads Enrolled in Study	166	Failure to Sense	1
Cumulative Months of Follow-Up	5,129		
Number of Leads Active in Study	95		




## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 0 Total

Number of Leads Enrolled in Study	148
Cumulative Months of Follow-Up	4,983
Number of Leads Active in Study	82



# Pacing Leads continued

## 4003, 4003M CapSure

### Product Characteristics


US Market Release	Jul-86	Serial Number Prefix	IH or LAX	<u>US Returned Product Analysis</u>	
Registered US Implants	38,000	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	24
Estimated Active US Implants	5,300	Polarity	Unipolar	Electrical Malfunction	60
Advisories	None	Steroid	Yes	Other	2

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications

Number of Leads Enrolled in Study	711	Conductor Fracture	1	Failure to Sense	1
Cumulative Months of Follow-Up	47,264	Extra Cardiac Stimulation	2	Oversensing	2
Number of Leads Active in Study	146	Failure to Capture	6		



Leads

# Pacing Leads continued

## 4004, 4004M CapSure

### Product Characteristics


US Market Release	Feb-89	Serial Number Prefix	PS or LAV	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	72,600	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	56
Estimated Active US Implants	5,800	Polarity	Bipolar	Electrical Malfunction	688
Advisories	1	Steroid	Yes	Other	19
See page 157 – 1993 Lead Survival Below Expectations					

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 277 Total

Number of Leads Enrolled in Study	1,640	Conductor Fracture	7	Insulation (ESC)	4
Cumulative Months of Follow-Up	71,653	Electrical Abandonment	1	Insulation (MIO)	4
Number of Leads Active in Study	4	Extra Cardiac Stimulation	2	Insulation (not further defined)	7
		Failure to Capture	131	Medical Judgment	1
		Failure to Sense	62	Oversensing	25
		Impedance Out of Range	32	Unspecified Clinical Failure	1



	1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	9 yr	10 yr <sup>at 129 mo</sup>
%	99.8	99.3	96.3	87.4	77.1	69.4	64.1	58.1	53.1	51.5
#	1,192	1,020	824	630	453	314	231	161	118	78
Effective Sample Size										

## 4011 Target Tip

### Product Characteristics


US Market Release	Nov-82	Serial Number Prefix	IB	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	58,400	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	29
Estimated Active US Implants	6,200	Polarity	Unipolar	Electrical Malfunction	152
Advisories	None	Steroid	No	Other	5

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 25 Total

Number of Leads Enrolled in Study	851	Conductor Fracture	1	Insulation (not further defined)	10
Cumulative Months of Follow-Up	54,409	Extra Cardiac Stimulation	4	Oversensing	1
Number of Leads Active in Study	2	Failure to Capture	9		



	1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	9 yr	10 yr	11 yr	12 yr	13 yr	14 yr	15 yr <sup>at 183 mo</sup>
%	99.4	99.2	99.1	98.8	97.6	96.4	96.0	96.0	96.0	95.0	93.6	92.8	91.9	91.9	91.9
#	626	556	475	414	353	299	250	219	189	165	134	109	81	71	55
Effective Sample Size															

# Pacing Leads continued

## 4012 Target Tip

### Product Characteristics


US Market Release	Jul-83	Serial Number Prefix	HQ	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	93,700	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	50
Estimated Active US Implants	6,500	Polarity	Bipolar	Electrical Malfunction	827
Advisories	1	Steroid	No	Other	34
See page 158 – 1991 Lead Survival Below Expectations					

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 316 Total

Number of Leads Enrolled in Study	2,543	Conductor Fracture	6	Insulation (ESC)	9
Cumulative Months of Follow-Up	151,162	Extra Cardiac Stimulation	3	Insulation (MIO)	4
Number of Leads Active in Study	10	Failure to Capture	126	Insulation (not further defined)	16
		Failure to Sense	77	Medical Judgment	1
		Impedance Out of Range	26	Oversensing	48



Years After Implant

	1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	9 yr	10 yr	11 yr	12 yr	13 yr	14 yr	15 yr	at 189 mo
%	99.6	99.1	98.4	95.9	92.6	88.1	83.9	77.8	73.6	69.8	66.9	65.7	65.0	62.9	62.2	62.2
#	1,935	1,714	1,528	1,310	1,084	888	698	522	400	307	243	200	144	98	69	51

Effective Sample Size

Leads

## 4023 CapSure SP

### Product Characteristics


US Market Release	Aug-91	Serial Number Prefix	LAK	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	41,200	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	47
Estimated Active US Implants	9,600	Polarity	Unipolar	Electrical Malfunction	21
Advisories	None	Steroid	Yes	Other	6

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 20 Total

Number of Leads Enrolled in Study	1,158	Extra Cardiac Stimulation	1	Insulation (not further defined)	1
Cumulative Months of Follow-Up	65,853	Failure to Capture	15	Lead Dislodgement	2
Number of Leads Active in Study	341	Impedance Out of Range	1		



Years After Implant

	1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	9 yr	10 yr	11 yr	12 yr	at 147 mo
%	99.9	99.3	98.8	98.6	98.6	98.2	97.2	96.8	96.8	96.3	95.0	95.0	95.0
#	886	765	681	602	513	435	329	238	171	111	65	53	45

Effective Sample Size

# Pacing Leads continued

## 4024 CapSure SP

### Product Characteristics

US Market Release	Oct-91	Serial Number Prefix	LAJ	<u>US Returned Product Analysis</u>	
Registered US Implants	222,100	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	264
Estimated Active US Implants	56,100	Polarity	Bipolar	Electrical Malfunction	135
Advisories	None	Steroid	Yes	Other	34


### Ventricular Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	1,214
Cumulative Months of Follow-Up	52,403
Number of Leads Active in Study	23

#### Qualifying Complications 3 Total

Failure to Capture	3
--------------------	---



## 4033 CapSure Z

### Product Characteristics

US Market Release	Not US released	Serial Number Prefix	LCA	<u>US Returned Product Analysis</u>	
Registered US Implants	NA	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	2
Estimated Active US Implants	NA	Polarity	Unipolar	Electrical Malfunction	0
Advisories	None	Steroid	Yes	Other	0


### Ventricular Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	539
Cumulative Months of Follow-Up	28,545
Number of Leads Active in Study	39

#### Qualifying Complications 10 Total

Conductor Fracture	1
Failure to Capture	8
Impedance Out of Range	1



# Pacing Leads continued

## 4057, 4057M Screw-In

### Product Characteristics


US Market Release		Serial Number Prefix		XQ or LAN		US Returned Product Analysis	
Registered US Implants	10,100	Type and/or Fixation	Transvenous, V or A, Screw-in		Implant Damage	39	
Estimated Active US Implants	1,800	Polarity	Unipolar		Electrical Malfunction	6	
Advisories	None	Steroid	No		Other	4	

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 7 Total

Number of Leads Enrolled in Study	259	Conductor Fracture	2	Failure to Sense	1
Cumulative Months of Follow-Up	15,351	Extra Cardiac Stimulation	2		
Number of Leads Active in Study	9	Failure to Capture	2		





# Pacing Leads continued

## 4058, 4058M Screw-In

### Product Characteristics


US Market Release	Jan-89	Serial Number Prefix	ZY or LAW	<u>US Returned Product Analysis</u>	
Registered US Implants	101,900	Type and/or Fixation	Transvenous, V or A, Screw-in	Implant Damage	388
Estimated Active US Implants	20,800	Polarity	Bipolar	Electrical Malfunction	261
Advisories	None	Steroid	No	Other	23

## Atrial Placement

### System Longevity Study Results

### Qualifying Complications 33 Total

Number of Leads Enrolled in Study	2,364	Extra Cardiac Stimulation	1	Insulation (not further defined)	1
Cumulative Months of Follow-Up	131,441	Failure to Capture	15	Lead Dislodgement	3
Number of Leads Active in Study	44	Failure to Sense	7	Oversensing	1
		Impedance Out of Range	5		




## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 53 Total

Number of Leads Enrolled in Study	1,690	Conductor Fracture	5	Impedance Out of Range	7
Cumulative Months of Follow-Up	77,493	Extra Cardiac Stimulation	3	Insulation (not further defined)	4
Number of Leads Active in Study	49	Failure to Capture	22	Lead Dislodgement	1
		Failure to Sense	10	Oversensing	1



# Pacing Leads continued

## 4067 CapSureFix

### Product Characteristics


				US Returned Product Analysis	
US Market Release	Jan-97	Serial Number Prefix	LCV	Implant Damage	3
Registered US Implants	1,000	Type and/or Fixation	Transvenous, V or A, Screw-in	Electrical Malfunction	1
Estimated Active US Implants	300	Polarity	Unipolar	Other	1
Advisories	None	Steroid	Yes		

## Atrial Placement

### System Longevity Study Results

### Qualifying Complications 6 Total

Number of Leads Enrolled in Study	109	Failure to Capture	5
Cumulative Months of Follow-Up	6,490	Oversensing	1
Number of Leads Active in Study	39		



# Pacing Leads continued

## 4068 CapSureFix

### Product Characteristics

US Market Release	Mar-96	Serial Number Prefix	LCE	<u>US Returned Product Analysis</u>	
Registered US Implants	124,800	Type and/or Fixation	Transvenous, V or A, Screw-in	Implant Damage	406
Estimated Active US Implants	43,900	Polarity	Bipolar	Electrical Malfunction	111
Advisories	None	Steroid	Yes	Other	11


### Atrial Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	2,411
Cumulative Months of Follow-Up	124,306
Number of Leads Active in Study	552

#### Qualifying Complications 60 Total

Conductor Fracture	2	Insulation (ESC)	2
Extra Cardiac Stimulation	1	Insulation (not further defined)	1
Failure to Capture	19	Lead Dislodgement	8
Failure to Sense	11	Oversensing	7
Impedance Out of Range	6	Unspecified Clinical Failure	3




### Ventricular Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	1,799
Cumulative Months of Follow-Up	89,241
Number of Leads Active in Study	457

#### Qualifying Complications 38 Total

Conductor Fracture	2	Impedance Out of Range	5
Extra Cardiac Stimulation	2	Oversensing	4
Failure to Capture	20	Unspecified Clinical Failure	2
Failure to Sense	3		



# Pacing Leads continued

## 4073 CapSure Sense

### Product Characteristics


Product Characteristics				US Returned Product Analysis	
US Market Release	Jun-02	Serial Number Prefix	BBF	Implant Damage	1
Registered US Implants	600	Type and/or Fixation	Transvenous, Vent., Tines	Electrical Malfunction	0
Estimated Active US Implants	400	Polarity	Unipolar	Other	0
Advisories	None	Steroid	Yes		

## Atrial Placement

### System Longevity Study Results

Qualifying Complications 0 Total

Number of Leads Enrolled in Study	1
Cumulative Months of Follow-Up	52
Number of Leads Active in Study	1




## Ventricular Placement

### System Longevity Study Results

Qualifying Complications 0 Total

Number of Leads Enrolled in Study	100
Cumulative Months of Follow-Up	4,770
Number of Leads Active in Study	83



# Pacing Leads continued

## 4074 CapSure Sense

### Product Characteristics


				US Returned Product Analysis	
US Market Release	Jun-02	Serial Number Prefix	BBD	Implant Damage	13
Registered US Implants	63,000	Type and/or Fixation	Transvenous, Vent., Tines	Electrical Malfunction	8
Estimated Active US Implants	45,200	Polarity	Bipolar	Other	1
Advisories	None	Steroid	Yes		

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications

		5	Total		
Number of Leads Enrolled in Study	622	Failure to Capture	1	Lead Dislodgement	2
Cumulative Months of Follow-Up	24,361	Failure to Sense	1		
Number of Leads Active in Study	516	Impedance Out of Range	1		



# Pacing Leads continued

## 4076 CapSureFix Novus

### Product Characteristics

US Market Release	Feb-04	Serial Number Prefix	BBL	<b>US Returned Product Analysis</b>	
Registered US Implants	199,200	Type and/or Fixation	Transvenous, V or A, Screw-in	Implant Damage	96
Estimated Active US Implants	165,500	Polarity	Bipolar	Electrical Malfunction	13
Advisories	None	Steroid	Yes	Other	8


### Atrial Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	745
Cumulative Months of Follow-Up	17,688
Number of Leads Active in Study	618

#### Qualifying Complications 3 Total

Failure to Capture	1
Lead Dislodgement	2




### Ventricular Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	668
Cumulative Months of Follow-Up	17,812
Number of Leads Active in Study	555

#### Qualifying Complications 2 Total

Failure to Capture	2
--------------------	---



Leads

# Pacing Leads continued

## 4081 Target Tip

### Product Characteristics


US Market Release	Jul-89	Serial Number Prefix	LAC	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	3,900	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	4
Estimated Active US Implants	800	Polarity	Unipolar	Electrical Malfunction	5
Advisories	None	Steroid	No	Other	0

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 3 Total

Number of Leads Enrolled in Study	260	Conductor Fracture	1
Cumulative Months of Follow-Up	9,940	Failure to Sense	2
Number of Leads Active in Study	9		



## 4092 CapSure SP Novus

### Product Characteristics


US Market Release	Sep-98	Serial Number Prefix	LEP	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	151,000	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	39
Estimated Active US Implants	83,100	Polarity	Bipolar	Electrical Malfunction	19
Advisories	None	Steroid	Yes	Other	5

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 17 Total

Number of Leads Enrolled in Study	1,145	Conductor Fracture	3	Impedance Out of Range	1
Cumulative Months of Follow-Up	55,839	Extra Cardiac Stimulation	1	Lead Dislodgement	4
Number of Leads Active in Study	512	Failure to Capture	8		



# Pacing Leads continued

## 4503, 4503M CapSure

### Product Characteristics


US Market Release	Jul-86	Serial Number Prefix	MQ, LAY	<u>US Returned Product Analysis</u>	
Registered US Implants	8,000	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	2
Estimated Active US Implants	1,500	Polarity	Unipolar	Electrical Malfunction	12
Advisories	None	Steroid	Yes	Other	0

### Atrial Placement

#### System Longevity Study Results

#### Qualifying Complications 1 Total

Number of Leads Enrolled in Study	59	Failure to Sense	1
Cumulative Months of Follow-Up	3,278		
Number of Leads Active in Study	6		



## 4504, 4504M CapSure

### Product Characteristics


US Market Release	Mar-90	Serial Number Prefix	QM or LBA	<u>US Returned Product Analysis</u>	
Registered US Implants	15,400	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	5
Estimated Active US Implants	1,700	Polarity	Bipolar	Electrical Malfunction	172
Advisories	1	Steroid	Yes	Other	4
See page 156 – 1996 Lead Survival Below Expectations					

### Atrial Placement

#### System Longevity Study Results

#### Qualifying Complications 48 Total

Number of Leads Enrolled in Study	368	Electrical Abandonment	3	Impedance Out of Range	9
Cumulative Months of Follow-Up	19,879	Extra Cardiac Stimulation	1	Insulation (MIO)	1
Number of Leads Active in Study	1	Failure to Capture	14	Lead Dislodgement	1
		Failure to Sense	16	Oversensing	3



	1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	at 105 mo
%	100.0	100.0	99.1	98.2	90.3	82.2	73.0	69.9	66.1
#	294	260	220	186	145	109	74	59	47
Effective Sample Size									



# Pacing Leads continued

## 4512 Target Tip

### Product Characteristics


US Market Release	Jul-83	Serial Number Prefix	PF	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	10,300	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	4
Estimated Active US Implants	1,200	Polarity	Bipolar	Electrical Malfunction	85
Advisories	None	Steroid	No	Other	8

## Atrial Placement

### System Longevity Study Results

### Qualifying Complications

Number of Leads Enrolled in Study	600	Electrical Abandonment	1	Insulation (MIO)	4
Cumulative Months of Follow-Up	39,833	Failure to Capture	6	Insulation (not further defined)	2
Number of Leads Active in Study	4	Failure to Sense	14	Lead Dislodgement	1
		Impedance Out of Range	3	Oversensing	2
		Insulation (ESC)	2		



## 4523 CapSure SP

### Product Characteristics


US Market Release	Aug-91	Serial Number Prefix	ZE	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	11,200	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	5
Estimated Active US Implants	3,200	Polarity	Unipolar	Electrical Malfunction	2
Advisories	None	Steroid	Yes	Other	1

## Atrial Placement

### System Longevity Study Results

### Qualifying Complications

Number of Leads Enrolled in Study	121	Impedance Out of Range	1
Cumulative Months of Follow-Up	7,194	Lead Dislodgement	2
Number of Leads Active in Study	15	Oversensing	1



# Pacing Leads continued

## 4524 CapSure SP

### Product Characteristics

US Market Release	Oct-91	Serial Number Prefix	LAR	<u>US Returned Product Analysis</u>	
Registered US Implants	101,700	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	47
Estimated Active US Implants	32,300	Polarity	Bipolar	Electrical Malfunction	31
Advisories	None	Steroid	Yes	Other	8

### Atrial Placement


#### System Longevity Study Results

Number of Leads Enrolled in Study	911
Cumulative Months of Follow-Up	40,567
Number of Leads Active in Study	57

#### Qualifying Complications

6 Total

Failure to Capture	3
Failure to Sense	2
Lead Dislodgement	1



Leads

## 4533 CapSure Z

### Product Characteristics

US Market Release	Not US released	Serial Number Prefix	LCB	<u>US Returned Product Analysis</u>	
Registered US Implants	NA	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	0
Estimated Active US Implants	NA	Polarity	Unipolar	Electrical Malfunction	0
Advisories	None	Steroid	Yes	Other	0

### Atrial Placement


#### System Longevity Study Results

Number of Leads Enrolled in Study	206
Cumulative Months of Follow-Up	11,286
Number of Leads Active in Study	16

#### Qualifying Complications

4 Total

Failure to Capture	1	Oversensing	1
Failure to Sense	1		
Lead Dislodgement	1		



# Pacing Leads continued

## 4557, 4557M Screw-In

### Product Characteristics


US Market Release	Aug-88	Serial Number Prefix	VQ or LAM	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	19,700	Type and/or Fixation	Transvenous, Atrial-J, Screw-in	Implant Damage	53
Estimated Active US Implants	4,400	Polarity	Unipolar	Electrical Malfunction	14
Advisories	None	Steroid	No	Other	4

### Atrial Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	294	Extra Cardiac Stimulation	1	Oversensing	1
Cumulative Months of Follow-Up	18,465	Failure to Capture	3		
Number of Leads Active in Study	10	Failure to Sense	1		

#### Qualifying Complications



## 4558M Screw-In

### Product Characteristics


US Market Release	Nov-94	Serial Number Prefix	LDC	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	20,000	Type and/or Fixation	Transvenous, Atrial-J, Screw-in	Implant Damage	111
Estimated Active US Implants	5,400	Polarity	Bipolar	Electrical Malfunction	12
Advisories	None	Steroid	No	Other	1

### Atrial Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	539	Electrical Abandonment	1	Impedance Out of Range	2
Cumulative Months of Follow-Up	22,441	Failure to Capture	3	Insulation (not further defined)	1
Number of Leads Active in Study	25	Failure to Sense	2	Oversensing	2

#### Qualifying Complications



# Pacing Leads continued

## 4568 CapSureFix

### Product Characteristics


US Market Release	Jan-97	Serial Number Prefix	LDD	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	69,700	Type and/or Fixation	Transvenous, Atrial-J, Screw-in	Implant Damage	198
Estimated Active US Implants	31,700	Polarity	Bipolar	Electrical Malfunction	18
Advisories	None	Steroid	Yes	Other	4

### Atrial Placement

#### System Longevity Study Results

#### Qualifying Complications 33 Total

Number of Leads Enrolled in Study	587	Failure to Capture	18	Lead Dislodgement	9
Cumulative Months of Follow-Up	24,835	Failure to Sense	3	Medical Judgment	1
Number of Leads Active in Study	213	Impedance Out of Range	2		



Leads

## 4574 CapSure Sense

### Product Characteristics


US Market Release	Jun-02	Serial Number Prefix	BBE	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	39,900	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	7
Estimated Active US Implants	30,300	Polarity	Bipolar	Electrical Malfunction	2
Advisories	None	Steroid	Yes	Other	0

### Atrial Placement

#### System Longevity Study Results

#### Qualifying Complications 0 Total

Number of Leads Enrolled in Study	18
Cumulative Months of Follow-Up	325
Number of Leads Active in Study	16



# Pacing Leads continued

## 4592 CapSure SP Novus

### Product Characteristics


US Market Release	Oct-98	Serial Number Prefix	LER	<u>US Returned Product Analysis</u>	
Registered US Implants	74,400	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	13
Estimated Active US Implants	43,500	Polarity	Bipolar	Electrical Malfunction	5
Advisories	None	Steroid	Yes	Other	0

### Atrial Placement

#### System Longevity Study Results

#### Qualifying Complications

Number of Leads Enrolled in Study	244	Failure to Capture	2
Cumulative Months of Follow-Up	11,373	Failure to Sense	1
Number of Leads Active in Study	84	Lead Dislodgement	2



## 5023, 5023M CapSure SP

### Product Characteristics


US Market Release	Nov-88	Serial Number Prefix	SX or LAS	<u>US Returned Product Analysis</u>	
Registered US Implants	9,900	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	15
Estimated Active US Implants	2,300	Polarity	Unipolar	Electrical Malfunction	7
Advisories	None	Steroid	Yes	Other	0

### Ventricular Placement

#### System Longevity Study Results

#### Qualifying Complications

Number of Leads Enrolled in Study	1,353	Conductor Fracture	2	Impedance Out of Range	1
Cumulative Months of Follow-Up	70,587	Extra Cardiac Stimulation	4		
Number of Leads Active in Study	470	Failure to Capture	8		



# Pacing Leads continued

## 5024, 5024M CapSure SP

### Product Characteristics

US Market Release	Mar-90	Serial Number Prefix	SY or LAT	<u>US Returned Product Analysis</u>	
Registered US Implants	201,500	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	723
Estimated Active US Implants	56,200	Polarity	Bipolar	Electrical Malfunction	115
Advisories	None	Steroid	Yes	Other	29


### Ventricular Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	8,140
Cumulative Months of Follow-Up	431,840
Number of Leads Active in Study	594

#### Qualifying Complications

Conductor Fracture	3	Insulation (ESC)	1
Extra Cardiac Stimulation	2	Insulation (not further defined)	5
Failure to Capture	25	Lead Dislodgement	5
Failure to Sense	2	Oversensing	1
Impedance Out of Range	3	Unspecified Clinical Failure	1



Leads

## 5026 CapSure

### Product Characteristics

US Market Release	Feb-88	Serial Number Prefix	RZ	<u>US Returned Product Analysis</u>	
Registered US Implants	7,400	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	60
Estimated Active US Implants	1,100	Polarity	Bipolar	Electrical Malfunction	7
Advisories	None	Steroid	Yes	Other	1


### Ventricular Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	168
Cumulative Months of Follow-Up	9,611
Number of Leads Active in Study	5

#### Qualifying Complications

Electrical Abandonment	1
Failure to Capture	3



# Pacing Leads continued

## 5033 CapSure Z

### Product Characteristics


US Market Release	Feb-96	Serial Number Prefix	LDK	<u>US Returned Product Analysis</u>	
Registered US Implants	2,400	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	6
Estimated Active US Implants	700	Polarity	Unipolar	Electrical Malfunction	1
Advisories	None	Steroid	Yes	Other	3

### Ventricular Placement

#### System Longevity Study Results

#### Qualifying Complications

Number of Leads Enrolled in Study	1,899	Cardiac Perforation	1	Impedance Out of Range	4
Cumulative Months of Follow-Up	96,195	Conductor Fracture	7	Insulation (not further defined)	1
Number of Leads Active in Study	263	Failure to Capture	11	Lead Dislodgement	2



## 5034 CapSure Z

### Product Characteristics


US Market Release	Feb-96	Serial Number Prefix	LDF	<u>US Returned Product Analysis</u>	
Registered US Implants	56,300	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	85
Estimated Active US Implants	17,000	Polarity	Bipolar	Electrical Malfunction	31
Advisories	None	Steroid	Yes	Other	11

### Ventricular Placement

#### System Longevity Study Results

#### Qualifying Complications

Number of Leads Enrolled in Study	1,594	Conductor Fracture	1	Impedance Out of Range	1
Cumulative Months of Follow-Up	82,761	Failure to Capture	9	Lead Dislodgement	1
Number of Leads Active in Study	228	Failure to Sense	1		



# Pacing Leads continued

## 5054 CapSure Z Novus

### Product Characteristics


US Market Release	Jun-98	Serial Number Prefix	LEH	<u>US Returned Product Analysis</u>	
Registered US Implants	85,400	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	43
Estimated Active US Implants	44,900	Polarity	Bipolar	Electrical Malfunction	16
Advisories	None	Steroid	Yes	Other	6

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 11 Total

Number of Leads Enrolled in Study	1,392	Failure to Capture	7	Lead Dislodgement	2
Cumulative Months of Follow-Up	63,013	Failure to Sense	1		
Number of Leads Active in Study	463	Impedance Out of Range	1		



Leads



# Pacing Leads continued

## 5068 CapSureFix

### Product Characteristics

US Market Release	Jan-97	Serial Number Prefix	LDJ	<u>US Returned Product Analysis</u>	
Registered US Implants	103,100	Type and/or Fixation	Transvenous, V or A, Screw-in	Implant Damage	455
Estimated Active US Implants	40,600	Polarity	Bipolar	Electrical Malfunction	75
Advisories	None	Steroid	Yes	Other	15


### Atrial Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	968
Cumulative Months of Follow-Up	32,929
Number of Leads Active in Study	65

#### Qualifying Complications

Failure to Capture	2	Oversensing	1
Impedance Out of Range	2		
Lead Dislodgement	1		




### Ventricular Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	1,359
Cumulative Months of Follow-Up	37,307
Number of Leads Active in Study	119

#### Qualifying Complications

Conductor Fracture	1	Lead Dislodgement	1
Failure to Capture	3		
Insulation (not further defined)	1		



# Pacing Leads continued

## 5072 SureFix

### Product Characteristics


				US Returned Product Analysis	
US Market Release	Jun-98	Serial Number Prefix	LEM	Implant Damage	28
Registered US Implants	8,900	Type and/or Fixation	Transvenous, V or A, Screw-in	Electrical Malfunction	5
Estimated Active US Implants	4,500	Polarity	Bipolar	Other	1
Advisories	None	Steroid	Yes		

## Atrial Placement

### System Longevity Study Results

### Qualifying Complications 2 Total

Number of Leads Enrolled in Study	451	Cardiac Perforation	1
Cumulative Months of Follow-Up	21,911	Failure to Capture	1
Number of Leads Active in Study	114		



# Pacing Leads continued

## 5076 CapSureFix Novus

### Product Characteristics


US Market Release	Aug-00	Serial Number Prefix	PJN	<u>US Returned Product Analysis</u>	
Registered US Implants	1,025,600	Type and/or Fixation	Transvenous, V or A, Screw-in	Implant Damage	949
Estimated Active US Implants	723,800	Polarity	Bipolar	Electrical Malfunction	276
Advisories	None	Steroid	Yes	Other	84

### Atrial Placement

#### System Longevity Study Results

#### Qualifying Complications 16 Total

Number of Leads Enrolled in Study	2,678	Cardiac Perforation	1	Impedance Out of Range	2
Cumulative Months of Follow-Up	104,689	Conductor Fracture	1	Insulation (not further defined)	1
Number of Leads Active in Study	1,207	Extra Cardiac Stimulation	2	Lead Dislodgement	4
		Failure to Capture	4	Oversensing	1




### Ventricular Placement

#### System Longevity Study Results

#### Qualifying Complications 9 Total

Number of Leads Enrolled in Study	1,521	Cardiac Perforation	1	Failure to Sense	1
Cumulative Months of Follow-Up	55,080	Conductor Fracture	1	Impedance Out of Range	1
Number of Leads Active in Study	556	Failure to Capture	3	Lead Dislodgement	2



# Pacing Leads continued

## 5092 CapSure SP Novus

### Product Characteristics


US Market Release	Jun-98	Serial Number Prefix	LET	<u>US Returned Product Analysis</u>	
Registered US Implants	111,800	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	48
Estimated Active US Implants	61,600	Polarity	Bipolar	Electrical Malfunction	28
Advisories	None	Steroid	Yes	Other	11

### Ventricular Placement

#### System Longevity Study Results

#### Qualifying Complications 8 Total

Number of Leads Enrolled in Study	1,171	Extra Cardiac Stimulation	1
Cumulative Months of Follow-Up	41,628	Failure to Capture	2
Number of Leads Active in Study	227	Lead Dislodgement	5



## 5524, 5524M CapSure SP

### Product Characteristics


US Market Release	Mar-90	Serial Number Prefix	XV or LAV	<u>US Returned Product Analysis</u>	
Registered US Implants	60,600	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	67
Estimated Active US Implants	21,000	Polarity	Bipolar	Electrical Malfunction	25
Advisories	None	Steroid	Yes	Other	7

### Atrial Placement

#### System Longevity Study Results

#### Qualifying Complications 38 Total

Number of Leads Enrolled in Study	4,445	Conductor Fracture	1	Insulation (not further defined)	2
Cumulative Months of Follow-Up	241,326	Failure to Capture	22	Lead Dislodgement	4
Number of Leads Active in Study	482	Failure to Sense	4	Oversensing	4
		Impedance Out of Range	1		



# Pacing Leads continued

## 5534 CapSure Z

### Product Characteristics

US Market Release	Feb-96	Serial Number Prefix	LDG	<u>US Returned Product Analysis</u>	
Registered US Implants	26,300	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	29
Estimated Active US Implants	9,300	Polarity	Bipolar	Electrical Malfunction	8
Advisories	None	Steroid	Yes	Other	5

### Atrial Placement


#### System Longevity Study Results

Number of Leads Enrolled in Study	261
Cumulative Months of Follow-Up	12,779
Number of Leads Active in Study	23

#### Qualifying Complications

Failure to Capture	5
Impedance Out of Range	1

6 Total



## 5554 CapSure Z Novus

### Product Characteristics

US Market Release	Jun-98	Serial Number Prefix	LEJ	<u>US Returned Product Analysis</u>	
Registered US Implants	54,800	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	8
Estimated Active US Implants	31,500	Polarity	Bipolar	Electrical Malfunction	12
Advisories	None	Steroid	Yes	Other	4

### Atrial Placement


#### System Longevity Study Results

Number of Leads Enrolled in Study	338
Cumulative Months of Follow-Up	18,883
Number of Leads Active in Study	60

#### Qualifying Complications

Failure to Capture	1	Oversensing	1
Impedance Out of Range	1		
Lead Dislodgement	1		

4 Total



# Pacing Leads continued

## 5568 CapSureFix

### Product Characteristics

US Market Release	Jan-97	Serial Number Prefix	LDN	<u>US Returned Product Analysis</u>	
Registered US Implants	70,400	Type and/or Fixation	Transvenous, Atrial-J, Screw-in	Implant Damage	264
Estimated Active US Implants	45,000	Polarity	Bipolar	Electrical Malfunction	15
Advisories	None	Steroid	Yes	Other	12


### Atrial Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	920
Cumulative Months of Follow-Up	30,230
Number of Leads Active in Study	184

#### Qualifying Complications

Failure to Capture	2
Failure to Sense	1
Lead Dislodgement	1



## 5592 CapSure SP Novus

### Product Characteristics

US Market Release	Jun-98	Serial Number Prefix	LEU	<u>US Returned Product Analysis</u>	
Registered US Implants	27,800	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	6
Estimated Active US Implants	18,100	Polarity	Bipolar	Electrical Malfunction	3
Advisories	None	Steroid	Yes	Other	0


### Atrial Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	666
Cumulative Months of Follow-Up	26,129
Number of Leads Active in Study	179

#### Qualifying Complications

Failure to Capture	2
Lead Dislodgement	2



# Pacing Leads continued

## 5594 CapSure SP Novus

### Product Characteristics


US Market Release	Jun-01	Serial Number Prefix	LFD	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	10,900	Type and/or Fixation	Transvenous, Atrial-J, Tines	Implant Damage	0
Estimated Active US Implants	8,100	Polarity	Bipolar	Electrical Malfunction	4
Advisories	None	Steroid	Yes	Other	0

### Atrial Placement

#### System Longevity Study Results

Qualifying Complications 0 Total

Number of Leads Enrolled in Study	18
Cumulative Months of Follow-Up	1,068
Number of Leads Active in Study	12



## 6940 CapSureFix

### Product Characteristics


US Market Release	Oct-98	Serial Number Prefix	TCP	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	25,500	Type and/or Fixation	Transvenous, A or V, Screw-in	Implant Damage	114
Estimated Active US Implants	10,100	Polarity	Bipolar	Electrical Malfunction	21
Advisories	None	Steroid	Yes	Other	3

### Atrial Placement

#### System Longevity Study Results

Qualifying Complications 7 Total

Number of Leads Enrolled in Study	818	Conductor Fracture	1	Oversensing	3
Cumulative Months of Follow-Up	38,362	Failure to Sense	2		
Number of Leads Active in Study	151	Lead Dislodgement	1		



# Pacing Leads continued

## 6957 Spectraflex

### Product Characteristics


US Market Release	Jul-79	Serial Number Prefix	VC	<u>US Returned Product Analysis</u>	
Registered US Implants	29,100	Type and/or Fixation	Transvenous, V or A, Screw-in	Implant Damage	85
Estimated Active US Implants	2,000	Polarity	Unipolar	Electrical Malfunction	39
Advisories	None	Steroid	No	Other	25

### Atrial Placement

#### System Longevity Study Results

#### Qualifying Complications 10 Total

Number of Leads Enrolled in Study	673	Extra Cardiac Stimulation	1	Oversensing	1
Cumulative Months of Follow-Up	24,255	Failure to Capture	3		
Number of Leads Active in Study	1	Failure to Sense	5		




### Ventricular Placement

#### System Longevity Study Results

#### Qualifying Complications 42 Total

Number of Leads Enrolled in Study	1,853	Conductor Fracture	14	Impedance Out of Range	1
Cumulative Months of Follow-Up	96,335	Extra Cardiac Stimulation	2	Insulation (not further defined)	1
Number of Leads Active in Study	17	Failure to Capture	18	Oversensing	4
		Failure to Sense	2		



Leads



# Pacing Leads continued

## 6957J Spectraflex

### Product Characteristics

US Market Release	Sep-80	Serial Number Prefix	GG	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	30,000	Type and/or Fixation	Transvenous, Atrial-J, Screw-in	Implant Damage	74
Estimated Active US Implants	2,100	Polarity	Unipolar	Electrical Malfunction	30
Advisories	None	Steroid	No	Other	30


### Atrial Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	2,348
Cumulative Months of Follow-Up	160,477
Number of Leads Active in Study	22

#### Qualifying Complications

Conductor Fracture	13	Insulation (ESC)	1
Extra Cardiac Stimulation	3	Insulation (not further defined)	3
Failure to Capture	48	Lead Dislodgement	2
Failure to Sense	14	Oversensing	3
Impedance Out of Range	1		



## 6961 Tenax

### Product Characteristics

US Market Release	Jan-78	Serial Number Prefix	TB	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	44,700	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	103
Estimated Active US Implants	2,100	Polarity	Unipolar	Electrical Malfunction	27
Advisories	None	Steroid	No	Other	0


### Ventricular Placement

#### System Longevity Study Results

Number of Leads Enrolled in Study	627
Cumulative Months of Follow-Up	42,864
Number of Leads Active in Study	0

#### Qualifying Complications

Extra Cardiac Stimulation	4	Insulation (not further defined)	2
Failure to Capture	8	Lead Dislodgement	1
Failure to Sense	6	Oversensing	1



# Pacing Leads continued

## 6962 Tenax

### Product Characteristics


US Market Release	Jan-78	Serial Number Prefix	UB	<u>US Returned Product Analysis</u>	
Registered US Implants	70,600	Type and/or Fixation	Transvenous, Vent., Tines	Implant Damage	170
Estimated Active US Implants	3,200	Polarity	Bipolar	Electrical Malfunction	84
Advisories	None	Steroid	No	Other	0

## Ventricular Placement

### System Longevity Study Results

### Qualifying Complications 52 Total

Number of Leads Enrolled in Study	1,483	Conductor Fracture	5	Impedance Out of Range	3
Cumulative Months of Follow-Up	109,942	Extra Cardiac Stimulation	1	Insulation (not further defined)	2
Number of Leads Active in Study	2	Failure to Capture	27	Lead Dislodgement	1
		Failure to Sense	10	Oversensing	3



Leads









## US Returned Product Analysis Summary

Model Number	Family	US Market Release	Estimated US Implants	Estimated US Active	Implant Damage	Electrical Malfunction	Other
3830	SelectSecure	Aug-05	12,200	10,600	21	5	1
4003, 4003M	CapSure	Jul-86	38,000	5,300	24	60	2
4004, 4004M	CapSure	Feb-89	72,600	5,800	56	688	19
4011	Target Tip	Nov-82	58,400	6,200	29	152	5
4012	Target Tip	Jul-83	93,700	6,500	50	827	34
4023	CapSure SP	Aug-91	41,200	9,600	47	21	6
4024	CapSure SP	Oct-91	222,100	56,100	264	135	34
4033	CapSure Z	Not US released	NA	NA	2	0	0
4057, 4057M	Screw-in	Aug-88	10,100	1,800	39	6	4
4058, 4058M	Screw-in	Jan-89	101,900	20,800	388	261	23
4067	CapSureFix	Jan-97	1,000	300	3	1	1
4068	CapSureFix	Mar-96	124,800	43,900	406	111	11
4073	CapSure Sense	Jun-02	600	400	1	0	0
4074	CapSure Sense	Jun-02	63,000	45,200	13	8	1
4076	CapSureFix Novus	Feb-04	199,200	165,500	96	13	8
4081	Target Tip	Jul-89	3,900	800	4	5	0
4092	CapSure SP Novus	Sep-98	151,000	83,100	39	19	5
4503, 4503M	CapSure	Jul-86	8,000	1,500	2	12	0
4504, 4504M	CapSure	Mar-90	15,400	1,700	5	172	4
4512	Target Tip	Jul-83	10,300	1,200	4	85	8
4523	CapSure SP	Aug-91	11,200	3,200	5	2	1
4524	CapSure SP	Oct-91	101,700	32,300	47	31	8
4533	CapSure Z	Not US released	NA	NA	0	0	0
4557, 4557M	Screw-in	Aug-88	19,700	4,400	53	14	4
4558M	Screw-in	Nov-94	20,000	5,400	111	12	1
4568	CapSureFix	Jan-97	69,700	31,700	198	18	4
4574	CapSure Sense	Jun-02	39,900	30,300	7	2	0
4592	CapSure SP Novus	Oct-98	74,400	43,500	13	5	0
5023, 5023M	CapSure SP	Nov-88	9,900	2,300	15	7	0
5024, 5024M	CapSure SP	Mar-90	201,500	56,200	723	115	29
5026	CapSure	Feb-88	7,400	1,100	60	7	1
5033	CapSure Z	Feb-96	2,400	700	6	1	3
5034	CapSure Z	Feb-96	56,300	17,000	85	31	11
5054	CapSure Z Novus	Jun-98	85,400	44,900	43	16	6
5068	CapSureFix	Jan-97	103,100	40,600	455	75	15
5072	SureFix	Jun-98	8,900	4,500	28	5	1
5076	CapSureFix Novus	Aug-00	1,025,600	723,800	949	276	84
5092	CapSure SP Novus	Jun-98	111,800	61,600	48	28	11
5524, 5524M	CapSure SP	Mar-90	60,600	21,000	67	25	7
5534	CapSure Z	Feb-96	26,300	9,300	29	8	5
5554	CapSure Z Novus	Jun-98	54,800	31,500	8	12	4
5568	CapSureFix	Jan-97	70,400	45,000	264	15	12

continued

## US Returned Product Analysis Summary continued

Model Number	Family	US Market Release	Estimated US Implants	Estimated US Active	Implant Damage	Electrical Malfunction	Other
5592	CapSure SP Novus	Jun-98	27,800	18,100	6	3	0
5594	CapSure SP Novus	Jun-01	10,908	1,700	0	4	0
6940	CapSureFix	Oct-98	25,500	10,100	114	21	3
6957	Spectraflex	Jul-79	29,100	2,000	85	39	25
6957J	Spectraflex	Sep-80	30,000	2,100	74	30	30
6961	Tenax	Jan-78	44,700	2,100	103	27	0
6962	Tenax	Jan-78	70,600	3,200	170	84	0



## Reference Chart

Model Number	Family	Type	Insulation	Conductor Material	Tip Electrode	Connector Type
3830	SelectSecure	Transvenous V or A Screw-In	Polyurethane/Silicone (55D,4719)	MP35N 5 Filars/ Cable	1.8 mm Helix/Steroid	IS-1 BI
4003, 4003M	CapSure	Transvenous Ventricular Tines	Polyurethane (80A)	MP35N 4 Filars	Porous/Steroid	5 mm (4003) IS-1 UNI (4003M)
4004, 4004M	CapSure	Transvenous Ventricular Tines	Polyurethane (80A)	MP35N 6/4 Filars	Porous/Steroid	3.2 mm Low Profile (4004) IS-1 BI (4004M)
4011	Target Tip	Transvenous Ventricular Tines	Polyurethane (80A)	MP35N 4 Filars	Target Tip Concentric Grooves	5 mm
4012	Target Tip	Transvenous Ventricular Tines	Polyurethane (80A)	MP35N 6/4 Filars	Target Tip Concentric Grooves	3.2 mm Low Profile
4023	CapSure SP	Transvenous Ventricular Tines	Polyurethane (55D)	MP35N 4 Filars	Porous Platinized/ Steroid	IS-1 UNI
4024	CapSure SP	Transvenous Ventricular Tines	Polyurethane (55D)	MP35N 4/5 Filars	Porous Platinized/ Steroid	IS-1 BI
4033	CapSure Z	Transvenous Ventricular Tines	Polyurethane (55D)	MP35N 2 Filars	CapSure Z Platinized/ Steroid	IS-1 UNI
4057, 4057M	Screw-In	Transvenous V or A Screw-In	Polyurethane (80A)	MP35N 1 Filar	2.0 mm Helix	5 mm (4057) IS-1 UNI (4057M)
4058, 4058M	Screw-In	Transvenous V or A Screw-In	Polyurethane (80A/55D)	MP35N 4/1 Filars	2.0 mm Helix	3.2 mm Low Profile (4058) IS-1 BI (4058M)
4067	CapSureFix	Transvenous V or A Screw-In	Polyurethane (80A)	MP35N 3 Filars	1.8 mm Helix/Steroid	IS-1 UNI
4068	CapSureFix	Transvenous V or A Screw-In	Polyurethane (80A/55D)	MP35N 4/3 Filars	1.8 mm Helix/Steroid	IS-1 BI
4073	CapSure Sense	Transvenous Ventricular Tines	Polyurethane (55D)	MP35N 5 Filars	TiN Coated Platinum Iridium/ Steroid	IS-1 UNI
4074	CapSure Sense	Transvenous Ventricular Tines	Polyurethane/Silicone (55D, 4719)	MP35N 5/5 Filars	TiN Coated Platinum Iridium/ Steroid	IS-1 BI
4076	CapSureFix Novus	Transvenous V or A Screw-In	Polyurethane/Silicone (55D, 4719)	MP35N 4/6 Filars	TiN Coated Platinum Alloy/ Steroid	IS-1 BI
4081	Target Tip	Transvenous Ventricular Tines	Polyurethane (80A)	MP35N 4 Filars	Target Tip Concentric Grooves	IS-1 UNI w/Removable 5 mm Sleeve
4092	CapSure SP Novus	Transvenous Ventricular Tines	Polyurethane/Silicone (55D/4719)	MP35N 6/4 Filars	Porous Platinized/ Steroid	IS-1 BI
4503, 4503M	CapSure	Transvenous Atrial-J Tines	Polyurethane (80A)	MP35N 4 Filars	Porous/Steroid	5 mm (4503) IS-1 UNI (4503M)
4504, 4504M	CapSure	Transvenous Atrial-J Tines	Polyurethane (80A)	MP35N 3/4 Filars	Porous/Steroid	3.2 mm Low Profile (4504) IS-1 BI (4504M)
4512	Target Tip	Transvenous Atrial-J Tines	Polyurethane (80A)	MP35N 2/4 Filars	Target Tip Concentric Grooves	3.2 mm Low Profile
4523	CapSure SP	Transvenous Atrial-J Tines	Polyurethane (55D)	MP35N 2 Filars	Porous Platinized/ Steroid	IS-1 UNI
4524	CapSure SP	Transvenous Atrial-J Tines	Polyurethane (55D)	MP35N 4/5 Filars	Porous Platinized/ Steroid	IS-1 BI
4533	CapSure Z	Transvenous Atrial-J Tines	Polyurethane (55D)	MP35N 2 Filars	CapSure Z Platinized/ Steroid	IS-1 UNI
4557, 4557M	Screw-In	Transvenous Atrial-J Screw-In	Polyurethane (80A)	MP35N 1 Filar	1.5 mm Helix	5 mm (4557) IS-1 UNI (4557M)
4558M	Screw-In	Transvenous Atrial-J Screw-In	Polyurethane (80A/55D)	MP35N 6/3 Filars	1.8 mm Helix/Steroid	IS-1 BI
4568	CapSureFix	Transvenous Atrial-J Screw-In	Polyurethane (80A/55D)	MP35N 6/3 Filars	1.8 mm Helix/Steroid	IS-1 BI
4574	CapSure Sense	Transvenous Atrial-J Tines	Polyurethane/Silicone (55D,4719)	MP35N 5/5 Filars	TiN Coated Platinum Iridium	IS-1 BI
4592	CapSure SP Novus	Transvenous Atrial-J Tines	Polyurethane/Silicone (55D/4719)	MP35N 6/3 Filars	Porous Platinized/ Steroid	IS-1 BI
5023, 5023M	CapSure SP	Transvenous Ventricular Tines	Silicone	MP35N 4 Filars	Porous Platinized/ Steroid	5 mm (5023) IS-1 UNI (5023M)

continued

# Pacing Leads continued

## Reference Chart continued

Model Number	Family	Type	Insulation	Conductor Material	Tip Electrode	Connector Type
5024, 5024M	CapSure SP	Transvenous Ventricular Tines	Silicone	MP35N 4/5 Filars	Porous Platinized/ Steroid	3.2 mm Low Profile (5024) IS-1 BI (5024M)
5026	CapSure	Transvenous Ventricular Tines	Silicone	MP35N 6/4 Filars	Porous Platinized/ Steroid	3.2 mm Low Profile
5033	CapSure Z	Transvenous Ventricular Tines	Silicone	MP35N 4 Filars	CapSure Z Platinized/Steroid	IS-1 UNI
5034	CapSure Z	Transvenous Ventricular Tines	Silicone	MP35N 4/5 Filars	CapSure Z Platinized/Steroid	IS-1 BI
5054	CapSure Z Novus	Transvenous Ventricular Tines	Silicone (4719)	MP35N 5/5 Filars	CapSure Z Porous/Platinized/ Steroid	IS-1 BI
5068	CapSureFix	Transvenous V or A Screw-In	Silicone	MP35N 4/3 Filars	1.8 mm Helix/Steroid	IS-1 BI
5072	SureFix	Transvenous V or A Screw-In	Silicone	MP35N 4/5 Filars	1.8 mm Helix/Steroid	IS-1 BI
5076	CapSureFix Novus	Transvenous V or A Screw-In	Silicone (4719)	MP35N 4/6 Filars	Porous Platinized/ Steroid	IS-1 BI
5092	CapSure SP Novus	Transvenous Ventricular Tines	Silicone (4719)	MP35N 5/5 Filars	Porous Platinized/ Steroid	IS-1 BI
5524, 5524M	CapSure SP	Transvenous Atrial-J Tines	Silicone	MP35N 6/5 Filars	Porous Platinized/ Steroid	3.2 mm Low Profile (5524) IS-1 BI (5524M)
5534	CapSure Z	Transvenous Atrial-J Tines	Silicone	MP35N 4/5 Filars	CapSure Z Platinized/Steroid	IS-1 BI
5554	CapSure Z Novus	Transvenous Atrial-J Tines	Silicone (4719)	MP35N 6/5 Filars	CapSure Z Porous Platinized/ Steroid	IS-1 BI
5568	CapSureFix	Transvenous Atrial-J Screw-In	Silicone	MP35N 6/3 Filars	1.8 mm Helix/Steroid	IS-1 BI
5592	CapSure SP Novus	Transvenous Atrial-J Tines	Silicone (4719)	MP35N 6/5 Filars	Porous Platinized/ Steroid	IS-1 BI
5594	CapSure SP Novus	Transvenous Atrial-J Tines	Silicone (4719)	MP35N 6/5 Filars	Platinized Platinum/ Steroid	IS-1 BI
6940	CapSureFix	Transvenous A or V Screw-In	Silicone	MP35N 3/6 Filars	Platinum Alloy	IS-1 BI
6957	Spectraflex	Transvenous V or A Screw-In	Polyurethane (80A)	MP35N 1 Filar	2.0 mm Helix	5 mm
6957J	Spectraflex	Transvenous Atrial-J Screw-In	Polyurethane (80A)	MP35N 1 Filar	1.5 mm Helix	5 mm
6961	Tenax	Transvenous Ventricular Tines	Silicone	MP35N 3 Filars	Ring Tip	5 mm
6962	Tenax	Transvenous Ventricular Tines	Silicone	MP35N 4 Filars	Ring Tip	5 mm Bifurcated

# Epi/Myocardial Pacing Leads

## 4951, 4951M Spectraflex


### Product Characteristics

US Market Release	Oct-81	Serial Number Prefix	TF or LBJ	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	23,100	Type and/or Fixation	Myocardial Stab-in, V or A, Peds	Implant Damage	15
Estimated Active US Implants	2,500	Polarity	Unipolar	Electrical Malfunction	97
Advisories	None	Steroid	No	Other	28

### System Longevity Study Results

### Qualifying Complications 10 Total

Number of Leads Enrolled in Study	179	Failure to Capture	4	Insulation (ESC)	1
Cumulative Months of Follow-Up	6,512	Failure to Sense	3	Insulation (not further defined)	1
Number of Leads Active in Study	4	Impedance Out of Range	1		



## 4965 CapSure Epi


### Product Characteristics

US Market Release	Sep-96	Serial Number Prefix	LBT	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	17,700	Type and/or Fixation	Epicardial Suture-On V or A	Implant Damage	8
Estimated Active US Implants	8,900	Polarity	Unipolar	Electrical Malfunction	115
Advisories	None	Steroid	Yes	Other	2

### System Longevity Study Results

### Qualifying Complications 8 Total

Number of Leads Enrolled in Study	188	Conductor Fracture	3	Oversensing	2
Cumulative Months of Follow-Up	5,451	Failure to Capture	2		
Number of Leads Active in Study	28	Failure to Sense	1		



# Epi/Myocardial Pacing Leads continued

## 4968 CapSure Epi


### Product Characteristics

US Market Release	Sep-99	Serial Number Prefix	LEN	<u>US Returned Product Analysis</u>	
Registered US Implants	16,000	Type and/or Fixation	Epicardial Suture-On V or A	Implant Damage	2
Estimated Active US Implants	10,000	Polarity	Bipolar	Electrical Malfunction	11
Advisories	None	Steroid	Yes	Other	0

### System Longevity Study Results

### Qualifying Complications 34 Total

Number of Leads Enrolled in Study	543	Conductor Fracture	7	Impedance Out of Range	4
Cumulative Months of Follow-Up	22,604	Failure to Capture	14	Insulation (not further defined)	2
Number of Leads Active in Study	326	Failure to Sense	3	Oversensing	4



Leads

## 5071


### Product Characteristics

US Market Release	Dec-92	Serial Number Prefix	LAQ	<u>US Returned Product Analysis</u>	
Registered US Implants	33,900	Type and/or Fixation	Myocardial Screw-in Vent.	Implant Damage	29
Estimated Active US Implants	13,400	Polarity	Unipolar	Electrical Malfunction	7
Advisories	None	Steroid	No	Other	2

### System Longevity Study Results

### Qualifying Complications 11 Total

Number of Leads Enrolled in Study	229	Failure to Capture	9
Cumulative Months of Follow-Up	7,092	Oversensing	2
Number of Leads Active in Study	30		



# Epi/Myocardial Pacing Leads continued

## 6917, 6917A Tenax


### Product Characteristics

US Market Release	Jun-73	Serial Number Prefix	WV or WC	<u>US Returned Product Analysis</u>	
Registered US Implants	180,100	Type and/or Fixation	Myocardial Screw-in Vent.	Implant Damage	115
Estimated Active US Implants	5,000	Polarity	Unipolar	Electrical Malfunction	42
Advisories	None	Steroid	No	Other	1

### System Longevity Study Results

### Qualifying Complications 69 Total

Number of Leads Enrolled in Study	985	Conductor Fracture	6	Impedance Out of Range	2
Cumulative Months of Follow-Up	47,481	Extra Cardiac Stimulation	1	Insulation (MIO)	1
Number of Leads Active in Study	4	Failure to Capture	30	Oversensing	18
		Failure to Sense	11		



# Epi/Myocardial Pacing Leads continued

## Lead Survival Summary (95% Confidence Interval)

Model Number	Family	US Market Release	Leads Enrolled	Leads Active in Study	Qualifying Complications	Cumulative Months of Follow-Up in Study	Device Survival Probability (%)																
							Years After Implant																
							1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	10 yr	12 yr	14 yr	16 yr					
<b>4951, 4951M</b>	<b>Spectraflex</b>	Oct-81	179	4	10	6,512	977 +1.6/-4.8	96.5 +2.2/-5.8	96.5 +2.2/-5.8	93.4 +3.7/-8.1	93.4 +3.7/-8.1 at 57 mo												
<b>4965</b>	<b>CapSure Epi</b>	Sep-96	188	28	8	5,451	98.8 +0.9/-3.6	97.7 +1.6/-5.0	96.5 +2.2/-6.1	91.9 +4.3/-8.9 at 45 mo													
<b>4968</b>	<b>CapSure Epi</b>	Sep-99	543	326	34	22,604	99.6 +0.3/-1.4	97.4 +1.2/-2.3	95.9 +1.7/-2.7	94.5 +2.1/-3.4	93.2 +2.5/-4.0	91.7 +3.0/-4.7	91.7 +3.0/-4.7	91.7 +3.0/-4.7									
<b>5071</b>	(No brand name)	Dec-92	229	30	11	7,092	96.3 +2.1/-4.4	94.5 +2.8/-5.5	91.7 +4.0/-7.3	91.7 +4.0/-7.3	91.7 +4.0/-7.3	91.7 +4.0/-7.3 at 66 mo											
<b>6917, 6917A</b>	<b>Tenax</b>	Jun-73	985	4	69	47,481	99.0 +0.5/-1.0	97.6 +0.9/-1.5	95.9 +1.3/-2.0	93.6 +1.8/-2.5	92.8 +2.0/-2.6	91.1 +2.3/-3.1	89.5 +2.7/-3.4	87.2 +3.1/-4.1	84.2 +3.9/-4.8	83.2 +4.2/-5.3	80.4 +5.2/-6.8 at 150 mo						

# Epi/Myocardial Pacing Leads continued

## US Returned Product Analysis Summary

Model Number	Family	US Market Release	Estimated US Implants	Estimated US Active	Implant Damage	Electrical Malfunction	Other
4951, 4951M	Spectraflex	Oct-81	23,100	2,500	15	97	28
4965	CapSure Epi	Sep-96	17,700	8,900	8	115	2
4968	CapSure Epi	Sep-99	16,000	10,000	2	11	0
5071	(No brand name)	Dec-92	33,900	13,400	29	7	2
6917, 6917A	Tenax	Jun-73	180,100	5,000	115	42	1

Source: Returned Product Analysis  
Data as of January 31, 2009

## Reference Chart

Model Number	Family	Type	Insulation	Conductor Material	Tip Electrode	Connector Type
4951, 4951M	Spectraflex	Myocardial Stab-In V or A/Peds	Polyurethane (80A)	MP35N 4 Filars	Barb	5 mm (4951) IS-1 UNI (4951M)
4965	CapSure Epi	Epicardial Suture-On V or A	Silicone	MP35N 5 Filars	Porous Platinized/ Steroid	IS-1 UNI
4968	CapSure Epi	Epicardial Suture V or A	Silicone	MP35N 5 Filars	Porous Platinized/ Steroid	IS-1 B1
5069	(No brand name)	Myocardial Screw-In	Silicone	MP35N Multifilars	3-Turn Helix	IS-1 UNI
5071	(No brand name)	Myocardial Screw-In Ventricular	Silicone	MP35N Multifilars	2-Turn Helix	IS-1 UNI
6917	Tenax	Myocardial Screw-In Ventricular	Silicone	Pt Ir Tinsel Wire	3-Turn Helix	5 mm
6917A	Tenax	Myocardial Screw-In Ventricular	Silicone	Pt Ir Tinsel Wire	2-Turn Helix	5 mm

# VDD Single Pass Pacing Leads

## 5032 CapSure VDD


### Product Characteristics

US Market Release	Mar-96	Serial Number Prefix	LCL, LCN, LCM	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	5,400	Type and/or Fixation	Transvenous, Atr-Vent.,Tines	Implant Damage	24
Estimated Active US Implants	1,600	Polarity	Quadripolar	Electrical Malfunction	12
Advisories	None	Steroid	Yes	Other	0

### System Longevity Study Results

### Qualifying Complications

Number of Leads Enrolled in Study	38	Failure to Sense	1	Total	1
Cumulative Months of Follow-Up	2,011				
Number of Leads Active in Study	0				



## 5038 CapSure VDD-2


### Product Characteristics

US Market Release	Sep-98	Serial Number Prefix	LEE, LEG, or LEF	<a href="#">US Returned Product Analysis</a>	
Registered US Implants	8,200	Type and/or Fixation	Transvenous, Atr-Vent.,Tines	Implant Damage	6
Estimated Active US Implants	3,800	Polarity	Quadripolar	Electrical Malfunction	5
Advisories	None	Steroid	Yes	Other	1

### System Longevity Study Results

### Qualifying Complications

Number of Leads Enrolled in Study	552	Conductor Fracture	2	Total	5
Cumulative Months of Follow-Up	21,740	Failure to Capture	1		
Number of Leads Active in Study	92	Failure to Sense	2		



Leads



# VDD Single Pass Pacing Leads continued

## Lead Survival Summary (95% Confidence Interval)

Model Number	Family	US Market Release	Leads Enrolled	Leads Active in Study	Qualifying Complications	Cumulative Months of Follow-Up in Study	Device Survival Probability (%)									
							Years After Implant									
							1 yr	2 yr	3 yr	4 yr	5 yr	6 yr	7 yr	8 yr	9 yr	10 yr
5032	CapSure VDD	Mar-96	38	0	1	2,011	Survival estimate not available due to insufficient sample size									
5038	CapSure VDD-2	Sep-98	552	92	5	21,740	99.8 +0.2/-1.4	99.8 +0.2/-1.4	99.8 +0.2/-1.4	98.8 +0.8/-2.6	98.2 +1.2/-3.1	98.2 +1.2/-3.1	97.1 +1.8/-4.5	97.1 +1.8/-4.5 at 93 mo		

Source: System Longevity Study  
Data as of January 31, 2009

## US Returned Product Analysis Summary

Model Number	Family	US Market Release	Estimated US Implants	Estimated US Active	Implant Damage	Electrical Malfunction	Other
5032	CapSure VDD	Mar-96	5,400	1,600	24	12	0
5038	CapSure VDD-2	Sep-98	8,200	3,800	6	5	1

Source: Returned Product Analysis  
Data as of January 31, 2009

## Reference Chart

Model Number	Family	Type	Insulation	Conductor Material	Tip Electrode	Connector Type
5032	CapSure VDD	Transvenous V and A Tines	Silicone	MP35N 5/6/1 Filars	Porous Platinized/ Steroid	Atr. IS-1 BI, Vent. IS-1 BI
5038	CapSure VDD-2	Transvenous Ventricular Tines	Silicone	MP35N	Porous Platinized/ Steroid	Atr. IS-1 BI, Vent. IS-1 BI

# ICD and CRT-D Charge Time Performance

Medtronic continues its commitment to providing updated information on charge time performance.

## Introduction

Information on charge time performance of Medtronic is presented in this section of the CRDM Product Performance Report. Medtronic implemented the collection of charge time data on July 1, 1999. The data are collected via our ongoing active clinical study of long-term system performance called the System Longevity Study. The study protocol requests device data be routinely taken and sent to Medtronic at no more than 6-month intervals.

In our analysis performed for this report, only charge times resulting from full energy charges are considered. To ensure consistent reporting across devices, the charge time reported at implant represents the last charge time available from date of implant. When more than one charge time is available in a 6-month interval, a conservative approach has been adopted whereby only the maximum charge in each 6-month interval is reported. As charge time is directly proportional to the time elapsed since the last capacitor reformation, charges occurring within 15 days of a previous charge are excluded. This precludes the reporting of overly optimistic results.



Data from over 20,000 devices contribute to the charge time data in this report. By tracking and reporting this charge time data, Medtronic is able to ascertain the actual performance of its charging circuitry. The insight gained through this information is applied to Medtronic's ongoing efforts to provide charge times that are short and consistent over the life of the product.

## Data Presentation

Charge time data for ICD and CRT-D models are presented using boxplots at 6-month intervals. The shaded box on the plots represents the middle half of the data – the Interquartile Range (IQR). The white line in the middle of each box is the median charge time. The top of the box representing the IQR is the 75th percentile (i.e., 75% of all charge times fall below this line), whereas the bottom of the box represents the 25th percentile. The brackets around the box represent the variance in charge times (analogous to a confidence interval), and are calculated as a function of the median and IQR. Individual values falling outside the brackets are labeled as outliers.


## Results

As shown in the graph below, the performance of Medtronic ICD and CRT-D devices has improved. This graph shows the overall maximums, minimums, and medians for Medtronic ICD and CRT-D products, beginning with the 7221 Micro Jewel. A progression toward shorter mean charge times and less variation has occurred between 1996 and 2002. Models released after 2002 have limited experience but appear to be continuing this performance.




# ICD and CRT-D Charge Time Performance continued


## 7227 GEM Charge Time




## 7229 GEM II VR Charge Time




## 7230 Marquis VR Charge Time




## 7231 GEM III VR Charge Time



## 7232 Maximo VR Charge Time




## 7271 GEM DR Charge Time




# ICD and CRT-D Charge Time Performance continued


## 7272 InSync ICD Charge Time




## 7274 Marquis DR Charge Time




## 7275 GEM III DR Charge Time




## 7276 GEM III AT Charge Time



## 7277 InSync Marquis Charge Time




## 7278 Maximo DR Charge Time




ICD Charge Times

# ICD and CRT-D Charge Time Performance continued


## 7279 InSync III Marquis Charge Time




## 7285 InSync III Protect Charge Time




## 7288 Intrinsic Charge Time




## 7289 InSync II Marquis Charge Time



## 7290 Onyx Charge Time




## 7297 InSync Sentry Charge Time




# ICD and CRT-D Charge Time Performance continued


## 7299 InSync Sentry Charge Time




## 7303 InSync Maximo Charge Time




## 7304 InSync Maximo Charge Time




## C154DWK Concerto Charge Time



## D153ATG, D153DRG EnTrust Charge Time




## D154ATG, D154DRG EnTrust Charge Time




ICD Charge Times

# ICD and CRT-D Charge Time Performance continued


## D154VRC EnTrust Charge Time



## D154AWG Virtuoso Charge Time



## D154VWC Virtuoso Charge Time



## 6930, 6931, 6948, 6949 Sprint Fidelis Defibrillation Leads

Original Date of Advisory: October 2007

### Potential Conductor Wire Fracture

#### Product

All Model 6930, 6931, 6948, and 6949 implantable defibrillation leads

#### Advisory

There are two primary locations where chronic conductor fractures have occurred on Sprint Fidelis leads: 1) the distal portion of the lead, affecting the anode (ring electrode) and 2) near the anchoring sleeve tie-down, predominantly affecting the cathode (helix tip electrode), and occasionally the high voltage conductor. These two locations account for approximately 90% of the chronic fractures identified in Returned Product Analysis (RPA). The remaining 10% of chronic fractures occurred in the DF-1 connector leg and the proximal portion of the RV coil. High voltage conductor fractures could result in the inability to deliver defibrillation therapy. Anode or cathode conductor fractures (at either location) may present clinically as increased impedance, oversensing, increased interval counts, multiple inappropriate shocks, and/or loss of pacing output.

#### Patient Management Recommendations

Medtronic recommends you consider the following as part of routine follow-up for each patient:

- To reduce the risk of inappropriate detection and therapy due to oversensing, program VF detection for initial Number of Intervals to Detect (NID) to nominal settings 8/24) or longer at physician discretion and Redetect NID to nominal settings (12/16)
- Turn ON Patient Alert for RV Pacing, RV Defibrillation, and SVC Defibrillation impedance. For Concerto, Virtuoso, Consulta, Secura and Maximo II devices enrolled on the Medtronic CareLink Network, turn ON the Medtronic CareAlert Notifications for these same parameters.
- To optimize effectiveness of the lead impedance alert:
  - Review V. Pacing Lead Performance Trend to determine typical chronic impedance value for the patient (typical values for Fidelis leads should be 350-1,000 ohms)
  - Program lead impedance alert threshold for RV Pacing to 1,000 ohms, if the typical chronic impedance for the patient is  $\leq 700$  ohms, or
  - Program lead impedance alert threshold for RV Pacing to 1,500 ohms, if the typical chronic impedance for the patient is  $> 700$  ohms
  - Program lead impedance alert threshold for RV Defibrillation and SVC Defibrillation to 100 ohms

#### Status Update

Sprint Fidelis lead performance continues to be in line with the information provided in the October 2007, May 2008 and March 2009 advisory communications. In consultation with the Independent Physician Quality Panel, our patient management recommendations are as follows:

- When a lead fracture is suspected or confirmed, we strongly recommend prompt patient attention. Patients should contact their physician without delay if they experience unexpected shocks.
- The Lead Integrity Alert (LIA)<sup>†</sup> is expected to provide 3 days advance notice prior to inappropriate therapy to 76% of the patients with lead fractures. As a result, we strongly recommend that all Sprint Fidelis patients who have the ability to upgrade to Lead Integrity Alert do so promptly.

- The risk of prophylactic intervention appears to be greater than the risk of serious injury resulting from lead fracture even for pacemaker dependent patients, except in select individual patient circumstances as determined by the physician
- Special circumstances may apply to device change-out or upgrade procedures when a lead fracture has not occurred. At least four options are available, each of which carries risks and benefits that should be taken into consideration:
  - Leave a properly performing lead intact; this is likely to be the best choice for the majority of patients
  - Place a new ICD lead without extraction of the existing lead
  - Place a pace sense lead without the extraction of the existing lead. This option reflects the observation that approximately 90% of Fidelis failures are related to fractures in the pace sense circuit. It is unknown what the failure rate of the high voltage conductor would be should a pace sense conductor failure occur in the existing Sprint Fidelis lead.
  - Unusual patient circumstances may warrant extracting and implanting a new ICD lead. Factors to consider when making this decision include patient life expectancy, age, and comorbidities, number of implanted leads and duration of implant, and patient preference. **Medtronic's Independent Physician Quality Panel recommends that if a lead requires removal, the procedure be performed by a physician with extensive lead extraction experience.** (A new HRS consensus document on lead extraction is expected to be available in May 2009.)

Out of the initial implant population of 204,000 in the United States, approximately 150,100 remain implanted. According to System Longevity Study results, lead survival is estimated to be 93.6% (+2.1/-3.1) at 48 months. As the implanted population ages and the sample size increases for each time interval, the accuracy of the estimated survival probability will increase as shown by tighter confidence intervals.

As part of our commitment to keep you informed about Sprint Fidelis lead performance, Medtronic publishes the quarterly System Longevity Study's all-cause lead survival curve and the CareLink dataset lead survival curve for the Model 6949 lead at [www.medtronic.com/fidelis](http://www.medtronic.com/fidelis). Semi-annual updates will also continue to be provided in the Product Performance Report. Additional information about the Sprint Fidelis lead is available at [www.medtronic.com/fidelis](http://www.medtronic.com/fidelis).

#### Lead Integrity Alert<sup>†</sup>

Medtronic has released Lead Integrity Alert (LIA) software. LIA was designed to provide patients more advance notice via an audible sound of a potential lead fracture that could result in an unnecessary shock.

Data show that with LIA, approximately 76% of the patients with Sprint Fidelis leads are expected to receive 3 or more days advance warning of a potential lead fracture that could result in an unnecessary shock.

Upon hearing the alert, patients should contact their physician without delay.

LIA can be downloaded into nearly all Medtronic implantable cardioverter defibrillators (ICDs) and cardiac resynchronization therapy defibrillators (CRT-Ds) implanted worldwide.



## Sigma Implantable Pulse Generators

Original Date of Advisory: November 2005

### Potential Separation of Interconnect Wires

#### Product

A specific subset of Sigma series pacemakers may fail due to separation of interconnect wires from the hybrid circuit. Specific model and serial numbers of affected devices are available online at <http://SigmaSNList.medtronic.com>.

#### Advisory

This subset of Sigma series pacemakers that may fail due to separation of interconnect wires from the hybrid circuit may present clinically as loss of rate response, premature battery depletion, intermittent or total loss of telemetry, or no output.

Separation of redundant interconnect wires has been observed on hybrid terminal blocks. Device failure occurs only where both interconnect wires separate from a hybrid terminal block. In October 2005, testing and analysis identified the root cause of these failures and the affected population. Hybrid circuits used in this subset of devices were cleaned during manufacturing with a particular cleaning solvent that could potentially reduce the strength of the interconnect wire bond over time.

Our modeling predicts a failure rate from 0.17% to 0.30% over the remaining lifetime of these pacemakers. No provocative testing can predict which devices may fail.

#### Patient Management Recommendations

To assist physicians in their patient care and after discussion with physician consultants, Medtronic offers the following recommendations:

- Medtronic does not recommend replacement of these devices prior to normal elective replacement (ERI), based on the low probability of occurrence of a serious event in this population
- Continue routine follow-up in accordance with standard practice
- Advise patients to seek attention immediately if they experience return of symptoms (e.g., syncope or light-headedness)
- Determine whether device replacement is warranted on a case-by-case basis based upon consultation with patients, review of the individual

patient's medical history, and consideration of the relative risks of an invasive procedure

#### Status Update

**The Sigma Family device performance related to the interconnect wires separation mechanism continues to be within Medtronic's engineering projections.** As of January 31, 2009, 213 devices out of approximately 40,000 devices worldwide have been confirmed as having experienced interconnect wire separation. Fifty-three (53) of these devices were returned from the United States.

One hundred fifty-eight (158) of the 213 devices (0.40%) were returned with information indicating a problem with the patient's pacing system prior to explant. The remaining 55 devices (0.14%) were returned with no information indicating a potential malfunction while implanted or with insufficient information to determine the state of the device at explant. Lacking definite information indicating proper operation until explant, these devices are conservatively categorized as having experienced interconnect wire separation while implanted.

Implant duration for the 213 devices confirmed as having experienced interconnect wire separation has ranged between 17 and 83 months, with an average of 61.3 months.

Out of the initial advisory population of 40,000 worldwide, approximately 16,200 remain implanted. Approximately 3,900 of these are in the United States.

The Patient Management Recommendations set forth in the advisory remain unchanged.

7274 Marquis DR      7278 Maximo DR      7277 InSync Marquis      7279 InSync III Marquis  
 7230 Marquis VR      7232 Maximo VR      7289 InSync II Marquis      7285 InSync III Protect  
 Original Date of Advisory: February 2005

## Potential Premature Battery Depletion Due to Battery Short

### Product

The specific subset of Marquis family ICD and CRT-D devices having batteries manufactured prior to December 2003 is affected. Devices manufactured with batteries produced after December 2003 are not affected. Specific model and serial numbers of affected devices are available online at <http://MarquisSNList.medtronic.com>.

### Advisory

Medtronic Marquis family of ICD and CRT-D devices having batteries manufactured prior to December 2003 may experience rapid battery depletion due to a specific internal battery short mechanism. Battery design changes were implemented in December 2003 that eliminate the possibility of this internal shorting mechanism.

Highly accelerated bench testing indicated the rate of this shorting mechanism may increase as the battery is depleted. As of February 2005, the rate of shorting was approximately 1 in 10,000 (0.01%); bench test data indicated the rate may increase to between 0.2% and 1.5% over the second half of device life.

No provocative testing can predict which of these devices will experience this issue. Once a short occurs, battery depletion can take place within a few hours to a few days. After depletion the device ceases to function. It is also possible that as the battery depletes quickly, patients may experience temporary warmth in the area surrounding the ICD.

### Patient Management Recommendations

We recommend you consider the following patient management options:

- Conduct quarterly (i.e., every 3 months) follow-up procedures
- Inform patients that should they experience warmth in the area surrounding the ICD to seek follow-up care promptly
- Program Low Battery Voltage ERI Patient Alert to “On-High.” This will result in an audible, alternating tone in the limited circumstances where a battery depletes slowly over a number of days. Data indicates most shorts will occur rapidly and will not be detected by this feature.

- Provide a hand-held magnet to patients to check device status and program the Low Battery Voltage ERI Patient Alert to “On-High.” Device operation may be monitored periodically (e.g., daily) by patients placing the magnet over the device for 1-2 seconds. If the device is functional, a steady tone will sound for approximately 20 seconds. If no tone is heard, follow-up care should be sought promptly.

### Status Update

**The Marquis Family device performance related to the battery shorting mechanism continues to be within Medtronic’s engineering projections.** As of January 31, 2009, 129 Marquis Family devices have been confirmed as having this internal battery shorting mechanism. Sixty-six (66) of these devices were returned from the United States.

Of the 129 returns, 40 have been identified by patients reporting warmth in the ICD pocket, 44 by a regularly scheduled follow-up or during a non-device-related hospital visit, 17 by hand-held magnet test or CareLink attempt, nine by return of bradycardia symptoms, five by the Patient Alert sounding, and 14 unknown.

Implant duration for the 129 devices ranged between 11 to 70 months, with an average of 41 months.

Consistent with Medtronic projections, the observed rate of shorting is higher in the second half of device life than in the first half of device life. Of the devices that have exhibited shorting in the last half of device life, 48% occurred in the last quarter of device life and 31% in the last 10% of device life.

Out of the initial advisory population of 87,000 worldwide, approximately 23,900 remain implanted. Approximately 20,700 of these are in the United States.

The Patient Management Recommendations set forth in the advisory remain unchanged.

### Warranty

All Marquis devices are subject to a Limited Warranty. Devices returned to Medtronic and determined to be malfunctioning will be replaced or credited in accordance with the warranty terms. In addition, for devices subject to the Marquis advisory, should a physician decide to replace a device for a patient who is pacemaker dependent or who receives frequent VT/VF therapy, Medtronic will, upon receipt of a written statement from the physician setting forth the basis for early replacement of the device, provide a replacement device at no cost.

## Kappa 600, 700 Dual Chamber (D, DR, and VDD) IPGs

Original Date of Advisory: March 15, 2002

### Potential Fractured Power Supply Wires

#### Product

A specific subset of Kappa 700/600 dual chamber (D, DR, and VDD) implantable pulse generators has been identified by serial numbers. Hospitals and Physicians were notified. Specific model and serial numbers of affected devices are available by calling US Technical Services at 1 (800) 505-4636.

#### Advisory

As of March 15, 2002, Medtronic observed 53 related failures (0.02%) in over 255,000 Kappa 700/600 dual chamber (D, DR, and VDD) series devices sold worldwide. Medtronic voluntarily communicated this information to physicians because these failures had been observed in patients having submuscular implants.

These devices have presented with an electrical reset, intermittent output, or no output. Our investigation identified the root cause as fractured wires supplying power to the pacemaker. This has been directly correlated to submuscular placement of these devices. Submuscular implant locations (e.g., subpectoral, abdominal, etc.) can result in additional stress and repetitive flexing on the implanted device causing excessive fatigue on these wires. Of the estimated 4,000 devices implanted submuscular, approximately 200 (5%) may experience this failure. These stresses on the implanted device are unique to submuscular implant sites and do not exist with subcutaneous implants.

#### Patient Management Recommendations

While there is no provocative testing or time dependency that will predict which submuscular placed device will fail, certain electrical resets may be an indicator that a wire fracture has occurred. Normal electrical resets can occur as a result of electro-surgical procedures such as cautery and ablation or from defibrillation therapy. If none of the normal causes of electrical reset can be confirmed, or if a device serial number presents as "000000" following an electrical reset, this may be an indicator of a wire fracture.

For patients who have submuscular implants of devices within the designated serial number range and who are pacemaker dependent with no underlying rhythm, replacement of the device should be considered. Medtronic will provide the replacement device free of charge under the terms of its warranty program if a device is replaced in these patients.

For patients having subcutaneous implants, no change to your current patient care and follow-up is advised.

#### Status

**Device performance related to this advisory continues to be within Medtronic's engineering projections. Patient management recommendations remain unchanged.** As of January 31, 2009, 308 out of approximately 180,000 distributed (0.17% incidence) Kappa family devices worldwide have been confirmed as having fractured power supply wires. One hundred sixty-three (163) of these devices were returned from the United States. Out of the initial implant population of 121,000 in the United States, approximately 24,700 remain implanted.

## 7227Cx GEM 7229Cx GEM II VR

Original Date of Advisory: October 15, 1999

### Potential Circuit Overload

#### Product

Model 7227Cx and Model 7229Cx implantable cardioverter defibrillators supplied before October 15, 1999, with serial numbers ending in an “H.” For example, PIPxxxxxxH or PJJxxxxxxH, where x is a variable numeric, may be affected. Specific model and serial numbers of affected devices are available by calling US Technical Services at 1 (800) 723-4636.

#### IMPORTANT REMINDER:

Medtronic strongly advises physicians who have patients under their care affected by this issue to reprogram the Patient Alert feature “ON” without delay.

#### Advisory

Manufacturing error in a small percentage of devices may cause circuit overload when  $AX \geq B$  High Voltage energy is delivered via an integrated bipolar lead. GEM Model 7227Cx and GEM II VR Model 7229Cx devices with dedicated bipolar sensing leads are not affected by this issue. Devices affected may not be able to subsequently charge to full energy and experience “charge circuit timeout.”

#### Patient Management Recommendations

- Assessment of all patients with the potentially affected devices implanted **AND** an integrated bipolar ICD lead such as the Models 6942 and 6945 should take place without delay
- Reprogram polarity pathway to  $B \geq AX$  for all cardioversion and defibrillation therapies
- Confirm correct device function:
  - Perform a full energy charging sequence
  - If “charge circuit timeout” is observed, contact your Medtronic representative
  - If device charges normally, it has not been damaged and will function appropriately with polarity programmed  $B \geq AX$

Recent studies have demonstrated that DFTs are similar or lower in a  $B \geq AX$  polarity pathway when compared to  $AX \geq B$ .

Devices implanted with functional dedicated bipolar leads such as the Models 6932, 6934S, 6936, 6943, and 6966 are not affected.

#### Status

**Device performance related to this advisory continues to be within Medtronic’s engineering projections. Patient management recommendations remain unchanged.** Out of the initial implant population of 10,000 in the United States, approximately 1,200 remain implanted. The devices affected by this advisory are nearing the end of their expected battery longevity.

## 4504, 4504M CapSure Atrial Lead

## 4582 Target Tip Atrial Lead

Original Date of Advisory: October 4, 1996

### Lead Survival Below Expectations

#### Product

All Model 4504, 4504M, and 4582 implantable pacing leads

#### Advisory

Lead survival probability is below expectations and is primarily associated with insulation degradation due to Metal Ion Oxidation (MIO).

#### Patient Management Recommendations

- Follow patients in accordance with Medicare Guidelines
- Avoid the use of the AAI or AOO mode
- During patient evaluation, give careful attention to lead performance such as:
  - Review patient ECG for indications of transient sensing and/or capture abnormalities
  - Monitor in clinic for impedance less than 250 ohms or a decrease of more than 30% from implant values (or an established baseline using telemetry), which would suggest lead failure
- Consider the use of unipolar if the pulse generator has this capability
- At the time of pacemaker system revision (e.g., normal pulse generator or ventricular lead revision), carefully evaluate lead integrity and patient status before choosing to reuse

#### Status

**Patient management recommendations remain unchanged. Laboratory analysis trends and engineering conclusions remain unchanged.**

Out of the initial implant population of 16,600 in the United States, approximately 1,700 remain implanted. According to System Longevity Study results, lead survival is estimated to be 66.1% at 8 years, 9 months.

## 4004, 4004M CapSure Ventricular Lead

## 4082 Target Tip Ventricular Lead

Original Date of Advisory: October 8, 1993

### Lead Survival Below Expectations

#### Product

All Model 4004/4004M and 4082 implantable pacing leads

#### Advisory

Lead survival probability is below expectations due primarily to polyurethane insulation failure (MIO) and conductor fracture (associated with “subclavian crush”).

#### Patient Management Recommendations

- Increase, as appropriate, the frequency of patient evaluation through in-clinic visits supplemented with transtelephonic and/or ambulatory monitoring; for example, consistent with Guideline I under Medicare Pacemaker Monitoring Guidelines (50-1 Cardiac Pacemaker Evaluation Services)
- During patient evaluations, give careful attention to lead performance such as:
  - Reviewing patient ECGs carefully for indications of transient sensing and/or capture abnormalities
  - Monitoring in-clinic for impedances less than 300 ohms or a decrease of more than 30% from implant values (or an established baseline using telemetry), which would suggest lead failure
  - Eliciting and thoroughly investigating any patient complaints suggestive of lead failure
- Consider whether prophylactic replacement would be appropriate, especially in patients at high risk, such as pacemaker dependent patients
- Carefully evaluate lead integrity when performing routine pulse generator replacements. Replace lead if:
  - Insulation breaches are observed
  - Lead impedance is less than 300 ohms or has decreased by more than 30% from implant values
  - Impedance or voltage threshold measurements vary significantly when multiple readings are taken
  - If the risk of continued use outweighs the risk associated with implanting a new lead
- As always, individual circumstances and medical judgment dictate patient care and frequency of follow-up
- Consider lead replacement during normal pulse generator change-out. Carefully evaluate lead integrity and patient status before choosing to reuse.

#### Status

**Patient management recommendations remain unchanged. Laboratory analysis trends and engineering conclusions remain unchanged.**

Out of the initial implant population of 77,000 in the United States, approximately 6,000 remain implanted. According to System Longevity Study results, lead survival is 50.6% at 10 years, 9 months.

## 4012 Target Tip Ventricular Lead

Original Date of Advisory: September 26, 1991

### Lead Survival Below Expectations

#### Product

All Model 4012 implantable pacing leads

#### Advisory

Lead survival probability beyond 5 years is below expectations due primarily to polyurethane insulation failure (due to ESC and/or MIO) and conductor fracture (associated with “subclavian crush”).

#### Patient Management Recommendations

Consider increasing frequency of monitoring (e.g., from quarterly to bimonthly or monthly). Consider the following activities as part of normal follow-up procedures:

- Monitor for significant changes in impedance, which could be an indication of impending failure (pulse generator must have impedance telemetry capabilities)
- Review patient ECGs carefully for indications of transient sensing and/or capture abnormalities. This can be done using transtelephonic or in-clinic monitoring and/or using ambulatory monitoring.
- Elicit any patient complaints suggestive of lead failure and investigate thoroughly lead integrity/performance characteristics following reports of patient complaints or symptoms using the above techniques

- Consider whether prophylactic replacement would be appropriate in patients at high risk, such as pacemaker dependent patients
- Evaluate carefully the integrity of the lead during routine pulse generator replacement before choosing to reuse. Specifically, Medtronic recommends placement of a new lead if:
  - Insulation breaches are observed
  - Lead impedance is less than 300 ohms or has decreased by more than 30% from implant values
  - Electrical properties such as impedance and threshold vary significantly when multiple readings are taken

As always, medical judgment must be used to establish the appropriate schedule and course of care for every individual, particularly pacemaker dependent or other patients at higher risk.

#### Status

**Patient management recommendations remain unchanged. Laboratory analysis trends and engineering conclusions remain unchanged.**

Out of the initial implant population of 96,800 in the United States, approximately 6,500 remain implanted. The System Longevity Study results show 62.2% lead survival at 15 years, 9 months.

## Minix, Minix ST, Micro Minix IPGs

Original Date of Advisory: May 6, 1991

### Potential Delayed Restoration of Permanent Settings

#### Product

All Models of the Minix, Minix ST, and Micro Minix families of implantable pulse generators

#### Advisory

Possibility of delayed restoration of permanent pacing mode and parameters, after the magnet or programming head is removed under certain conditions.

#### Patient Management Recommendations

To eliminate any potential risk associated with temporary programming, depress the INTERROGATE key and verify successful interrogation before moving the programming head away from the pulse generator. The delay condition can also be terminated by repositioning the programming head and depressing the EMERGENCY VVI key.

#### Status

**Device performance related to this advisory continues to be within Medtronic's engineering projections. Patient management recommendations remain unchanged.** Out of the initial implant population of 65,000 in the United States, approximately 4,100 remain implanted. The devices affected by this advisory are nearing the end of their expected longevity.



# Performance Notes

## Clinical Management of VCM near Elective Replacement

### Background

Medtronic Technical Services has received reports of devices going to ERI or end of life (EOL) sooner than expected after a normal follow-up in which the device longevity was projected to be approximately 18 months. It has been noted that these cases typically involve Kappa 700 devices where Ventricular Capture Management set the ventricular lead to high output (5 V, 1 ms), which occurs by device design when a high threshold is measured. It is important for physicians and allied professionals to understand VCM behavior as it relates to longevity so that they can, in turn, understand how this affects management of the device and follow-up visits as VCM equipped IPGs near the end of their expected longevity.

### Device Longevity and VCM Behavior

Ventricular Capture Management is a feature that uses evoked response sensing to determine the stimulation threshold needed to capture the ventricular chamber. Proper detection of the evoked response is crucial to the VCM algorithm determining an accurate capture threshold. There are rare conditions, however, during which the VCM algorithm will not be able to measure the evoked response accurately.<sup>1</sup> When this occurs, for safety reasons the VCM algorithm will reprogram the output to 5 V, 1 ms until the subsequent VCM measurement.

If the device has considerable remaining longevity, these occasional excursions to high output do not substantially affect remaining longevity. However, if the device has less than approximately 18 months remaining longevity, there is the possibility that the high output condition caused by the 5 V, 1 ms output will drain the battery and trigger ERI.

When ERI is declared by the device, VCM is disabled and the outputs are left at 5 V, 1 ms until the device is reprogrammed at an in-office follow-up. This increased current drain of a high output condition will speed depletion of the device, possibly resulting in the device getting to the EOL (battery voltage  $\leq 2.15$  V).

Please note that the following parameter changes occur when the device goes to ERI:

Table: IPG Therapy Parameter Changes at ERI

Parameter	Value
Pacing Mode	VVI
Lower Rate	65 bpm
Single Chamber Hysteresis	OFF
Sleep Function	OFF
Ventricular Capture Management	OFF
Atrial Sensing Assurance	OFF
Ventricular Sensing Assurance	OFF

Kappa 700 is Medtronic's first-generation VCM algorithm, which has a relatively higher incidence of evoked response undersensing compared to subsequent algorithms, resulting in more frequent high output conditions. Therefore, Kappa 700 products are the primary focus of this note. It should be noted that IPGs equipped with the second-generation VCM algorithm (Kappa 900, EnPulse, Adapta/Versa/Sensia, and Relia) have not been observed with evoked response undersensing in the general population, though the items listed in "Follow-Up Considerations" may also be used on these devices.

### Follow-Up Considerations

- Estimated longevity in the event the device goes to high output can be determined by the following steps. This allows the clinician to determine follow-up frequency if he or she is concerned the device may go to ERI due to high output.
  - Program the ventricular channel to 5 V, 1 ms
  - Navigate to Data/Battery and Lead Measurements
  - When the message stating "Warning – Old Data" is displayed, select "Yes" to measure battery voltage and lead impedance at the new ventricular outputs
  - An updated remaining longevity estimate will be calculated on the elevated outputs. Note the "Minimum Remaining Longevity." Clinical decisions can be based on this value.
  - Program the Amplitude and Pulse Widths back to their original values before leaving the session
- If the capture trends and lead impedance trends are stable, VCM can be programmed to "Monitor Only" for the remaining device life. This should be considered only if remaining longevity is 18 months or less.
- Follow-up frequency can be increased for those patients who do not have stable capture or lead impedance trends. This can be done via a CareLink Home Monitor, or in-office.

<sup>1</sup> Medtronic, Inc. (2001). Medtronic Kappa 700/600 Series Pacemaker Reference Guide (Chapter 4, p. 27). Can be retrieved from <http://manuals.medtronic.com>.

## Ensuring the Accuracy of Battery Longevity Estimates

### Purpose of This Information

This article is intended to help the clinician understand how Medtronic estimates CRT-D, ICD, and IPG device longevity and Medtronic’s performance against these estimates.

### Device Longevity and Battery Depletion

The device service life ends when the usable battery capacity is depleted. The time to battery depletion depends on three factors:

- The amount of electrical energy expended in providing therapy to the patient
- The amount of energy consumed by the electronic circuitry to perform the functions of the device (e.g., operating the microprocessor, telemetry, memory, and charging component)
- The energy capacity of the battery

Medtronic has developed a statistical model for device longevity that accounts for each of these factors, and has validated the model with real time clinical performance. During the development of its products, Medtronic engineers characterize device longevity using this model. Testing begins during development and continues after market release to ensure the accuracy of device longevity estimates.


### Using Survival Curves to Assess Longevity

The survival curves in the Product Performance Report represent the composite experience of thousands of devices over a wide range of programming options and patient use conditions. While the curves are useful for understanding the overall performance of a population of devices, they cannot be used to accurately predict the longevity of a specific device in a specific patient. To get a longevity prediction for a specific device, the longevity model must be used. The model is available by contacting Medtronic’s Technical Services Department.

Because the survival curves are an aggregate result, the Reference pages in the Product Performance Report include several longevity estimates for a range of use conditions. These longevity estimates, originally published in the device Technical Manual, are mean values calculated for the parameters given. This range of longevity estimates can be compared to the survival curve including normal battery depletion to assess the overall clinical performance of a device model against the original longevity estimates.

If most of a device model’s population is being used at nominal parameters and conditions, the time at which the survival curve including normal battery depletion equals 50% should approximate the midpoint in the range of longevity estimates.

If devices tend to be used at conditions that consume more or less energy than nominal, then the time at which the survival curve equals 50% should tend toward the lower or higher end of the range of longevity estimates, respectively.



## Interactions between Cardiac Pacing and Ventricular Arrhythmia Initiation

### Purpose of This Information

This article is intended to provide information for consideration when programming pacemaker operation in ICDs and pacemakers.

### Background

Right ventricular pacing has been associated with increased risk of appropriate therapy for ventricular tachycardia (VT) and ventricular fibrillation (VF) in ICD patients.<sup>1</sup> Abrupt changes in ventricular cycle lengths (short-long-short, S-L-S) may precede initiation of VT/VF in some instances. S-L-S sequences may be permitted in all forms of cardiac pacing. The pause lengths depend upon pacing mode and lower rate programming.<sup>2-4</sup> Because pauses may be associated with VT/VF initiation, pause suppression algorithms have been developed in ICDs. Although pause suppression may have utility in specific patients with repolarization abnormalities and pause dependent VT, it has not been shown to reduce arrhythmia incidence in the general ICD population.<sup>5</sup> Conversely, S-L-S sequences may occur with ventricular pacing in a variety of ways, including atrial tracking of premature atrial contractions (PACs) or by terminating pauses with ventricular paced beats.<sup>6</sup> In some patients, the ectopic depolarization pattern of a ventricular paced beat may be pro-arrhythmic, independent of pause timing. These observations have further enforced the desire to reduce unnecessary ventricular pacing.

### Clinical Trial Observations

Medtronic-sponsored clinical trials were retrospectively analyzed to further understand pause-mediated (i.e., S-L-S) scenarios prior to VT/VF. S-L-S onset scenarios were observed in a minority of patients in all pacing modes. Pacemaker interactions prior to VT/VF are dependent on patient conditions, as well as the technical aspects of pacing operation (i.e., pacing mode, lower rate, and AV interval). Because a very low frequency of ventricular pacing is observed during Managed Ventricular Pacing (MVP)<sup>7-9</sup> or VVI 40 pacing modes,<sup>10</sup> the long interval tended to terminate with a ventricular sense. In DDD mode, the long interval tended to be terminated by a ventricular pace. Long intervals of > 1,000 ms prior to VT/VF were rare in MVP mode. In these analyses, only an association between cardiac pacing and VT/VF initiation can be observed, causality cannot be established. The ongoing MVP (Managed Ventricular Pacing vs. VVI 40 Pacing) Trial, a 2-year, 1,000-patient prospective, randomized trial in ICD patients may offer more insight into the frequency of VT/VF across pacing modes.<sup>11</sup>

### Pacemaker Patients

In pacemaker patients, ventricular pacing has been associated with higher incidence of AT/AF and heart failure hospitalization.<sup>12,13</sup> MVP provides atrial rate support while dramatically reducing ventricular pacing in patients with sinus node dysfunction and transient AV block.<sup>9</sup> However, as stated in Medtronic reference manuals, depending upon the patient's intrinsic rhythm and conduction, MVP may allow ventricular cycle variation and occasional pauses of up to twice the lower rate.

DDD pacing with long AV intervals may reduce ventricular pacing and may decrease the potential length of pauses compared to MVP. However, DDD with long AV interval programming does not appear to be as effective as AAI-based pacing modes at reducing ventricular pacing,<sup>13,14</sup> may lead to endless loop tachycardia,<sup>14,15</sup> and does not completely eliminate pauses. Also, in DDD mode, a higher programmed lower rate or activation of rate response can lead to an increase in AV conduction times and a higher percentage of ventricular pacing. The potential benefits of reducing ventricular pacing must be weighed against the potential for longer ventricular pauses. Therefore, careful consideration should be given to pacemaker mode and lower rate programming, particularly in the setting of frequent AV block and repolarization abnormalities due to congenital Long QT, electrolyte imbalances, and some medications that prolong QT.

### Conclusion

Pacemaker operation may interact with VT/VF initiation in a variety of ways. The patient's heart failure status, arrhythmia substrate, medications, and the relative importance of maintaining ventricular synchrony versus ensuring ventricular rate support must be weighed when choosing optimal hardware (ICD vs. pacemaker) and pacemaker programming (pacing mode, lower rate, etc.).

### References

- Steinberg JS, Fischer A, Wang P, et al. The clinical implications of cumulative right ventricular pacing in the multicenter automatic defibrillator trial II. *J Cardiovasc Electrophysiol*. April 2005;16(4):359-365.
- Pinski SL, Eguia LE, Trohman RG. What is the minimal pacing rate that prevents torsades de pointes? Insights from patients with permanent pacemakers. *PACE*. November 2002; 25(11):1612-1615.
- Goldman DS, Levine PA. Pacemaker-mediated polymorphic ventricular tachycardia. *PACE*. October 1998; 21(10):1993-1995.
- Gray CJ, Basta M, Sapp JL, Parkash R, Gardner MJ. Inappropriate application of managed ventricular pacing in a patient with Brugada syndrome leading to polymorphic ventricular tachycardia, ventricular fibrillation and implantable cardioverter defibrillator shocks. *Heart Rhythm*. 2006, Abstract P1-89.
- Friedman PA, Jalal S, Kaufman S, et al. Effects of a rate smoothing algorithm for prevention of ventricular arrhythmias: results of the Ventricular Arrhythmia Suppression Trial (VAST). *Heart Rhythm*. May 2006;3(5):573-580.
- Himmrich E, Przibille O, Zellerhoff C, et al. Proarrhythmic effect of pacemaker stimulation in patients with implanted cardioverter-defibrillators. *Circulation*. July 15, 2003;108(2):192-197.
- Sweeney MO, Ellenbogen KA, Casavant D, et al. Multicenter, prospective, randomized trial of a new atrial-based Managed Ventricular Pacing Mode (MVP) in dual chamber ICDs. *J Cardiovasc Electrophysiol*. 2005;16:1-7.
- Sweeney MO, Shea JB, Fox V, et al. Randomized pilot study of a new atrial-based minimal ventricular pacing mode in dual-chamber implantable cardioverter-defibrillators. *Heart Rhythm*. July 2004;1(2):160-167.
- Gillis AM, Pürerfellner H, Israel CW, et al. Reducing unnecessary right ventricular pacing with the managed ventricular pacing mode in patients with sinus node disease and AV block. *PACE*. July 2006; 29(7):697-705.
- Dual-chamber pacing or ventricular backup pacing in patients with an implantable defibrillator: the Dual Chamber and VVI Implantable Defibrillator (DAVID) Trial. *JAMA*. December 25, 2002;288(24):3115-3123.
- Sweeney MO, Ellenbogen KA, Miller EH, Serfese L, Sheldon T, Whellan D. The Managed ventricular pacing versus VVI 40 Pacing (MVP) Trial: clinical background, rationale, design, and implementation. *J Cardiovasc Electrophysiol*. December 2006;17(12):1295-1298.
- Sweeney MO, Ellenbogen KA, et al. Adverse effect of ventricular pacing on heart failure and atrial fibrillation among patients with normal baseline QRS duration in a clinical trial of pacemaker therapy for sinus node dysfunction. *Circulation*. June 17, 2003;107(23):2932-2937.
- Nielsen JC, Kristensen L, Andersen HR, Mortensen PT, Pedersen OL, Pedersen AK. A randomized comparison of atrial and dual-chamber pacing in 177 consecutive patients with sick sinus syndrome: echocardiographic and clinical outcome. *J Am Coll Cardiol*. August 20, 2003;42(4):614-623.
- Nielsen JC, Pedersen AK, Mortensen PT, Andersen HR. Programming a fixed long atrioventricular delay is not effective in preventing ventricular pacing in patients with sick sinus syndrome. *Europace*. April 1999; 1(2):113-120.
- Dennis MJ, Sparks PB. Pacemaker mediated tachycardia as a complication of the autointrinsic conduction search function. *PACE*. June 2004;27(6 Pt 1):824-826.

## AT500 Pacing System Follow-Up Protocol

### Purpose of This Information

This article is intended to provide clinical guidance regarding follow-up practice and patient management when the AT500 battery voltage approaches the Elective Replacement Indicator (ERI) level of 2.6 volts.

### Background

Many AT500 pacing systems are now reaching their ERI voltage level (2.6 volts). This is expected since the battery used has an approximate longevity of 5-6 years under normal conditions (100% DDD pacing, 3 volts, 0.4 ms).

Technical Services has received reports of battery voltage levels below end of life (EOL of 2.2 volts) where EGM pre-storage is programmed ON, or higher outputs and/or pacing rates are necessary. It is important for physicians and allied professionals to understand battery depletion characteristics between ERI and EOL so that they, in turn, can understand how this affects management of follow-up visits for the AT500 as this device nears the end of its expected longevity.

### AT500 Battery and Longevity Information

In contrast to other IPGs, the AT500 does not change its mode, stimulation rate, or any other parameter when the battery voltage drops below the ERI level of 2.6 volts (with or without magnet applied). The Threshold Margin Test (TMT) is also not available.

Therefore, it is not possible to perform transtelephonic assessment of AT500 battery status. This must be done during an in-clinic follow-up session. A warning will be displayed on the Quick Look screen at the beginning of a programmer (follow-up) session when the ERI battery level

occurs. The measured battery voltage will also appear on the programmer display and on printouts.


Battery depletion curves are shown in Figure 1, with special focus on device longevity when programming EGM pre-storage ON or OFF.

Medtronic's review of ongoing AT500 battery life test data matches our original longevity modeling and so meets our expectations. However, when using longer durations between follow-up periods (> 3 months), clinicians should consider the following in setting their remaining longevity expectations.

- Enabling the "EGM Pre-storage On" capability will increase current and reduce device longevity by approximately 9 days for each month pre-storage is ON
- Longevity decreases with an increase in pacing rate, an increase in pacing amplitude or pulse width, a decrease in pacing impedance, a higher ratio of paced to sensed events, or extended use of the Atrial Preference Pacing, EGM pre-storage, or Holter Telemetry features

### Recommendations

Follow-up frequency should always be accelerated as devices reach ERI voltage levels to ensure device explant/replacement occurs prior to end of life voltage levels. With the wide variety of follow-up schedules being used, Medtronic recommends a 3-month follow-up frequency for the AT500 pacing systems. This is particularly important for patients in whom EGM pre-storage is programmed ON, or higher outputs and/or pacing rates are necessary.



**Figure 1**  
AT500 battery depletion curve for common parameter settings of DDDR, LR 70 ppm, UR 120 ppm, 100% pacing, Atrial – 2 V, 0.4 ms, 600 ohms, Ventricle – 2 V, 0.6 ms, 900 ohms, and EGM Pre-storage ON versus OFF.

## Insertion of the Lead into the Device

The implantable system consists of a pulse generator and at least one lead. The system operation depends on proper electrical and mechanical operation. With the advent of internationally recognized connector standards, the challenge of ensuring proper mechanical fit between the lead and device connectors has been simplified, although the international connector standard does not address all aspects of the procedure for connecting a lead to the device.

If the lead connector is not fully installed, oversensing may result as described in the connector problems section of the technical article, “Clinical Management of High Voltage Lead System Oversensing.”

Performing the following steps can be used for each lead connection during the implant procedure:



**1** Insert the torque wrench into the appropriate setscrew. For easier lead insertion, insert the lead closest to the device first.

- 2** Look down the connector port to verify that the port is not obstructed. If the port is obstructed, retract the setscrew to clear the bore. Take care not to disengage the setscrew from the connector block.
- 3** Push the lead into the connector port until the lead pin is clearly visible beyond the setscrew block.
- 4** Hold the lead in position while tightening the setscrew until the torque wrench clicks.
- 5** Tug gently on the lead to confirm a secure fit.


Current publications may provide additional information on implant procedures used by others, e.g., radiographic evaluation of the terminal pin beyond the terminal post.<sup>1</sup>

<sup>1</sup> Pickett RA III, Saavedra P, Ali MF, Darbar D, Rottman JN. Implantable cardioverter-defibrillator malfunction due to mechanical failure of the header connection. *J Cardiovasc Electrophysiol*. September 2004;15(9):1095-1099.

**Connector module before lead insertion**



**Cross-section of connector module after lead fully installed**



**X-ray image of connector module after lead fully installed**



## GEM II DR/VR and GEM III DR/VR/AT ICD Battery Discharge Behavior

Medtronic manufactures and utilizes a unique lithium/silver vanadium oxide battery in the GEM II/III family of ICDs. This battery has a distinctive voltage discharge with two regions of constant voltage at 3.2 volts and 2.6 volts.

The battery discharge curve (see curve below) is characterized by a significant decrease in the battery voltage approaching middle of life (MOL), followed by a plateau (MOL to ERI) where the battery voltage remains around 2.6 volts. The transition to the plateau could be easily misinterpreted as the battery rapidly approaches ERI, which occurs at 2.55 volts, when the battery may in fact have several years remaining until ERI.

It is important to understand that this battery voltage decrease in the GEM II/III family of ICDs is a normal


characteristic of the battery function in these devices and should not create a need for additional follow-up or monitoring.

As a general rule of thumb, the longevity from implant to MOL = MOL to ERI.

The design of the battery in subsequently released models has been modified to present a more linear battery discharge curve.

If you are concerned about early ERI in your patient's device, you can utilize the battery trend measurements stored in the save-to-disk file, which can be accessed and interpreted through the Medtronic Technical Services at 1 (800) 723-4636.

**GEM II/III Battery Discharge Curve**



## General Follow-Up and Replacement of ICD Leads

Implanted leads operate in the challenging biochemical environment of the human body and the body's response to foreign objects. Implanted leads are also subject to mechanical stresses associated with heart motion, body motion, and patient anatomy.

In this environment, pacemaker and defibrillation leads cannot be expected to last forever. Unlike implantable cardioverter defibrillators (ICDs), a lead's longevity cannot be predicted nor are there simple indicators that a lead is approaching the end of its service life. The determination that a lead may be approaching end of service life requires follow-up of the chronically implanted lead and thorough evaluation of lead integrity at ICD replacement.

### Follow-Up of Chronically Implanted Leads

The frequency of follow-up for ICD patients will depend on a number of factors including the patient's medical condition, ICD system implant time, hospital/clinic follow-up practice, and Medicare guidelines. In all cases, it is important to assess the functionality of the ICD system and the integrity. For newly implanted leads, it is beneficial to establish a baseline of chronic performance parameters once the lead has stabilized, generally within 6 to 12 months after implant. These performance parameters should include pacing and sensing thresholds and impedance. During routine patient follow-up, these procedures can be used to evaluate lead integrity.

- Measure pacing and sensing threshold and compare to the chronic baseline. Significant increases or decreases may be indicative of lead failure, dislodgement, perforation, exit block, etc.
- Measure pacing impedance where possible and compare to the chronic baseline. Decreases of 30% or more or pacing impedances below 200-250 ohms may be indicative of insulation failure. Sudden and significant increases in pacing impedance may be indicative of conductor fracture.
- High voltage lead circuit impedance should be between 10-75 ohms at system implant. Chronic measurements below 10 and above 200 ohms may be indicative of high voltage lead circuit failure.
- Carefully review ECGs or the nonsustained detection log on Medtronic ICDs for indications of pacing and/or sensing abnormalities such as oversensing, undersensing, and loss of capture
- Elicit and investigate any patient complaints/symptoms that may be suggestive of potential lead failure

Where routine follow-up indicates, additional tools should be used to further evaluate performance. Tools include radiographic data, ICD electrograms, ICD Patient Alert and performance information from the System Longevity Study (SLS).

The final decision on the functional integrity and continued use of an implanted lead must be a matter of medical judgment based on these factors as well as specific patient conditions.

### General Criteria for Lead Replacement

The evaluation of a chronically implanted lead is an important part of the decision to continue to use the lead with a new ICD. However, these results alone do not necessarily predict the future integrity of that lead. With the expected longevity of today's ICDs varying between approximately 5 and 10 years, a physician replacing a device should consider a number of factors, including those listed below.

Factors that should be considered in a decision to replace or continue to use include:

- Pacing and sensing thresholds should be evaluated for the potential to maintain acceptable levels
- Pacing impedance should be measured. Bear in mind that pacing impedance below 250 ohms results in excessive battery current drain, which may seriously compromise ICD longevity, regardless of lead integrity.
- The physical appearance of the lead should be examined for insulation cracks, breaches, or other indications of lead wear or degradation
- Medtronic System Longevity Study data should be referenced. Actuarial survival of the lead and the observed lead failure mechanisms are specific factors to consider. Use of a new lead should be considered if failure mechanisms suggest an increased time dependency as suggested in the shape of performance curve for the specific lead model.
- Current publications may provide additional information on the clinical management of leads.<sup>1-3</sup> Ultimately, the decision to replace an implanted lead involves medical judgment.

<sup>1</sup> Hauser RG, Cannom D, Hayes DL, et al. Long-term structural failure of coaxial polyurethane implantable cardioverter defibrillator leads. *PACE*. June 2002;25(6):879-882.

<sup>2</sup> Ellenbogen KA, Wood MA, Shepard RK, et al. Detection and management of an implantable cardioverter defibrillator lead failure: incidence and clinical implications. *J Am Coll Cardiol*. January 1, 2003;41(1):73-80.

<sup>3</sup> Hauser RG, Kallinen LM, Almquist AK, Gornick CC, Katsiyannis WT. Early failure of a small-diameter high-voltage implantable cardioverter-defibrillator lead. *Heart Rhythm*. July 2007;4(7):892-896.

## Clinical Management of High-Voltage Lead System Oversensing

Appropriate sensing by an ICD system refers to the sensing of cardiac events that may or may not require therapy delivery. ICD systems must sense relatively large QRS complexes while avoiding sensing of smaller T waves, yet continue to sense often small variable amplitude ventricular fibrillation. Thus, ICD systems attempt to dynamically adjust sensing of electrical events and discriminate between them based on detection algorithms and programmed settings.

Inappropriate sensing can occur when an ICD system classifies events of non-cardiac origin as QRS/VF events, or senses and counts T and far-field P waves as ventricular depolarizations. This is often referred to as “oversensing,” and may result in delivery of inappropriate high-voltage therapies. This is due, in part, to the desire to err on the side of delivering lifesaving high voltage therapy rather than withholding

it. Thus, an ICD system that is experiencing oversensing issues will continue to deliver therapeutic shocks as required, but may also subject the patient to unnecessary shocks.

Oversensing can be difficult to manage, in that the precipitating cause of the oversensing can be problematic to isolate. Oversensing can be caused by many factors, including myopotentials/far-field sensing, electromagnetic interference, T wave sensing, connector issues, incomplete or complete conductor fractures, and insulation breaches. While the individual physician must exercise medical judgment in determination of appropriate clinical management of ICD systems, the chart below may assist in the process of causal factor differentiation and possible intervention.

Phenomenon	Causal Factors	Characteristics	Management/Comments
Myopotentials/ Far-field sensing	Diaphragmatic muscle potentials in breathing, wide tip-to-ring (coil on integrated bipolar leads) spacing	Nonphysiological sensed event on EGM, which may confuse detection potentially resulting in false positive shocks	Check R waves for deterioration. Reprogram sensitivity. Try repositioning lead. Consider change-out to true bipolar lead, or if true bipolar lead in use, one with closer tip-to-ring spacing than current lead.
EMI (Electro-Magnetic Interference)	Arc welders, electrical generators, store walk-through security scanners, poorly insulated electrical equipment	Multiple and consecutive short intervals (< 140 ms) independent of underlying sinus beats. Associated with proximity to the EMI source.	Avoid EMI areas. True bipolar leads less susceptible.
T-wave sensing	Drugs, ischemic tissue, exercise, Long QT syndrome, electrolyte imbalance	Sense markers seen on EGM related to T wave. False positive detection.	Check for R wave deterioration and characteristics. If R wave > 3.0 mV, reprogram sensitivity. If R wave < 3.0 mV, reposition/replace lead. Address causal factor (e.g., drugs [if appropriate/medically viable]).
Connector problems	Loose setscrew, cross-threaded setscrew, incomplete lead insertion into header	This is an acute phenomenon seen within 6 months of implant (usually sooner)	Requires invasive check of connections. May be reproducible with pocket manipulation.
Incomplete conductor fracture	One or more filars of a multifilar conductor fracturing while leaving enough filars intact to provide a conduction circuit	Characterized by chaotic oversensing related to motion of the fracture site	Check EGMs and x-rays. Manipulate lead at suspected fracture site if possible as a provocative test. If confirmed, replace lead.
Lead insulation breach	Cuts, tears, metal ion oxidization, abrasion, cold flow, environmental stress cracking	Characterized by cyclical and/or erratic, intermittent, spontaneous oversensing; often post-pace or post-shock can cause false positives	Replace lead. If acute, usually secondary to implant damage/replacement damage. If late, material characteristic.
Oversensing during interrogation with programming head (not wireless telemetry) with complete lead fracture	Interrogation with a programming head in combination with complete lead fracture that creates an open circuit can induce noise on the sensing circuitry inside the ICD can	Nonphysiologic sensed event on EGM. If detection is enabled during interrogation, oversensing may result in inappropriate therapy.	Quickly remove the programming head. CANCEL the interrupted interrogation and manually load the software for the specific device model. Reposition the programmer head over the device and immediately select SUSPEND. Device will resume detection when programming head is removed, or when RESUME is selected. Replace lead.

Technical Services is available at all times to advise clinicians in the troubleshooting and management of Medtronic products. For assistance in the United States, please call 1 (800) 723-4636. In other countries, please contact your local Medtronic representative.



## Tests and Observations for Clinical Assessment of Chronic Pacing Leads

Test/Observation	Possible Insulation Failure	Possible Conductor Failure	Possible Other System Failure	Effect on Test/Observation
<b>Pacing Impedance</b> (Telemetered or Measured Invasively)	Sudden and Significant Decrease	Sudden and Significant Increase	Dislodgement. . . . . Perforation. . . . . Electrolyte Imbalance. . . . . Improper IPG/Lead Connection. . .	Decrease Increase or Decrease Increase or Decrease Increase or Decrease
<b>Pacing Thresholds</b> (Telemetered/Programmed or Measured Invasively)	Sudden and Significant Increase, Especially in Bipolar System	Sudden and Significant Increase	Dislodgement. . . . . Exit Block. . . . . Infarct at Electrode Site. . . . . Perforation. . . . . Improper IPG/Lead Connection. . .	Increase Increase Increase Increase Increase
<b>Electrograms</b> (Telemetered or Measured Invasively)	Sudden and Significant Decrease in Amplitudes and/or Slew Rates for P and/or R Waves	Sudden and Significant Decrease or Disappearance of Amplitudes and/or Slew Rates for P and/or R Waves	Dislodgement. . . . . Perforation . . . . . Infarct at Electrode Site. . . . . Electrolyte Imbalance. . . . . Improper IPG/Lead Connection. . .	Decrease Decrease Decrease Decrease Decrease
<b>Waveform Analysis</b> (Oscillographs of Pacer Artifact from ECG Electrodes)	Sudden Increase in Ratios of Leading-Edge Voltages to Trailing-Edge Voltages (i.e., over 25% increase)	Intermittent or No Pacer Artifacts (Even in Asynchronous Mode)	Improper IPG/Lead Connection	Intermittent or No Pacer Artifacts (Even in Asynchronous Mode)
<b>Radiographs</b> (Post-Implant, Recent, Current)	Not Discernible	Visual Observation of Conductor/Connector/ Electrode Fracture (Sometimes Discernible)	Dislodgement or Perforation. Improper IPG/Lead Connection.	Sometimes Discernible
<b>Visual Inspection</b> (Invasive)	Insulation Breach and/or Degradation, or Ligature Cut-Through	Not Easily Discernible	Connector Defect or Connector Pulled Apart. Improper IPG/Lead Connection.	Sometimes Discernible
<b>Pectoral Muscle Stimulation</b>	Sudden Onset, Especially in Bipolar System		Connector Defect in Bipolar or Unipolar. Hypersensitivity to Unipolar Pulse Generator Can. Anti-Stim Coating or Protection Deficient.	
<b>Phrenic Nerve/ Diaphragmatic Stimulation</b>	Sudden Onset in Bipolar or Unipolar Systems		Perforation or Displacement of Atrial Lead (Phrenic Nerve)	
<b>Pacemaker ECG Stimulus</b> Artifact Size and Morphology Change (May Not Be Possible with Digital ECG)	Sudden Onset and Significant Change, Especially in Bipolar System (Increase in Size)	Sudden Changes, Usually a Decrease in Size	Perforation or Dislodgement. Connector Defect. Improper IPG/Lead Connection.	Sometimes Discernible
<b>Oversensing</b> (Intermittent or Continuous)	Sudden Onset, Especially in Bipolar Systems		Physical Contact between the Electrode(s) on the Lead and that of Another Lead. Inappropriate IPG Parameter Setting. Improper IPG/Lead Connection.	Sometimes Discernible
<b>Undersensing</b> (Intermittent or Continuous)	Sudden Onset in Either Unipolar or Bipolar Systems	Sudden Onset in Either Unipolar or Bipolar Systems	Dislodgement or Perforation. Infarct at Electrode Site. Electrolyte Imbalance. Inappropriate IPG Parameter Setting. Improper IPG/Lead Connection.	Sometimes Discernible
<b>Loss of Capture</b>	See "Pacing Thresholds" Above	See "Pacing Thresholds" Above	See "Pacing Thresholds" Above	

## By Model Number

2187 83, 86	6917A 140, 141, 142	8085 58, 73, 78	E2DR33 44, 70, 77
2188 83, 86	6921 87, 95, 96, 97	8086 58, 73, 78	E2SR01 45, 70, 77
3830 98, 130, 134, 136	6930 87, 95, 96, 97, 151	8088 58, 73, 78	E2SR03 45, 70, 77
4003 99, 130, 134, 136	6931 88, 95, 96, 97, 151	8089 58, 73, 78	E2SR06 45, 70, 77
4003M 99, 130, 134, 136	6932 88, 95, 96, 97, 155	8158 60, 74, 79	E2VDD01 45, 70, 77
4004 100, 130, 134, 136, 157	6933 89, 95, 96, 97	8160 60, 74, 79	KD701 49, 71, 78, 154
4004M 100, 130, 134, 136, 157	6934S 97, 155	8161 60, 74, 79	KD703 49, 71, 78, 154
4011 100, 130, 134, 136	6935 89, 95, 96, 97	8162 60, 74, 79	KD706 49, 71, 78, 154
4012 101, 130, 134, 136, 158	6936 90, 95, 96, 97, 155	8164 60, 74, 79	KDR401 47, 70, 77
4023 101, 130, 134, 136	6937 89, 95, 96, 97	8165 60, 74, 79	KDR403 47, 70, 77
4024 102, 130, 134, 136	6937A 89, 95, 96, 97	8166 60, 74, 79	KDR601 48, 71, 77, 154
4033 102, 130, 134, 136	6939 90, 95, 96, 97	8330 56, 73, 78, 159	KDR603 48, 71, 77, 154
4057 103, 130, 134, 136	6940 126, 133, 135, 137	8331 56, 73, 78, 159	KDR606 48, 71, 77, 154
4057M 103, 130, 134, 136	6942 91, 95, 96, 97, 155	8331M 56, 73, 78, 159	KDR651 49, 71, 77, 154
4058 104, 130, 134, 136	6943 91, 95, 96, 97, 155	8340 56, 73, 78, 159	KDR653 49, 71, 77, 154
4058M 104, 130, 134, 136	6944 92, 95, 96, 97	8341 56, 73, 78, 159	KDR701 50, 71, 78, 154
4067 105, 130, 134, 136	6945 92, 95, 96, 97, 155	8341M 56, 73, 78, 159	KDR703 50, 71, 78, 154
4068 106, 130, 134, 136	6947 93, 95, 96, 97	8342 56, 73, 78, 159	KDR706 50, 71, 78, 154
4073 107, 131, 134, 136	6948 93, 95, 96, 97, 151	8416 55, 72, 78	KDR721 51, 71, 78, 154
4074 108, 131, 134, 136	6949 94, 95, 96, 97, 151	8417 55, 72, 78	KDR801 53, 72, 78
4076 109, 131, 134, 136	6957 127, 133, 135, 137	8417M 55, 72, 78	KDR803 53, 72, 78
4081 110, 131, 134, 136	6957J 128, 133, 135, 137	8418 55, 72, 78	KDR901 53, 72, 78
4082 157	6961 128, 133, 135, 137	8419 55, 72, 78	KDR903 53, 72, 78
4092 110, 131, 134, 136	6962 129, 133, 135, 137	8424 56, 72, 78	KDR906 53, 72, 78
4193 84, 86	6963 89, 95, 96, 97	8426 56, 72, 78	KDR921 54, 72, 78
4194 84, 86	6966 90, 95, 96, 97, 155	8427 56, 72, 78	KSR401 48, 70, 77
4195 85, 86	6996 94, 95, 96, 97	8960i 67, 76, 79	KSR403 48, 70, 77
4196 85, 86	6999 90, 95, 96, 97	8961i 67, 76, 79	KSR701 52, 71, 78
4503 111, 131, 134, 136	7088 57, 73, 78	8962i 67, 76, 79	KSR703 52, 71, 78
4503M 111, 131, 134, 136	7089 57, 73, 78	8964i 66, 75, 79	KSR706 52, 71, 78
4504 111, 131, 134, 136, 156	7107 57, 73, 78	8965i 66, 75, 79	KSR901 54, 72, 78
4504M 111, 131, 134, 136, 156	7108 57, 73, 78	8966i 66, 75, 79	KSR903 54, 72, 78
4512 112, 131, 134, 136	7223 37	8968i 67, 76, 79	KSR906 54, 72, 78
4523 112, 131, 134, 136	7227 23, 34, 37, 146, 155	ADD01 40, 69, 77	KVDD701 52, 71, 78, 154
4524 113, 131, 134, 136	7229 23, 34, 37, 146, 155, 165	ADDR01 40, 69, 77	KVDD901 55, 72, 78
4533 113, 131, 134, 136	7230 24, 34, 37, 146, 153	ADDR03 40, 69, 77	P1501DR 46, 70, 77
4557 114, 131, 134, 136	7231 24, 34, 37, 146, 165	ADDR06 40, 69, 77	SDR203 62, 74, 79, 152
4557M 114, 131, 134, 136	7232 25, 34, 37, 146, 153	ADDRL1 40, 69, 77	SDR303 63, 75, 79, 152
4558M 114, 131, 134, 136	7253 42, 69, 77, 163	ADDRS1 41, 69, 77	SDR306 63, 75, 79, 152
4568 115, 132, 134, 136	7271 25, 34, 37, 146	ADSR01 41, 69, 77	SED01 61, 74, 79
4574 115, 132, 134, 136	7272 13, 20, 22, 37, 147	ADSR03 41, 69, 77	SED01 61, 74, 79
4582 156	7274 26, 34, 37, 147, 153	ADSR06 41, 69, 77	SEDRL1 79
4592 116, 132, 134, 136	7275 26, 35, 37, 147, 161, 165	ADVDD01 42, 69, 77	SES01 61, 74, 79
4951 138, 141, 142	7276 27, 35, 37, 147, 165	AT501 42, 69, 77, 163	SES01 61, 74, 79
4951M 138, 141, 142	7277 14, 20, 22, 37, 147, 153	C154DWK 18, 21, 22, 149	SS103 62, 74, 79, 152
4965 138, 141, 142	7278 27, 35, 37, 147, 153	C164AWK 18, 21, 22	SS106 62, 74, 79, 152
4968 139, 141, 142	7279 148, 153	C174AWK 18, 21, 22	SSR203 63, 74, 79, 152
5023 116, 132, 134, 136	7285 148, 153	D153ATG 29, 35, 38, 149	SSR303 64, 75, 79, 152
5023M 116, 132, 134, 136	7287 37	D153DRG 29, 35, 38, 149	SSR306 64, 75, 79, 152
5024 117, 132, 134, 137	7288 28, 35, 37, 148	D153VRC 38	SVDD303 64, 75, 79, 152
5024M 117, 132, 134, 137	7289 15, 20, 22, 148, 153	D154ATG 29, 35, 38, 149	VEDR01 68, 76, 79
5026 117, 132, 134, 137	7290 28, 35, 37, 148	D154AWG 30, 36, 38, 150	
5032 143, 144	7295 22	D154DRG 29, 35, 38, 149	
5033 118, 132, 134, 137	7297 15, 20, 22, 148	D154VRC 30, 36, 38, 150	
5034 118, 132, 134, 137	7299 16, 20, 22, 149	D154VWC 31, 36, 38, 150	
5038 143, 144	7303 16, 20, 22, 149	D164AWG 30, 36, 38	
5054 119, 132, 134, 137	7304 17, 20, 22, 149	D164VWC 31, 36, 38	
5068 120, 132, 134, 137	7860 59, 73, 79	D224DRG 31, 36, 38	
5069 142	7861 59, 73, 79	D224TRK 19, 21, 22	
5071 139, 141, 142	7862 59, 73, 79	D224VRC 32, 36, 38	
5072 121, 132, 134, 137	7864 59, 73, 79	D284DRG 32, 36, 38	
5076 122, 132, 134, 137	7865 59, 73, 79	D284TRK 19, 21, 22	
5092 123, 132, 134, 137	7866 59, 73, 79	D284VRC 33, 36, 38	
5524 123, 132, 134, 137	7960i 65, 75, 79	E1DR01 43, 69, 77	
5524M 123, 132, 134, 137	7961i 65, 75, 79	E1DR03 43, 69, 77	
5534 124, 133, 134, 137	7962i 65, 75, 79	E1DR06 43, 69, 77	
5554 124, 133, 134, 137	7964i 65, 75, 79	E1DR21 43, 69, 77	
5568 125, 133, 134, 137	7965i 65, 75, 79	E2DR01 44, 70, 77	
5592 125, 133, 135, 137	7966i 65, 75, 79	E2DR03 44, 70, 77	
5594 126, 133, 135, 137	7968i 66, 75, 79	E2DR06 44, 70, 77	
6721 87, 95, 96, 97	8040 17, 21, 22	E2DR21 44, 70, 77	
6917 140, 141, 142	8042 18, 21, 22	E2DR31 44, 70, 77	

## By Family

## A

Adapta DR  
 ADD01 40, 69, 77  
 ADDR01 40, 69, 77  
 ADDR03 40, 69, 77  
 ADDR06 40, 69, 77  
 ADDR1 40, 69, 77  
 ADDRS1 41, 69, 77  
 Adapta SR  
 ADSR01 41, 69, 77  
 ADSR03 41, 69, 77  
 ADSR06 41, 69, 77  
 Adapta VDD  
 ADVDD01 42, 69, 77  
 AT500  
 7253 42, 69, 77, 163  
 AT501 42, 69, 77, 163  
 Attain  
 2187 83, 86  
 2188 83, 86  
 4193 84, 86  
 4194 84, 86  
 4196 85, 86

## C

CapSure  
 4003 99, 130, 134, 136  
 4003M 99, 130, 134, 136  
 4004 100, 130, 134, 136,  
 157  
 4004M 100, 130, 134, 136,  
 157  
 4503 111, 131, 134, 136  
 4503M 111, 131, 134, 136  
 4504 111, 131, 134, 136,  
 156  
 4504M 111, 131, 134, 136,  
 156  
 5026 117, 132, 134, 137  
 CapSureFix  
 4067 105, 130, 134, 136  
 4068 106, 130, 134, 136  
 4568 115, 132, 134, 136  
 5068 120, 132, 134, 137  
 5568 125, 133, 134, 137  
 6940 126, 133, 135, 137  
 CapSureFix Novus  
 4076 109, 131, 134, 136  
 5076 122, 132, 134, 137  
 CapSure Epi  
 4965 138, 141, 142  
 4968 139, 141, 142  
 CapSure Sense  
 4073 107, 131, 134, 136  
 4074 108, 131, 134, 136  
 4574 115, 132, 134, 136  
 CapSure SP  
 4023 101, 130, 134, 136  
 4024 102, 130, 134, 136  
 4523 112, 131, 134, 136  
 4524 113, 131, 134, 136  
 5023 116, 132, 134, 136  
 5023M 116, 132, 134, 136  
 5024 117, 132, 134, 137  
 5024M 117, 132, 134, 137  
 5524 123, 132, 134, 137  
 5524M 123, 132, 134, 137

CapSure SP Novus  
 4092 110, 131, 134, 136  
 4592 116, 132, 134, 136  
 5092 123, 132, 134, 137  
 5592 125, 133, 135, 137  
 5594 126, 133, 135, 137  
 CapSure VDD  
 5032 143, 144  
 CapSure VDD-2  
 5038 143, 144  
 CapSure Z  
 4033 102, 130, 134, 136  
 4533 113, 131, 134, 136  
 5033 118, 132, 134, 137  
 5034 118, 132, 134, 137  
 5534 124, 133, 134, 137  
 CapSure Z Novus  
 5054 119, 132, 134, 137  
 5554 124, 133, 134, 137  
 Concerto  
 C154DWK 18, 21, 22, 149  
 C164AWK 18, 21, 22  
 C174AWK 18, 21, 22  
 Consulta  
 D224TRK 19, 21, 22

## E

EnPulse DR  
 E1DR01 43, 69, 77  
 E1DR03 43, 69, 77  
 E1DR06 43, 69, 77  
 E1DR21 43, 69, 77  
 EnPulse 2 DR  
 E2DR01 44, 70, 77  
 E2DR03 44, 70, 77  
 E2DR06 44, 70, 77  
 E2DR21 44, 70, 77  
 E2DR31 44, 70, 77  
 E2DR33 44, 70, 77  
 EnPulse 2 SR  
 E2SR01 45, 70, 77  
 E2SR03 45, 70, 77  
 E2SR06 45, 70, 77  
 EnPulse 2 VDD  
 E2VDD01 45, 70, 77  
 EnRhythm DR  
 P1501DR 46, 70, 77  
 EnTrust  
 D153ATG 29, 35, 38, 149  
 D153DRG 29, 35, 38, 149  
 D153VRC 38  
 D154ATG 29, 35, 38, 149  
 D154DRG 29, 35, 38, 149  
 D154VRC 30, 36, 38, 150  
 Epicardial Patch  
 6721 87, 95, 96, 97  
 6921 87, 95, 96, 97

## G

GEM  
 7227 23, 34, 37, 146, 155  
 GEM DR  
 7271 25, 34, 37, 146  
 GEM II VR  
 7229 23, 34, 37, 146, 155,  
 165  
 GEM III AT  
 7276 27, 35, 37, 147, 165  
 GEM III DR  
 7275 26, 35, 37, 147, 161,  
 165  
 GEM III VR  
 7231 24, 34, 37, 146, 165

## I

InSync  
 8040 17, 21, 22  
 InSync ICD  
 7272 13, 20, 22, 37, 147  
 InSync Marquis  
 7277 14, 20, 22, 147, 153  
 InSync Maximo  
 7303 16, 20, 22, 149  
 7304 17, 20, 22, 149  
 InSync Sentry  
 7297 15, 20, 22, 148  
 7299 16, 20, 22, 149  
 InSync II Marquis  
 7289 15, 20, 22, 148, 153  
 InSync II Protect  
 7295 22  
 InSync III Protect  
 7285 148, 153  
 InSync III  
 8042 18, 21, 22  
 InSync III Marquis  
 7279 148, 153  
 Intrinsic  
 7288 28, 35, 37, 148  
 Intrinsic 30  
 7287 37

## K

Kappa 400 DR  
 KDR401 47, 70, 77  
 KDR403 47, 70, 77  
 Kappa 400 SR  
 KSR401 48, 70, 77  
 KSR403 48, 70, 77  
 Kappa 600 DR  
 KDR601 48, 71, 77, 154  
 KDR603 48, 71, 77, 154  
 KDR606 48, 71, 77, 154  
 KDR651 49, 71, 77, 154  
 KDR653 49, 71, 77, 154

Kappa 700 D  
 KD701 49, 71, 78, 154  
 KD703 49, 71, 78, 154  
 KD706 49, 71, 78, 154  
 Kappa 700 DR  
 KDR701 50, 71, 78, 154  
 KDR703 50, 71, 78, 154  
 KDR706 50, 71, 78, 154  
 KDR721 51, 71, 78, 154  
 Kappa 700 SR  
 KSR701 52, 71, 78  
 KSR703 52, 71, 78  
 KSR706 52, 71, 78  
 Kappa 700 VDD  
 KVDD701 52, 71, 78, 154  
 Kappa 800 DR  
 KDR801 53, 72, 78  
 KDR803 53, 72, 78  
 Kappa 900 DR  
 KDR901 53, 72, 78  
 KDR903 53, 72, 78  
 KDR906 53, 72, 78  
 Kappa 900 SR  
 KSR901 54, 72, 78  
 KSR903 54, 72, 78  
 KSR906 54, 72, 78  
 Kappa 900 VDD  
 KVDD901 55, 72, 78  
 Kappa 920 DR  
 KDR921 54, 72, 78

## L

Legend  
 8416 55, 72, 78  
 8417 55, 72, 78  
 8417M 55, 72, 78  
 8418 55, 72, 78  
 8419 55, 72, 78  
 Legend II  
 8424 56, 72, 78  
 8426 56, 72, 78  
 8427 56, 72, 78

## M

Marquis DR  
 7274 26, 34, 37, 147, 153  
 Marquis VR  
 7230 24, 34, 37, 146, 153  
 Maximo DR  
 7278 27, 35, 37, 147, 153  
 Maximo II  
 D284TRK 19, 21, 22  
 Maximo II DR  
 D284DRG 32, 36, 38  
 Maximo II VR  
 D284VRC 33, 36, 38  
 Maximo VR  
 7232 25, 34, 37, 146, 153  
 Minix/Minix ST  
 8330 56, 73, 78, 159  
 8331 56, 73, 78, 159  
 8331M 56, 73, 78, 159  
 8340 56, 73, 78, 159  
 8341 56, 73, 78, 159  
 8341M 56, 73, 78, 159  
 8342 56, 73, 78, 159

By Family continued

Minuet  
7107 57, 73, 78  
7108 57, 73, 78

O

Onyx  
7290 28, 35, 37, 148

P

Preva DR  
7088 57, 73, 78  
7089 57, 73, 78  
Preva SR  
8088 58, 73, 78  
8089 58, 73, 78  
Prevail S  
8085 58, 73, 78  
8086 58, 73, 78

Prodigy D  
7864 59, 73, 79  
7865 59, 73, 79  
7866 59, 73, 79

Prodigy DR  
7860 59, 73, 79  
7861 59, 73, 79  
7862 59, 73, 79

Prodigy S  
8164 60, 74, 79  
8165 60, 74, 79  
8166 60, 74, 79

Prodigy SR  
8158 60, 74, 79  
8160 60, 74, 79  
8161 60, 74, 79  
8162 60, 74, 79

S

Screw-In  
4057 103, 130, 134, 136  
4057M 103, 130, 134, 136  
4058 104, 130, 134, 136  
4058M 104, 130, 134, 136  
4557 114, 131, 134, 136  
4557M 114, 131, 134, 136  
4558M 114, 131, 134, 136

Secura DR  
D224DRG 31, 36, 38

SelectSecure  
3830 98, 130, 134, 136

Sensia DR  
SED01 61, 74, 79  
SEDR01 61, 74, 79  
SEDR1 79

Sensia SR  
SES01 61, 74, 79  
SESR01 61, 74, 79

Sigma 100 S  
SS103 62, 74, 79, 152  
SS106 62, 74, 79, 152  
Sigma 200 DR  
SDR203 62, 74, 79, 152

Sigma 200 SR  
SSR203 63, 74, 79, 152

Sigma 300 DR  
SDR303 63, 75, 79, 152  
SDR306 63, 75, 79, 152

Sigma 300 SR  
SSR303 64, 75, 79, 152  
SSR306 64, 75, 79, 152

Sigma 300 VDD  
SVDD303 64, 75, 79, 152

Spectraflex  
4951 138, 141, 142  
4951M 138, 141, 142  
6957 127, 133, 135, 137  
6957J 128, 133, 135, 137

Sprint  
6932 88, 95, 96, 97  
6942 91, 95, 96, 97  
6943 91, 95, 96, 97  
6945 92, 95, 96, 97

Sprint Fidelis  
6930 87, 95, 96, 97, 151  
6931 88, 95, 96, 97, 151  
6948 93, 95, 96, 97, 151  
6949 94, 95, 96, 97, 151

Sprint Quattro  
6944 92, 95, 96, 97

Sprint Quattro Secure  
6935 89, 95, 96, 97  
6947 93, 95, 96, 97

StarFix  
4195 85, 86

Sub-Q Lead  
6996 94, 95, 96, 97

Sub-Q Patch  
6939 90, 95, 96, 97  
6999 90, 95, 96, 97

SureFix  
5072 121, 132, 134, 137

SVC/CS  
6933 89, 95, 96, 97  
6937 89, 95, 96, 97  
6937A 89, 95, 96, 97  
6963 89, 95, 96, 97

T

Target Tip  
4011 100, 130, 134, 136  
4012 101, 130, 134, 136,  
158

4081 110, 131, 134, 136  
4082 157  
4512 112, 131, 134, 136  
4582 156

Tenax  
6917 140, 141, 142  
6917A 140, 141, 142  
6961 128, 133, 135, 137  
6962 128, 133, 135, 137

Thera-i D  
7964i 65, 75, 79  
7965i 65, 75, 79  
7966i 65, 75, 79

Thera-i DR  
7960i 65, 75, 79

7961i 65, 75, 79  
7962i 65, 75, 79  
7968i 66, 75, 79

Thera-i S  
8964i 66, 75, 79  
8965i 66, 75, 79  
8966i 66, 75, 79

Thera-i SR  
8960i 67, 76, 79  
8961i 67, 76, 79  
8962i 67, 76, 79

Thera-i VDD  
8968i 67, 76, 79

Transvene  
6934S 97  
6936 90, 95, 96, 97  
6966 90, 95, 96, 97

V

Versa DR  
VEDR01 68, 76, 79

Virtuoso  
D154AWG 30, 36, 38, 150  
D164AWG 30, 36, 38  
D154VWC 31, 36, 38, 150  
D164VWC 31, 36, 38

*If you are looking for a model number or family that is not included in this report, you may call US Technical Services (see page 2).*

# Mailer Kits Available for Returning Product

Medtronic urges all physicians to return explanted products and to notify Medtronic when a product is no longer in use, regardless of reason for explant or removal from use. The procedures for returning products vary by geographic location.

Mailer kits (pictured right) with prepaid US postage are available for use within the United States to send CRT, ICD, IPG, and leads to Medtronic's CRDM Returned Product Analysis Lab. These mailers are sized to accommodate the devices and leads from a single patient or clinical event and are designed to meet postal regulations for mailing biohazard materials.

If the product being returned is located outside the United States, please contact your local Medtronic representative for instructions.

Medtronic also requests the return of devices from non-clinical sources, such as funeral homes, and will assume responsibility for storage and disposal of the product once received. For return of larger quantities of explanted products than the mailer can accommodate, Medtronic has handling and shipping guidelines available upon request.

Both mailers and guidelines can be requested by contacting the Returned Product Lab.

CRDM Returned Product Analysis Laboratory  
Medtronic, Inc.  
7000 Central Avenue NE MS RCE172  
Minneapolis, MN 55432-3576 USA  
Phone: 1 (800) 328-2518, ext. 44800  
Email: [crdm.returnedproduct@medtronic.com](mailto:crdm.returnedproduct@medtronic.com)



[www.medtronic.com](http://www.medtronic.com)

## World Headquarters

Medtronic, Inc.  
710 Medtronic Parkway  
Minneapolis, MN 55432-5604  
USA  
Tel: (763) 514-4000  
Fax: (763) 514-4879

Medtronic USA, Inc.  
Toll-free: 1 (800) 328-2518  
(24-hour technical support for  
physicians and medical professionals)

## Europe

Medtronic International Trading Sàrl  
Route du Molliau 31  
CH-1131 Tolochenaz  
Switzerland  
Tel: (41 21) 802 7000  
Fax: (41 21) 802 7900

## Canada

Medtronic of Canada Ltd.  
6733 Kitimat Road  
Mississauga, Ontario L5N 1W3  
Canada  
Tel: (905) 826-6020  
Fax: (905) 826-6620  
Toll-free: 1 (800) 268-5346

## Asia Pacific

Medtronic International, Ltd.  
16/F Manulife Plaza  
The Lee Gardens, 33 Hysan Avenue  
Causeway Bay  
Hong Kong  
Tel: (852) 2891 4456  
Fax: (852) 2891 6830  
[enquiryap@medtronic.com](mailto:enquiryap@medtronic.com)

## Latin America

Medtronic USA, Inc.  
Doral Corporate Center II  
3750 NW 87th Avenue Suite 700  
Miami, FL 33178  
USA  
Tel: (305) 500-9328  
Fax: (786) 709-4244

